Package 'ShortRead'

September 24, 2012

Type Package
Title Classes and methods for high-throughput short-read sequencing data.
Version 1.14.4
Author Martin Morgan, Michael Lawrence, Simon Anders
Maintainer Bioconductor Package Maintainer <maintainer@bioconductor.org></maintainer@bioconductor.org>
Description Base classes, functions, and methods for representation of high-throughput, short-read sequencing data.
License Artistic-2.0
LazyLoad yes
Depends methods, BiocGenerics (>= 0.1.0), IRanges (>= 1.13.4), GenomicRanges (>= 1.7.43), Biostrings (>= 2.23.2), lattice,Rsamtools (>= 1.7.42), latticeExtra
Imports BiocGenerics, IRanges, GenomicRanges, Biostrings, Biobase, hwriter, Rsamtools, zlibbioc, lattice
Enhances Rmpi, parallel
Suggests biomaRt, RUnit, GenomicFeatures, yeastNagalakshmi
LinkingTo IRanges, Biostrings
biocViews DataImport, Sequencing, HighThroughputSequencing, QualityControl
R topics documented:
ShortReadBase-package 3 .QA-class 3 accessors 4 AlignedDataFrame 5 AlignedDataFrame-class 6 AlignedRead 7 AlignedRead-class 8 alphabetByCycle 10 alphabetScore 11 BAMQA-class 12

SowtieQA-class		13
lean		14
ountLines		14
eprecated		15
ustyScore		17
xperimentPath-class		18
astqFile-class		19
SappedReads-class		21
ntensity-class		22
MAQMapQA-class		23
a		24
QA-class		25
a2		27
QualityScore		30
QualityScore-class		31
eadAligned		33
eadBaseQuality		38
eadBfaToc		39
eadFasta		40
eadFastq		41
eadIntensities		43
eadPrb		45
eadQseq		46
eadXStringColumns		47
enewable		48
eport		49
cochePath-class		51
cocheSet-class		
taIntensity		
taIntensity-class		
hortRead-class		
hortRead-deprecated		
hortReadQ-class		58
hortReadQA-class		60
napshot-class		61
napshotFunction-class		64
olexaExportQA-class		66
olexaIntensity		67
olexaIntensity class		68
olexaPath-class		69
olexaSet-class		72
pTrellis-class		73
p ViewPerFeature		74
·		76
rapply		
rdistance		77
rduplicated		78
rFilter		80
RFilter-class		83
RFilterResult-class		84
RSet-class		86
RUtil-class		87
ıbles		89

ShortRe	eadB	ase-j	paci	kage

3

	Utilites	92
Index	•	94

ShortReadBase-package Base classes and methods for high-throughput short-read sequencing data.

Description

Base classes, functions, and methods for representation of high-throughput, short-read sequencing data.

Details

See packageDescription('ShortRead')

Author(s)

Maintainer: Martin Morgan <mtmorgan@fhcrc.org>

.QA-class

Virtual class for representing quality assessment results

Description

Classes derived from .QA-class represent results of quality assurance analyses. Details of derived class structure are found on the help pages of the derived classes.

Objects from the Class

Objects from the class are created by ShortRead functions, in particular qa.

Extends

```
Class ". ShortReadBase", directly.
```

Methods

Methods defined on this class include:

```
rbind signature(...="list"): rbind data frame objects in .... All objects of ... must be of the same class; the return value is an instance of that class.
```

show signature(object = "SolexaExportQA"): Display an overview of the object contents.

Author(s)

Martin Morgan <mtmmorgan@fhcrc.org>

4 accessors

See Also

Specific classes derived from .QA

Examples

```
getClass(".QA", where=getNamespace("ShortRead"))
```

accessors

Accessors for ShortRead classes

Description

These functions and generics define 'accessors' (to get and set values) for objects in the **ShortRead** package; methods defined in other packages may have additional meaning.

Usage

```
## SRVector
vclass(object, ...)
## ShortRead / ShortReadQ
sread(object, ...)
id(object, ...)
## AlignedRead
chromosome(object, ...)
position(object, ...)
alignQuality(object, ...)
alignData(object, ...)
## Solexa
experimentPath(object, ...)
dataPath(object, ...)
scanPath(object, ...)
imageAnalysisPath(object, ...)
baseCallPath(object, ...)
analysisPath(object, ...)
## SolexaSet
solexaPath(object, ...)
laneDescription(object, ...)
laneNames(object, ...)
```

Arguments

. . .

object An object derived from class ShortRead. See help pages for individual objects, e.g., ShortReadQ. The default is to extract the contents of a slot of the corre-

sponding name (e.g., slot sread) from object.

Additional arguments passed to the accessor. The default definitions do not

make use of additional arguments.

AlignedDataFrame 5

Value

Usually, the value of the corresponding slot, or other simple content described on the help page of object.

Author(s)

Martin Morgan

Examples

```
sp <- SolexaPath(system.file('extdata', package='ShortRead'))
experimentPath(sp)
basename(analysisPath(sp))</pre>
```

AlignedDataFrame

AlignedDataFrame constructor

Description

Construct an AlignedDataFrame from a data frame and its metadata

Usage

```
AlignedDataFrame(data, metadata, nrow = nrow(data))
```

Arguments

data A data frame containing alignment information.

metadata A data frame describing the columns of data, and with number of rows of

metadata corresponding to number of columns of data. The data frame must contain a column labelDescription providing a verbose description of each

column of data.

nrow An optional argument, to be used when data is not provided, to construct an

AlignedDataFrame with the specified number of rows.

Value

An object of AlignedDataFrame.

Author(s)

Martin Morgan mtmorgan@fhcrc.org

AlignedDataFrame-class

"AlignedDataFrame" representing alignment annotations as a data frame

Description

This class extends AnnotatedDataFrame. It is a data frame and associated metadata (describing the columns of the data frame). The main purpose of this class is to contain alignment data in addition to the central information of AlignedRead.

Objects from the Class

Objects from the class are created by calls to the AlignedDataFrame function.

Slots

data: Object of class "data.frame" containing the data. See AnnotatedDataFrame for details.

varMetadata: Object of class "data.frame" describing columns of data. See AnnotatedDataFrame for details.

dimLabels: Object of class character describing the dimensions of the AnnotatedDataFrame. Used internally; see AnnotatedDataFrame for details.

.__classVersion__: Object of class "Versions" describing the version of this object. Used internally; see AnnotatedDataFrame for details.

Extends

Class "AnnotatedDataFrame", directly. Class "Versioned", by class "AnnotatedDataFrame", distance 2.

Methods

This class inherits methods pData (to retrieve the underlying data frame) and varMetadata (to retrieve the metadata) from AnnotatedDataFrame.

Additional methods include:

append signature(x = "AlignedDataFrame", values = "AlignedDataFrame", length = "missing"):
 append values after x. varMetadata of x and y must be identical; pData and varMetadata
 are appended using rbind.

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

See Also

AnnotatedDataFrame

AlignedRead 7

AlignedRead Construct objects of class "A	AlignedRead"
---	--------------

Description

This function constructs objects of AlignedRead. It will often be more convenient to create AlignedRead objects using parsers such as readAligned.

Usage

Arguments

sread	An object of class DNAStringSet, containing the DNA sequences of the short reads.
id	An object of class BStringSet, containing the identifiers of the short reads. This object is the same length as sread.
quality	An object of class BStringSet, containing the ASCII-encoded quality scores of the short reads. This object is the same length as sread.
chromosome	A factor describing the particular sequence within a set of target sequences (e.g. chromosomes in a genome assembly) to which each short read aligns.
position	A integer vector describing the (base pair) position at which each short read begins its alignment.
strand	A factor describing the strand to which the short read aligns.
alignQuality	A numeric vector describing the alignment quality.
alignData	An AlignedDataFrame with number of rows equal to the length of sread, containing additional information about alignments.

Value

An object of class AlignedRead.

Author(s)

Martin Morgan mtmorgan@fhcrc.org

See Also

AlignedRead.

AlignedRead-class

"AlignedRead" class for aligned short reads

Description

This class represents and manipulates reads and their genomic alignments. Alignment information includes genomic position, strand, quality, and other data.

Objects from the Class

Objects of this class can be created from a call to the AlignedRead constructor, or more typically by parsing appropriate files (e.g., readAligned).

Slots

chromosome Object of class "factor" the particular sequence within a set of target sequences (e.g. chromosomes in a genome assembly) to which each short read aligns.

position Object of class "integer" the (base-pair) position in the genome to which the read is aligned. AlignedRead objects created by readAligned use 1-based indexing, with alignemnts reported in 'left-most' coordinates, as described in the vignette.

strand Object of class "factor" the strand of the alignment.

alignQuality Object of class "numeric" representing an alignment quality score.

alignData Object of class "AlignedDataFrame" additional alignment information.

quality Object of class "BStringSet" representing base-call read quality scores.

sread Object of class "DNAStringSet" DNA sequence of the read.

id Object of class "BStringSet" read identifier.

Extends

```
Class "ShortReadQ", directly. Class "ShortRead", by class "ShortReadQ", distance 2. Class ".ShortReadBase", by class "ShortReadQ", distance 3.
```

Methods

See accessors for additional functions to access slot content, and ShortReadQ, ShortRead for inherited methods. Additional methods include:

```
[ signature(x = "AlignedRead", i = "ANY", j = "missing"): This method creates a new AlignedRead object containing only those reads indexed by i. chromosome is recoded to contain only those levels in the new subset.
```

```
append signature(x = "AlignedRead", values = "AlignedRead", length = "missing"):
    append values after x. chromosome and strand must be factors with the same levels. See
    methods for ShortReadQ, AlignedDataFrame for details of how these components of x and y
    are appended.
```

```
coerce signature(from = "PairwiseAlignedXStringSet", to = "AlignedRead"):
    signature(from = "AlignedRead", to = "RangesList"): signature(from = "AlignedRead", to = "RangesList"): signature(from = "AlignedRead", to = "Gapped signature(from = "AlignedRead", to = "GappedReads"):
```

AlignedRead-class 9

Invoke these methods with, e.g., as(from, "AlignedRead") to coerce objects of class from to class "AlignedRead".

Coercion from AlignedRead to RangesList, RangedData or GRanges assumes that position(from) uses a 'leftmost' (see coverage on this page) coordinate system. Since Ranges objects cannot store NA values, reads with NA in the position, width, chromosome or (in the case of GRanges) strand vectors are dropped.

chromosome signature(object = "AlignedRead"): access the chromosome slot of object.

position signature(object = "AlignedRead"): access the position slot of object.

strand signature(object = "AlignedRead"): access the strand slot of object.

coverage signature(x = "AlignedRead", shift = OL, width = NULL, weight = 1L, ..., coords = c("le Calculate coverage across reads present in x.

shift must be either 0L or a named integer vector with names including all levels(chromosome(x)). It specifies how the reads in x should be (horizontally) shifted *before* the coverage is computed.

width must be either NULL or a named vector of non-negative integers with names including all levels(chromosome(x)). In the latter case, it specifies for each chromosome the end of the chromosome region over which coverage is to be calculated *after* the reads have been

shifted. Note that this region always starts at chromosome position 1. If width is NULL, it ends at the rightmost chromosome position covered by at least one read.

weight must be 1L for now (weighting the reads is not supported yet, sorry).

coords specifies the coordinate system used to record position. Both systems number base pairs from left to right on the 5' strand. leftmost indicates the eland convention, where position(x) is the left-most (minimum) base pair, regardless of strand. fiveprime is the MAQ convention, where position(x) is the coordinate of the 5' end of the aligned read.

extend indicates the number of base pairs to extend the read. Extension is in the 3' direction, measured from the 3' end of the aligned read.

The return value of coverage is a SimpleRleList object.

```
%in% signature(x = "AlignedRead", table = "RangesList"):
```

Return a length(x) logical vector indicating whether the chromosome, position, and width of x overlap (see IRanges overlap) with ranges in table. Reads for which chromosome(), position(), or width() return NA *never* overlap with table. This function assumes that positions are in 'leftmost' coordinates, as defined in coverage.

```
srorder signature(x = "AlignedRead", ..., withSread=TRUE):
srrank signature(x = "AlignedRead", ..., withSread=TRUE):
srsort signature(x = "AlignedRead", ..., withSread=TRUE):
srduplicated signature(x = "AlignedRead", ..., withSread=TRUE):
```

Order, rank, sort, and find duplicates in AlignedRead objects. Reads are sorted by chromosome, strand, position, and then (if withSread=TRUE) sread; less fine-grained sorting can be accomplished with, e.g., x[srorder(sread(x))]. srduplicated behaves like duplicated, i.e., the first copy of a duplicate is FALSE while the remaining copies are TRUE.

show signature(object = "AlignedRead"): provide a compact display of the AlignedRead
 content.

detail signature(x = "AlignedRead"): display alignData in more detail.

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

10 alphabetByCycle

See Also

```
readAligned
```

Examples

```
showMethods(class="AlignedRead", where=getNamespace("ShortRead"))
dirPath <- system.file('extdata', 'maq', package='ShortRead')
(aln <- readAligned(dirPath, 'out.aln.1.txt', type="MAQMapview"))
coverage(aln)[[1]]
cvg <- coverage(aln, shift=c(ChrA=10L))
## remove 0 coverage on left ends
ltrim0 <- function(x) {
    i <- !cumprod(runValue(x) == 0)
        Rle(runValue(x)[i], runLength(x)[i])
}
endoapply(cvg, ltrim0)
## demonstration of show() and detail() methods
show(aln)
detail(aln)</pre>
```

alphabetByCycle

Summarize short read nucleotide or quality scores by cycle

Description

alphabetByCycle summarizes short read nucleotides or qualities by cycle, e.g., returning the number of occurrences of each nucleotide A, T, G, C across all reads from 36 cycles of a Solexa lane.

Usage

```
alphabetByCycle(stringSet, alphabet, ...)
```

Arguments

stringSet A R object representing the collection of reads or quality scores to be summarized.

alphabet The alphabet (character vector of length 1 strings) from which the sequences in stringSet are composed. Methods often define an appropriate alphabet, so that the user does not have to provide one.

Additional arguments, perhaps used by methods defined on this generic.

Details

The default method requires that stringSet extends the XStringSet class of Biostrings.

The following method is defined, in addition to methods described in class-specific documentation:

alphabetByCycle signature(stringSet = "BStringSet"): this method uses an alphabet spanning all ASCII characters, codes 1:255. alphabetScore 11

Value

A matrix with number of rows equal to the length of alphabet and columns equal to the maximum width of reads or quality scores in the string set. Entries in the matrix are the number of times, over all reads of the set, that the corresponding letter of the alphabet (row) appeared at the specified cycle (column).

Author(s)

Martin Morgan

See Also

The IUPAC alphabet in Biostrings.

http://www.bioperl.org/wiki/FASTQ_sequence_format for the BioPerl definition of fastq. Solexa documentation 'Data analysis - documentation : Pipeline output and visualisation'.

Examples

```
showMethods("alphabetByCycle")

sp <- SolexaPath(system.file('extdata', package='ShortRead'))
rfq <- readFastq(analysisPath(sp), pattern="s_1_sequence.txt")
alphabetByCycle(sread(rfq))

abcq <- alphabetByCycle(quality(rfq))
dim(abcq)
## 'high' scores, first and last cycles
abcq[64:94,c(1:5, 32:36)]</pre>
```

alphabetScore

Efficiently calculate the sum of quality scores across bases

Description

This generic takes a QualityScore object and calculates, for each read, the sum of the encoded nucleotide probabilities.

Usage

```
alphabetScore(object, ...)
```

Arguments

object An object of class QualityScore.
... Additional arguments, currently unused.

Value

A vector of numeric values of length equal to the length of object.

Author(s)

Martin Morgan mtmorgan@fhcrc.org

12 BAMQA-class

BAMQA-class

Quality assessment from BAM files

Description

This class contains a list-like structure with summary descriptions derived from visiting one or more BAM files.

Objects from the Class

Objects of the class are usually produced by a qa method, with the argument type="BAM".

Slots

.srlist: Object of class "list", containing data frames or lists of data frames summarizing the results of qa.

Extends

```
Class "SRList", directly. Class ".QA", directly. Class ".SRUtil", by class "SRList", distance 2. Class ".ShortReadBase", by class ".QA", distance 2.
```

Methods

Accessor methods are inherited from the SRList class.

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

See Also

qa.

```
showClass("BAMQA")
```

BowtieQA-class 13

BowtieQA-class

Quality assessment summaries from Bowtie files

Description

This class contains a list-like structure with summary descriptions derived from visiting one or more Bowtie files.

Objects from the Class

Objects of the class are usually produced by a qa method, with the argument type="Bowtie".

Slots

.srlist: Object of class "list", containing data frames or lists of data frames summarizing the results of qa.

Extends

```
Class "SRList", directly. Class ".QA", directly. Class ".SRUtil", by class "SRList", distance 2. Class ".ShortReadBase", by class ".QA", distance 2.
```

Methods

Accessor methods are inherited from the SRList class.

```
report signature(x="BowtieQA", ..., dest=tempfile(), type="html"): produces
an html file summarizing the QA results.
```

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

See Also

qa.

```
showClass("BowtieQA")
```

14 countLines

Country	clo	ean	Remove classes	sequences	with	ambiguous	nucleotides	from	short	read
---------	-----	-----	-------------------	-----------	------	-----------	-------------	------	-------	------

Description

Short reads may contain ambiguous base calls (i.e., IUPAC symbols different from A, T, G, C). This generic removes all sequences containing 1 or more ambiguous bases.

Usage

```
clean(object, ...)
```

Arguments

object An object for which clean methods exist; see below to discover these methods.

Additional arguments, perhaps used by methods.

Details

The following method is defined, in addition to methods described in class-specific documentation:

Value

An instance of class(object), containing only sequences with non-redundant nucleotides.

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

Examples

```
showMethods('clean')
```

countLines

Count lines in all (text) files in a directory whose file name matches a pattern

Description

countLines visits all files in a directory path dirPath whose base (i.e., file) name matches pattern. Lines in the file are counted as the number of new line characters.

Usage

```
countLines(dirPath, pattern=character(0), ..., useFullName=FALSE)
```

deprecated 15

Arguments

dirPath	A character vector (or other object; see methods defined on this generic) giving the directory path (relative or absolute) of files whose lines are to be counted.
pattern	The (grep-style) pattern describing files whose lines are to be counted. The default (character(0)) results in line counts for all files in the directory.
	Additional arguments, passed internally to list.files. See list.files.
useFullName	A logical(1) indicating whether elements of the returned vector should be named with the base (file) name (default; useFullName=FALSE) or the full path name (useFullName=TRUE).

Value

A named integer vector of line counts. Names are paths to the files whose lines have been counted, excluding dirPath.

Author(s)

Martin Morgan

Examples

```
sp <- SolexaPath(system.file('extdata', package='ShortRead'))
countLines(analysisPath(sp))
countLines(experimentPath(sp), recursive=TRUE)
countLines(experimentPath(sp), recursive=TRUE, useFullName=TRUE)</pre>
```

deprecated

Deprecated and defunct functions

Description

These functions were introduced but are now deprecated or defunct.

Usage

Arguments

object	For basePath, and object of class ExperimentPath.
	Additional arguments.
start	A vector with the start positions of each read on the reference sequence. All reads must correspond to the same reference sequence.
fraglength	A vector of the same length as 'start' with the lengths of all the fragments. Alternatively, a single integer, specifying one constant length to assume for all tags.
chrlength	The length of the reference sequence. You may use the function readBfaToc to extract this information from the .bfa file.

16 deprecated

dir A factor with level "-" and "+" of the same length as 'start', specifying whether

the fragment extends to the right (towards higher index values, '+') or to the left (towards lower index values, '-') beyond the read. See below for more explana-

tion.

readlength The length of the reads, either as a vector of the same length as 'start' or as a

single number. This parameter makes sense only if 'dir' is used, too. If not

specified, read lengths and fragment lengths are taken to be the same.

offset The index of the first base pair in the result vector. The default is 1, i.e. assumes

that the 'start' positions are in 1-based chromosome coordinates.

Value

pileup an integer vector of length 'chrlength', each element counting how many frag-

ments map to this basepair.

Note

(the following refers to the pileup function)

- 1. This function is not suitable for paired-end reads.
- 2. If the arguments 'dir' and 'readlength' are not used, the fragments are assumed to start at the positions given in 'start' and extend to the right by the number of basepairs given in fraglength. If 'dir' and 'readlength' are supplied then the interval starting at 'start' and extending to the right by the number of base pairs given in 'readlength' marks the position of the read, which is one end of the fragment. If 'dir' ist '+', it is taken as the left end and the fragment will be extended to the right to have the total length given by 'fraglength'. If 'dir' is '-', the end is taken as the right end and is extended to the left. Note that in the latter case, the 'start' position does mark the border between read and rest of fragment, not an actual 'end' of the fragment. If you are confused now, look at the examples below.
- 3. Sorry for the inconsequent use of 'width' and 'length' in a seemingly interchangeable fashion.

Author(s)

Simon Anders, EMBL-EBI, <sanders@fs.tum.de>

dustyScore 17

```
## End(Not run)
```

dustyScore

Summarize low-complexity sequences

Description

dustyScore identifies low-complexity sequences, in a manner inspired by the dust implementation in BLAST.

Usage

```
dustyScore(x, batchSize=NA, ...)
```

Arguments

X	A DNAStringSet object, or object derived from ShortRead, containing a collection of reads to be summarized.
batchSize	NA or an integer(1) vector indicating the maximum number of reads to be processed at any one time.
	Additional arguments, not currently used.

Details

The following methods are defined:

```
dustyScore signature(x = "DNAStringSet"): operating on an object derived from class DNAStringSet.
dustyScore signature(x = "ShortRead"): operating on the sread of an object derived from class ShortRead.
```

The dust-like calculations used here are as implemented at https://stat.ethz.ch/pipermail/bioc-sig-sequencing/2009-February/000170.html. Scores range from 0 (all triplets unique) to the square of the width of the longest sequence (poly-A, -C, -G, or -T).

The batchSize argument can be used to reduce the memory requirements of the algorithm by processing the x argument in batches of the specified size. Smaller batch sizes use less memory, but are computationally less efficient.

Value

A vector of numeric scores, with length equal to the length of x.

Author(s)

```
Herve Pages (code); Martin Morgan
```

References

Morgulis, Getz, Schaffer and Agarwala, 2006. WindowMasker: window-based masker for sequenced genomes, Bioinformatics 22: 134-141.

18 ExperimentPath-class

See Also

The WindowMasker supplement defining dust ftp://ftp.ncbi.nlm.nih.gov/pub/agarwala/windowmasker/windowmasker_suppl.pdf

Examples

```
sp <- SolexaPath(system.file('extdata', package='ShortRead'))
rfq <- readFastq(analysisPath(sp), pattern="s_1_sequence.txt")
range(dustyScore(rfq))</pre>
```

ExperimentPath-class "ExperimentPath" class representing a file hierarchy of data files

Description

Short read technologies often produce a hierarchy of output files. The content of the hierarchy varies. This class represents the root of the file hierarchy. Specific classes (e.g., SolexaPath) represent different technologies.

Objects from the Class

Objects from the class are created by calls to the constructor:

ExperimentPath(experimentPath)

experimentPath character(1) object pointing to the top-level directory of the experiment; see specific technology classes for additional detail.

verbose=FALSE (optional) logical vector which, when TRUE results in warnings if paths do not exist.

All paths must be fully-specified.

Slots

ExperimentPath has one slot, containing a fully specified path to the corresponding directory (described above).

basePath See above.

The slot is accessed with experimentPath.

Extends

```
Class ". ShortReadBase", directly.
```

Methods

Methods include:

```
show signature(object = "ExperimentPath"): briefly summarize the file paths of object.
detail signature(x = "ExperimentPath"): summarize file paths of x.
```

FastqFile-class 19

Author(s)

Michael Lawrence

Examples

```
showClass("ExperimentPath")
```

FastqFile-class

Sampling and streaming records from fastq files

Description

FastqFile represents a path and connection to a fastq file. FastqFileList is a list of such connections.

FastqSampler draws a subsample from a fastq file. yield is the method used to extract the sample from the FastqSampler instance; a short illustration is in the example below. FastqSamplerList is a list of FastqSampler elements.

FastqStreamer draws successive subsets from a fastq file. Iterating over the stream will eventually yield the entire file; a short illustration is in the example below. FastqStreamerList is a list of FastqStreamer elements.

Usage

```
## FastqFile and FastqFileList
FastqFile(con, ...)
FastqFileList(..., class="FastqFile")
## S3 method for class 'ShortReadFile'
open(con, ...)
## S3 method for class 'ShortReadFile'
close(con, ...)
## S4 method for signature 'FastqFile'
readFastq(dirPath, pattern=character(), ...)
## FastqSampler and FastqStreamer
FastqSampler(con, n=1e6, readerBlockSize=1e8, verbose=FALSE)
FastqSamplerList(..., n=1e6, readerBlockSize=1e8, verbose=FALSE)
FastqStreamer(con, n=1e6, readerBlockSize=1e8, verbose=FALSE)
FastqStreamerList(..., n=1e6, readerBlockSize=1e8, verbose=FALSE)
yield(x, ...)
```

Arguments

 $\hbox{con, dirPath} \qquad \hbox{A character string naming a connection, or (for con) an R connection (e.g., \verb|file|, or con) an R connection (e.g., \verb|file|, or con) and R connection (e.g., or con) and R connection (e.g.,$

gzfile).

The size of the sample (number of records) to be drawn.

readerBlockSize

The number of bytes or characters to be read at one time; smaller readerBlockSize

reduces memory requirements but is less efficient.

verbose Display progress.

20 FastqFile-class

An instance from the FastqSampler or FastqStreamer class.

Additional arguments. For FastqFileList, FastqSamplerList, or FastqStreamerList, this can either be a single character vector of paths to fastq files, or several instances of the corresponding FastqFile, FastqSampler, or FastqStreamer objects.

pattern Ignored.

For developer use, to specify the underlying class contained in the FastqFileList.

Objects from the class

Available classes include:

FastqFile A file path and connection to a fastq file.

FastqFileList A list of FastqFile instances.

FastqSampler Uniformly sample records from a fastq file.

FastqStreamer Iterate over a fastq file, returning successive parts of the file.

Methods

The following methods are available to users:

```
readFastq, FastqFile-method: see also ?readFastq.
writeFastq, ShortReadQ, FastqFile-method: see also ?writeFastq, ?"writeFastq, ShortReadQ, FastqFile-method: yield: Draw a single sample from the instance. Operationally this requires that the underlying data (e.g., file) represented by the Sampler instance be visited; this may be time consuming.
```

See Also

readFastq, writeFastq, yield.

```
sp <- SolexaPath(system.file('extdata', package='ShortRead'))</pre>
fl <- file.path(analysisPath(sp), "s_1_sequence.txt")</pre>
f <- FastqFile(fl)
rfq <- readFastq(f)
f <- FastqSampler(f1, 50)</pre>
            # sample of size n=50
yield(f)
yield(f)
            # independent sample of size 50
close(f)
f <- FastqStreamer(f1, 50)</pre>
           # records 1 to 50
yield(f)
            # records 51 to 100
yield(f)
close(f)
## iterating over an entire file
f <- FastqStreamer(f1, 50)</pre>
while (length(fq <- yield(f))) {</pre>
    ## do work here
    print(length(fq))
```

GappedReads-class 21

```
}
close(f)

## Internal fields, methods, and help; for developers
ShortRead:::.FastqSampler_g$methods()
ShortRead:::.FastqSampler_g$fields()
ShortRead:::.FastqSampler_g$help("yield")
```

GappedReads-class

GappedReads objects

Description

The GappedReads class extends the GappedAlignments class defined in the GenomicRanges package. A GappedReads object contains all the information contained in a GappedAlignments object plus the sequences of the queries. Those sequences can be accessed via the qseq accessor.

Constructor

readGappedReads(file, format="BAM", use.names=FALSE, ...): Read a file as a GappedReads object. Like with the readGappedAlignments constructor, by default (i.e. use.names=FALSE) the resulting object has no names. If use.names is TRUE, then the names are constructed from the query template names (QNAME field in a SAM/BAM file).

Note that this function is just a front-end that delegates to the format-specific back-end function specified via the format argument. The use names argument and any extra argument are passed to the back-end function. Only the BAM format is supported for now. Its back-end is the readBamGappedReads function defined in the Rsamtools package. See ?readBamGappedReads for more information (you might need to install and load the Rsamtools package first).

Accessors

In the code snippets below, x is a GappedReads object.

qseq(x): Extracts the sequences of the queries as a DNAStringSet object.

Author(s)

H. Pages and P. Aboyoun

References

```
http://samtools.sourceforge.net/
```

See Also

GappedAlignments-class, readBamGappedReads

```
greads_file <- system.file("extdata", "ex1.bam", package="Rsamtools")
greads <- readGappedReads(greads_file)
greads
qseq(greads)</pre>
```

22 Intensity-class

short read image intensities	Intensity-class	"Intensity", "IntensityInfo", and "IntensityMeasure" base classes for short read image intensities
------------------------------	-----------------	--

Description

The Intensity, IntensityMeasure, and IntensityInfo classes represent and manipulate image intensity measures. Instances from the class may also contain information about measurement errors, and additional information about the reads from which the intensities are derived.

Intensity, and IntensityMeasure, are virtual classes, and cannot be created directly. Classes derived from IntensityMeasure (e.g., ArrayIntensity) and Intensity (e.g., SolexaIntensity) are used to represent specific technologies.

Objects from the Class

ArrayIntensity objects can be created with calls of the form ArrayIntensity(array(0, c(1,2,3))).

Objects of derived classes can be created from calls such as the SolexaIntensity constructor, or more typically by parsing appropriate files (e.g., readIntensities).

Slots

Class Intensity has slots:

readInfo: Object of class "IntensityInfo" containing columns for the lane, tile, x, and y coordinates of the read.

intensity: Object of class "IntensityMeasure" containing image intensity data for each read and cycle.

measurementError: Object of class "IntensityMeasure" containing measures of image intensity uncertainty for each read and cycle.

.hasMeasurementError: Length 1 logical variable indicating whether intensity standard errors are included (internal use only).

Classes IntensityInfo and IntensityMeasure are virtual classes, and have no slots.

Extends

These classes extend ".ShortReadBase", directly.

Methods

Methods and accessor functions for Intensity include:

```
readInfo signature(object = "Intensity"): access the readInfo slot of object.
```

intensity signature(object = "Intensity"): access the intensity slot of object.

measurementError signature(object = "Intensity"): access the nse slot of object, or signal an error if no standard errors are available.

dim signature(object = "Intensity"): return the dimensions (e.g., number of reads by number of cycles) represented by object.

show signature(object = "Intensity"): provide a compact representation of the object.

MAQMapQA-class 23

Subsetting "[" is available for the IntensityMeasure class; the drop argument to "[" is ignored.

Subsetting with "[[" is available for the ArrayIntensity class. The method accepts three arguments, corresponding to the read, base, and cycle(s) to be selected. The return value is the array (i.e., underlying data values) corresponding to the selected indices.

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

See Also

```
readIntensities
```

Examples

```
showMethods(class="Intensity", where=getNamespace("ShortRead"))
example(readIntensities)
```

MAQMapQA-class

Quality assessment summaries from MAQ map files

Description

This class contains a list-like structure with summary descriptions derived from visiting one or more MAQMap files.

Objects from the Class

Objects of the class are usually produced by a qa method.

Slots

.srlist: Object of class "list", containing data frames or lists of data frames summarizing the results of qa.

Extends

```
Class "SRList", directly. Class ".QA", directly. Class ".SRUtil", by class "SRList", distance 2. Class ".ShortReadBase", by class ".QA", distance 2.
```

Methods

Accessor methods are inherited from the SRList class.

```
report signature(x="MAQMapQA", ..., dest=tempfile(), type="html"): produces
    an html file summarizing the QA results.
```

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

24 qa

See Also

qa.

Examples

```
showClass("MAQMapQA")
```

qa

Perform quality assessment on short reads

Description

This function is a common interface to quality assessment functions available in ShortRead. Results from this function may be displayed in brief, or integrated into reports using, e.g., report.

Usage

Arguments

dirPath A character vector or other object (e.g., SolexaPath; see showMethods, below)

locating the data for which quality assessment is to be performed. See help pages for defined methods (by evaluating the example code, below) for details

of available methods.

pattern A character vector limiting the files in dirPath to be processed, as with list.files.

Care should be taken to specify pattern to avoid reading unintended files.

type The type of file being parsed; must be a character vector of length 1, selected

from one of the types enumerated in the parameter.

... Additional arguments used by methods.

fapply, reduce: Influence how evaluation occurs when this function is run with the **Rmpi** or **parallel** packages; see srapply.

Lpattern, Rpattern: A character vector or XString object to be matched to the left end of a sequence. If either Lpattern or Rpattern are provided, trimLRPatterns is invoked to produce a measure of adapter contamination. Mismatch rates are 0.1 on the left and 0.2 on the right, with a minimum overlap of 10 nt.

QA-class 25

Details

The following methods are defined, in addition to those on S4 formal classes documented elsewhere:

qa,character-method Quality assessment is performed on all files in directory dirPath whose file name matches pattern. The type of analysis performed is based on the type argument. Use SolexaExport when all files matching pattern are Solexa _export.txt files. Use SolexaRealign for Solexa _realign.txt files. Use Bowtie for Bowtie files. Use MAQMapShort for MAQ map files produced by MAQ versions below 0.70 and MAQMap for more recent output. Use fastq for collections of fastq-format files. Use BAM for .bam files; see the param argument in readAligned for information on how to influence reads included in the qa summary. Quality assessment details vary depending on data source.

qa,list-method dirPath is a list of objects, all of the same class and typically derived from ShortReadQ, on which quality assessment is performed. All elements of the list must have names, and these should be unique.

Value

An object derived from class .QA. Values contained in this object are meant for use by report

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

See Also

```
.QA, SolexaExportQA MAQMapQA FastqQA BAMQA
```

Examples

```
sp <- SolexaPath(system.file("extdata", package = "ShortRead"))
qa1 <- qa(sp)
fq <- file.path(analysisPath(sp), pattern="s_1_sequence.txt")
qa2 <- qa(readFastq(fq), basename(fq)) # name fastq files
showMethods("qa", where=getNamespace("ShortRead"))</pre>
```

QA-class

(Updated) classes for representing quality assessment results

Description

Classes derived from .QA-class represent results of quality assurance analyses.

Objects from the Class

Users create instances of many of these classes by calling the corresponding constructors, as documented on the help page for qa2. Classes constructed in this way include QACollate, QAFastqSource, QAAdapterContamination, QAFrequentSequence, QANucleotideByCycle, QANucleotideUse, QAQualityByCycle, QAQualityUse, QAReadQuality, and QASequenceUse.

The classes QASource, QAFiltered, QAFlagged and QASummary are generated internally, not by users.

26 QA-class

Extends

.QA2 extends class ".ShortReadBase", directly.

QASummary is a virtual class extending .QA2; all user-creatable classes extend QASummary.

QASource extends QASummary. All classes used to represent raw data input (QAFastqSource) extend QASource.

QAData is a reference class, used to contain a single instance of the fastq used in all QA Summary steps.

QACollate extends .QA2. It contains a SimpleList instance with zero or more QASummary elements.

QA extends .QA2, and contains a SimpleList of zero or more QASummary elements. This class represents the results of the qa2 analysis.

Methods

Methods defined on this class include:

qa2 signature(object="QACollate", state, ...,verbose=FALSE) creates a QA report from the elements of QACollate. Methods on qa2 for objects extending class QASummary summarize QA statistics for that class, e.g., qa2,QAFrequentSequences-method implements the calculations required to summarize frequently used sequences, using data in state.

report signature(x="QA", ...) creates an HTML report. Methods on report for objects extending class QASummary are responsible for creating the html snippet for that QA component.

flag signature(object=".QA2", ..., verbose=FALSE) implements criteria to flag individual lanes as failing quality assessment. NOTE: flag is not fully implemented.

rbind signature(...="QASummary"): rbind multiple summary elements of the same class, as when these have been created by separately calculating statistics on a number of fastq files.

show signature(object = "SolexaExportQA"): Display an overview of the object contents.

Author(s)

Martin Morgan <mtmmorgan@fhcrc.org>

See Also

Specific classes derived from .QA2

```
getClass(".QA2", where=getNamespace("ShortRead"))
```

qa2 27

qa2 (Updated) quality assessment reports on short reads

Description

This page summarizes an updated approach to quality assessment reports in ShortRead.

Usage

```
## Input source for short reads
QAFastqSource(con = character(), n = 1e+06, readerBlockSize = 1e+08,
    flagNSequencesRange = NA_integer_, ...,
    html = system.file("template", "QASources.html", package="ShortRead"))
QAData(seq = ShortReadQ(), filter = logical(length(seq)), ...)
## Possible QA elements
QAFrequentSequence(useFilter = TRUE, addFilter = TRUE,
    n = NA_integer_, a = NA_integer_, flagK=.8, reportSequences = FALSE,
QANucleotideByCycle(useFilter = TRUE, addFilter = TRUE, ...)
QANucleotideUse(useFilter = TRUE, addFilter = TRUE, ...)
QAQualityByCycle(useFilter = TRUE, addFilter = TRUE, ...)
QAQualityUse(useFilter = TRUE, addFilter = TRUE, ...)
QAReadQuality(useFilter = TRUE, addFilter = TRUE,
    flagK = 0.2, flagA = 30L, ...)
QASequenceUse(useFilter = TRUE, addFilter = TRUE, ...)
QAAdapterContamination(useFilter = TRUE, addFilter = TRUE,
   Lpattern = NA_character_, Rpattern = NA_character_,
   max.Lmismatch = 0.1, max.Rmismatch = 0.2, min.trim = 9L, ...)
## Order QA report elements
QACollate(src, ...)
## perform analysis
qa2(object, state, ..., verbose=FALSE)
## Outputs from qa2
QA(src, filtered, flagged, ...)
QAFiltered(useFilter = TRUE, addFilter = TRUE, ...)
QAFlagged(useFilter = TRUE, addFilter = TRUE, ...)
## Summarize results as html report
## S4 method for signature 'QA'
report(x, ..., dest = tempfile(), type = "html")
## additional methods; 'flag' is not fully implemented
flag(object, ..., verbose=FALSE)
## S4 method for signature 'QASummary'
rbind(..., deparse.level = 1)
```

28 qa2

Arguments

con character(1) file location of fastq input, as used by FastqSampler.

n integer(1) number of records to input, as used by FastqStreamer (QAFastqSource).

integer(1) number of sequences to tag as 'frequent' (QAFrequentSequence).

readerBlockSize

integer(1) number of bytes to input, as used by FastqStreamer.

flagNSequencesRange

integer(2) minimum and maximum reads above which source files will be

flagged as outliers.

html character(1) location of the HTML template for summarizing this report ele-

ment.

seq ShortReadQ representation of fastq data.

filter logical() vector with length equal to seq, indicating whether elements of seq

are filtered (TRUE) or not.

useFilter, addFilter

logical(1) indicating whether the QA element should be calculating using the filtered (useFilter=TRUE) or all reads, and whether reads failing the QA element should be added to the filter used by subsequent steps (addFilter = TRUE)

or not.

a integer (1) count of number of sequences above which a read will be consid-

ered 'frequent' (QAFrequentSequence).

flagK, flagA flagK numeric(1) between 0 and 1 indicating the fraction of frequent sequences

greater than or equal to n or a above which a fastq file will be flagged (QAFrequentSequence).

flagK numeric{1} between 0 and 1 and flagA integer(1) indicating that a run should be flagged when the fraction of reads with quality greater than or equal

to flagA falls below threshold flagK.

reportSequences

logical(1) indicating whether frequent sequences are to be reported.

Lpattern, Rpattern, max.Lmismatch, max.Rmismatch, min.trim

Parameters influencing adapter identification, see matchPattern.

src The source, e.g., QAFastqSource, on which the quality assessment report will

be based.

object An instance of class derived from QA on which quality metrics will be derived;

for end users, this is usually the result of QACollate..

state The data on which quality assessment will be performed; this is not usually

necessary for end-users.

verbose logical(1) indicating whether progress reports should be reported.

filtered, flagged

Primarily for internal use, instances of QAFiltered and QAFlagged.

x An instance of QA on which a report is to be generated.

dest character(1) providing the directory in which the report is to be generated.

type character(1) indicating the type of report to be generated; only "html" is sup-

ported.

deparse.level see rbind.

.. Additional arguments, e.g., html to specify the location of the html source to

use as a template for the report.

qa2 29

Details

Use QACollate to specify an order in which components of a QA report are to be assembled. The first argument is the data source (e.g., QAFastqSource).

Functions related to data input include:

QAFastqSource defines the location of fastq files to be included in the report. con is used to construct a FastqSampler instance, and records are processed using ga2, QAFastqSource-method.

QAData is a class for representing the data during the QA report generation pass; it is primarily for internal use.

Possible elements in a QA report are:

QAFrequentSequence identifies the most-commonly occuring sequences. One of n or a can be non-NA, and determine the number of frequent sequences reported. n specifies the number of most-frequent sequences to filter, e.g., n=10 would filter the top 10 most commonly occurring sequences; a provides a threshold frequency (count) above which reads are filtered. The sample is flagged when a fraction flagK of the reads are filtered.

reportSequences determines whether the most commonly occurring sequences, as determined by n or a, are printed in the html report.

QANucleotideByCycle reports nucleotide frequency as a function of cycle.

QAQualityByCycle reports average quality score as a function of cycle.

QAQualityUse summarizes overall nucleotide qualities.

QAReadQuality summarizes the distribution of read qualities.

QASequenceUse summarizes the cumulative distribution of reads occurring 1, 2, ... times.

QAAdapterContamination reports the occurrence of 'adapter' sequences on the left and / or right end of each read.

Value

An object derived from class .QA. Values contained in this object are meant for use by report

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

See Also

QA.

```
sp <- SolexaPath(system.file('extdata', package='ShortRead'))
fl <- file.path(analysisPath(sp), "s_1_sequence.txt")
fls <- c(fl, fl)

coll <- QACollate(QAFastqSource(fls), QAReadQuality(),
        QAAdapterContamination(), QANucleotideUse(),
        QAQualityUse(), QASequenceUse(),
        QAFrequentSequence(n=10), QANucleotideByCycle(),
        QAQualityByCycle())
x <- qa2(coll, verbose=TRUE)
res <- report(x)</pre>
```

30 QualityScore

```
if (interactive())
    browseURL(res)

## Not run: ## parallel evaluation
options(srapply_fapply="parallel")
coll <- QACollate(QAFastqSource(fls), qaFastqTemplate)
x2 <- qa2(coll, verbose=TRUE)
browseURL(res2 <- report(x2))

## End(Not run)</pre>
```

QualityScore

Construct objects indicating read or alignment quality

Description

Use these functions to construct quality indicators for reads or alignments. See QualityScore for details of object content and methods available for manipulating them.

Usage

```
NumericQuality(quality = numeric(0))
IntegerQuality(quality = integer(0))
MatrixQuality(quality = new("matrix"))
FastqQuality(quality, ...)
SFastqQuality(quality, ...)
```

Arguments

quality

An object used to initialize the data structure. Appropriate objects are indicated in the constructors above for Numeric, Integer, and Matrix qualities. For FastqQuality and SFastqQuality, methods are defined for BStringSet, character, and missing.

... Additional arguments, currently unused.

Value

Constructors return objects of the corresponding class derived from QualityScore.

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

See Also

```
QualityScore, readFastq, readAligned
```

```
nq <- NumericQuality(rnorm(20))
nq
quality(nq)
quality(nq[10:1])</pre>
```

QualityScore-class 31

QualityScore-class

Quality scores for short reads and their alignments

Description

This class hierarchy represents quality scores for short reads. QualityScore is a virtual base class, with derived classes offering different ways of representing qualities. Methods defined on QualityScore are implemented in all derived classes.

Objects from the Class

Objects from the class are created using constructors (e.g., NumericQuality) named after the class name.

Defined classes are as follows:

QualityScore Virtual base class; instances cannot be instantiated.

NumericQuality A single numeric vector, where values represent quality scores on an arbitrary scale.

IntegerQuality A integer numeric vector, where values represent quality scores on an arbitrary scale.

MatrixQuality A rectangular matrix of quality scores, with rows representing reads and columns cycles. The content and interpretation of row and column entries is arbitrary; the rectangular nature implies quality scores from equal-length reads.

FastqQuality 'fastq' encoded quality scores stored in a BStringSet instance. Base qualities of a single read are represented as an ASCII character string. The integer-valued quality score of a single base is encoded as its ASCII equivalent plus 33. The precise definition of the integer-valued quality score is unspecified, but is usually a Phred score; the meaning can be determined from the source of the quality scores. Multiple reads are stored as a BStringSet, and so can be of varying lengths.

SolexaQuality As with FastqQuality, but with integer qualities encoded as ASCII equivalent plus 64.

Extends

```
Class ". ShortReadBase", directly.
```

Methods

The following methods are defined on all QualityScore and derived classes:

```
[ signature(x = "QualityScore", i = "ANY", j = "missing")
[ signature(x = "MatrixQuality", i = "ANY", j = "missing"):
```

Subset the object, with index i indicating the reads for which quality scores are to be extracted. The class of the result is the same as the class of x. It is an error to provide any argument other than i.

```
[[ signature(x = "QualityScore", i = "ANY", j = "ANY"):
```

Subset the object, returning the quality score (e.g., numeric value) of the ith read.

```
[[ signature(x = "MatrixQuality", i = "ANY", j = "ANY"):
```

Returns the vector of quality scores associated with the ith read.

32 QualityScore-class

```
dim signature(x = "MatrixQuality"):
     The integer(2) dimension (e.g., number of reads, read width) represented by the quality score.
length signature(x = "QualityScore"):
length signature(x = "MatrixQuality"):
     The integer(1) length (e.g., number of reads) represented by the quality score. Note that
     length of MatrixQuailty is the number of rows of the corresponding matrix, and not the
     length of the corresponding numeric vector.
append signature(x = "QualityScore", values = "QualityScore", length = "missing"):
     append values after x.
width signature(x = "QualityScore"):
width signature(x = "NumericQuality"):
width signature(x = "MatrixQuality"):
width signature(x = "FastqQuality"):
     A numeric vector with length equal to the number of quality scores, and value equal to the
     number of quality scores for each read. For instance, a FastqQuality will have widths equal
     to the number of nucleotides in the underlying short read.
show signature(object = "QualityScore"):
show signature(object = "NumericQuality"):
show signature(object = "FastqQuality"):
     provide a brief summary of the object content.
detail signature(x = "QualityScore"):
     provide a more detailed view of object content.
The following methods are defined on specific classes:
alphabet signature(x = "FastqQuality", ...): Return a character vector of valid quality
     characters.
alphabetFrequency signature(stringSet = "FastqQuality"):
     Apply alphabetFrequency to quality scores, returning a matrix as described in alphabetFrequency.
alphabetByCycle signature(stringSet = "FastqQuality"):
     Apply alphabetByCycle to quality scores, returning a matrix as described in alphabetByCycle.
alphabetScore signature(object = "FastqQuality"):
alphabetScore signature(object = "SFastqQuality"):
     Apply alphabetScore (i.e., summed base quality, per read) to object.
coerce signature(from = "FastqQuality", to = "numeric"):
coerce signature(from = "FastqQuality", to = "matrix"):
coerce signature(from = "FastqQuality", to = "PhredQuality"):
coerce signature(from = "SFastqQuality", to = "matrix"):
coerce signature(from = "SFastqQuality", to = "SolexaQuality"):
     Use like as(from, "matrix")) to coerce objects of class from to class to, using the quality
     encoding implied by the class. When to is "matrix", all quality scores must be of the same
     width, and the result is a matrix of type integer. The result always represents the integer
     encoding of the corresponding quality string.
narrow signature(x = "FastqQuality", start = NA, end = NA, width = NA, use.names = TRUE):
     'narrow' quality so that scores are between start and end bases, according to narrow in the
     IRanges package.
```

```
trimTailw signature(object="FastqQuality", k="integer",a="character", halfwidth="integer", ..., trim trailing nucleotides when a window of width 2 * halfwidth + 1 contains k or more quality scores falling at or below a.
```

trimTails signature(object="FastqQuality", k="integer",a="character", successive=FALSE, ..., rang trim trailing scores if k scores fall below the quality encoded by a. If successive=FALSE, the k'th failing score and all subsequent scores are trimmed. If successive=TRUE, failing scores must occur successively; the sequence is trimmed from the first of the successive failing score.

```
srorder signature(x = "FastqQuality"):
srrank signature(x = "FastqQuality"):
srduplicated signature(x = "FastqQuality"):
```

Apply srsort, srorder, srrank, and srduplicated to quality scores, returning objects as described on the appropriate help page.

Integer representations of SFastqQuality and FastqQuality can be obtained with as(x, "matrix").

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

See Also

NumericQuality and other constructors.

Examples

```
names(slot(getClass("QualityScore"), "subclasses"))
```

readAligned

Read aligned reads and their quality scores into R representations

Description

Import files containing aligned reads into an internal representation of the alignments, sequences, and quality scores. Most methods (see 'details' for exceptions) read all files into a single R object.

Usage

```
readAligned(dirPath, pattern=character(0), ...)
```

Arguments

dirPath	A character vector (or other object; see methods defined on this generic) giving the directory path (relative or absolute; some methods also accept a character vector of file names) of aligned read files to be input.
pattern	The (grep-style) pattern describing file names to be read. The default (character(0)) results in (attempted) input of all files in the directory.
	Additional arguments, used by methods. When dirPath is a character vector, the argument type must be provided. Possible values for type and their meaning are described below. Most methods implement filter=srFilter(), allowing objects of SRFilter to selectively returns aligned reads.

Details

There is no standard aligned read file format; methods parse particular file types.

The readAligned, character-method interprets file types based on an additional type argument. Supported types are:

type="SolexaExport" This type parses .*_export.txt files following the documentation in the Solexa Genome Alignment software manual, version 0.3.0. These files consist of the following columns; consult Solexa documentation for precise descriptions. If parsed, values can be retrieved from AlignedRead as follows:

Machine see below

Run number stored in alignData

Lane stored in alignData

Tile stored in alignData

 ${f X}$ stored in alignData

Y stored in alignData

Multiplex index see below

Paired read number see below

Read sread

Quality quality

Match chromosome chromosome

Match contig alignData

Match position position

Match strand strand

Match description Ignored

Single-read alignment score alignQuality

Paired-read alignment score Ignored

Partner chromosome Ignored

Partner contig Ignored

Partner offset Ignored

Partner strand Ignored

Filtering alignData

The following optional arguments, set to FALSE by default, influence data input

withMultiplexIndex When TRUE, include the multiplex index as a column multiplexIndex in alignData.

withPairedReadNumber When TRUE, include the paired read number as a column pairedReadNumber in alignData.

withId When TRUE, construct an identifier string as 'Machine_Run:Lane:Tile:X:Y#multiplexIndex/pairedReadNum The substrings '#multiplexIndex' and '/pairedReadNumber' are not present if withMultiplexIndex=FALSE or withPairedReadNumber=FALSE.

withAll A convencience which, when TRUE, sets all with* values to TRUE.

Note that not all paired read columns are interpreted. Different interfaces to reading alignment files are described in SolexaPath and SolexaSet.

type="SolexaPrealign" See SolexaRealign
type="SolexaAlign" See SolexaRealign

type="SolexaRealign" These types parse s_L_TTTT_prealign.txt, s_L_TTTT_align.txt or s_L_TTTT_realign.txt files produced by default and eland analyses. From the Solexa documentation, align corresponds to unfiltered first-pass alignments, prealign adjusts alignments for error rates (when available), realign filters alignments to exclude clusters failing to pass quality criteria.

Because base quality scores are not stored with alignments, the object returned by readAligned scores all base qualities as -32.

If parsed, values can be retrieved from AlignedRead as follows:

Sequence stored in sread

Best score stored in alignQuality
Number of hits stored in alignData
Target position stored in position

Strand stored in strand

Target sequence Ignored; parse using readXStringColumns

Next best score stored in alignData

type="SolexaResult" This parses s_L_eland_results.txt files, an intermediate format that does not contain read or alignment quality scores.

Because base quality scores are not stored with alignments, the object returned by readAligned scores all base qualities as -32.

Columns of this file type can be retrieved from AlignedRead as follows (description of columns is from Table 19, Genome Analyzer Pipeline Software User Guide, Revision A, January 2008):

Id Not parsed

Sequence stored in sread

Type of match code Stored in alignData as matchCode. Codes are (from the Eland manual): NM (no match); QC (no match due to quality control failure); RM (no match due to repeat masking); U0 (best match was unique and exact); U1 (best match was unique, with 1 mismatch); U2 (best match was unique, with 2 mismatches); R0 (multiple exact matches found); R1 (multiple 1 mismatch matches found, no exact matches); R2 (multiple 2 mismatch matches found, no exact or 1-mismatch matches).

Number of exact matches stored in alignData as nExactMatch

Number of 1-error mismatches stored in alignData as nOneMismatch

Number of 2-error mismatches stored in alignData as nTwoMismatch

Genome file of match stored in chromosome

Position stored in position

Strand (direction of match) stored in strand

'N' treatment stored in alignData, as NCharacterTreatment. '.' indicates treatment of 'N' was not applicable; 'D' indicates treatment as deletion; 'l' indicates treatment as insertion

Substitution error stored in alignData as mismatchDetailOne and mismatchDetailTwo. Present only for unique inexact matches at one or two positions. Position and type of first substitution error, e.g., 11A represents 11 matches with 12th base an A in reference but not read. The reference manual cited below lists only one field (mismatchDetailOne), but two are present in files seen in the wild.

- type="MAQMap", records=-1L Parse binary map files produced by MAQ. See details in the next section. The records option determines how many lines are read; -1L (the default) means that all records are input. For type="MAQMap", dir and pattern must match a single file.
- type="MAQMapShort", records=-1L The same as type="MAQMap" but for map files made with Maq prior to version 0.7.0. (These files use a different maximum read length [64 instead of 128], and are hence incompatible with newer Maq map files.). For type="MAQMapShort", dir and pattern must match a single file.

type="MAQMapview" Parse alignment files created by MAQ's 'mapiew' command. Interpretation of columns is based on the description in the MAQ manual, specifically

...each line consists of read name, chromosome, position, strand, insert size from the outer coordinates of a pair, paired flag, mapping quality, single-end mapping quality, alternative mapping quality, number of mismatches of the best hit, sum of qualities of mismatched bases of the best hit, number of 0-mismatch hits of the first 24bp, number of 1-mismatch hits of the first 24bp on the reference, length of the read, read sequence and its quality.

The read name, read sequence, and quality are read as XStringSet objects. Chromosome and strand are read as factors. Position is numeric, while mapping quality is numeric. These fields are mapped to their corresponding representation in AlignedRead objects.

Number of mismatches of the best hit, sum of qualities of mismatched bases of the best hit, number of 0-mismatch hits of the first 24bp, number of 1-mismatch hits of the first 24bp are represented in the AlignedRead object as components of alignData.

Remaining fields are currently ignored.

type="Bowtie" Parse alignment files created with the Bowtie alignment algorithm. Parsed columns can be retrieved from AlignedRead as follows:

Identifier id
Strand strand

Chromosome chromosome

Position position; see comment below

Read sread; see comment below

Read quality quality; see comments below

Similar alignments alignData, 'similar' column; Bowtie v. 0.9.9.3 (12 May, 2009) documents this as the number of other instances where the same read aligns against the same reference characters as were aligned against in this alignment. Previous versions marked this as 'Reserved'

Alignment mismatch locations alignData 'mismatch', column

NOTE: the default quality encoding changes to FastqQuality with **ShortRead** version 1.3.24. This method includes the argument qualityType to specify how quality scores are encoded. Bowtie quality scores are 'Phred'-like by default, with qualityType='FastqQuality', but

can be specified as 'Solexa'-like, with qualityType='SFastqQuality'.

Bowtie outputs positions that are 0-offset from the left-most end of the + strand. ShortRead parses position information to be 1-offset from the left-most end of the + strand.

Bowtie outputs reads aligned to the - strand as their reverse complement, and reverses the quality score string of these reads. ShortRead parses these to their original sequence and orientation.

type="SOAP" Parse alignment files created with the SOAP alignment algorithm. Parsed columns can be retrieved from AlignedRead as follows:

id id
seq sread; see comment below
qual quality; see comment below
number of hits alignData
a/b alignData (pairedEnd)

readAligned 37

```
length alignData(alignedLength)
+/- strand
chr chromosome
location position; see comment below
types alignData(typeOfHit: integer portion; hitDetail: text portion)
```

This method includes the argument qualityType to specify how quality scores are encoded. It is unclear from SOAP documentation what the quality score is; the default is 'Solexa'-like, with qualityType='SFastqQuality', but can be specified as 'Phred'-like, with qualityType='FastqQuality'.

SOAP outputs positions that are 1-offset from the left-most end of the + strand. ShortRead preserves this representation.

SOAP reads aligned to the - strand are reported by SOAP as their reverse complement, with the quality string of these reads reversed. ShortRead parses these to their original sequence and orientation.

type="BAM" Parse BAM files produced by samtools and other third party programs. This method includes the argument param=ScanBamParam(). The ScanBamParam object is recycled for all files. The which and flag arguments to ScanBamParam() can be used to influence which reads in the BAM file are parsed; see ScanBamParam. The following values override user settings (issuing a warning if contradictory values are provided):

simpleCigar=TRUE Reads aligned with indels are ignored; this is required for representation in AlignedRead.

reverseComplement=TRUE By default, BAM stores reads as they are aligned to the reference genome, whereas AiignedRead stores them as they are prior to alignment; this flag converts reads from the BAM to AlignedRead format.

what=c("qname", "flag", "rname", "strand", "pos", "mapq", "seq", "qual") These BAM fields are mapped to corresponding fields in AlignedRead.

BAM fields are mapped to AlignedRead as:

```
qname id
seq sread
qual quality
strand strand
rname chromosome
pos position
mapq alignQuality
flag alignData
```

Value

A single R object (e.g., AlignedRead) containing alignments, sequences and qualities of all files in dirPath matching pattern. There is no guarantee of order in which files are read.

Author(s)

Martin Morgan mtmorgan@fhcrc.org, Simon Anders anders@ebi.ac.uk (MAQ map)

See Also

The AlignedRead class.

Genome Analyzer Pipeline Software User Guide, Revision A, January 2008.

38 readBaseQuality

```
The MAQ reference manual, http://maq.sourceforge.net/maq-manpage.shtml#5, 3 May, 2008. The Bowtie reference manual, http://bowtie-bio.sourceforge.net, 28 October, 2008. The SOAP reference manual, http://soap.genomics.org.cn/soap1, 16 December, 2008. The BAM file format specification, http://samtools.sourceforge.net.
```

Examples

```
sp <- SolexaPath(system.file("extdata", package="ShortRead"))</pre>
ap <- analysisPath(sp)</pre>
## ELAND_EXTENDED
(aln0 <- readAligned(ap, "s_2_export.txt", "SolexaExport"))</pre>
## PhageAlign
(aln1 <- readAligned(ap, "s_5_.*_realign.txt", "SolexaRealign"))</pre>
## MAO
dirPath <- system.file('extdata', 'maq', package='ShortRead')</pre>
list.files(dirPath)
## First line
readLines(list.files(dirPath, full.names=TRUE)[[1]], 1)
countLines(dirPath)
## two files collapse into one
(aln2 <- readAligned(dirPath, type="MAQMapview"))</pre>
## select only chr1-5.fa, '+' strand
filt <- compose(chromosomeFilter("chr[1-5].fa"),</pre>
                 strandFilter("+"))
(aln3 <- readAligned(sp, "s_2_export.txt", filter=filt))</pre>
```

readBaseQuality

Read short reads and their quality scores into R representations

Description

readBaseQuality reads all base call files in a directory dirPath whose file name matches seqPattern and all quality score files whose name matches prbPattern, returning a compact internal representation of the sequences, and quality scores in the files. Methods read all files into a single R object.

Usage

```
readBaseQuality(dirPath, ...)
## S4 method for signature 'character'
readBaseQuality(dirPath, seqPattern=character(0),
prbPattern=character(0), type=c("Solexa"), ...)
```

Arguments

dirPath A character vector (or other object; see methods defined on this generic) giving

the directory path (relative or absolute) of files to be input.

seqPattern The (grep-style) pattern describing base call file names to be read. The default

(character(0)) results in (attempted) input of all files in the directory.

readBfaToc 39

prbPattern	The (grep-style) pattern describing quality score file names to be read. The
	default (character(0)) results in (attempted) input of all files in the directory.
type	The type of file to be parsed. Supported types include: Solexa: parse reads and their qualities from _seq.txt and _prb.txt-formatted files, respectively.
	Additional arguments, perhaps used by methods.

Value

A single R object (e.g., ShortReadQ) containing sequences and qualities of all files in dirPath matching seqPattern and prbPattern respectively. There is no guarantee of order in which files are read.

Author(s)

See Also

```
A ShortReadQ object.
readXStringColumns, readPrb
```

Examples

```
sp <- SolexaPath(system.file("extdata", package="ShortRead"))
readBaseQuality(sp, seqPattern="s_1.*_seq.txt", prbPattern="s_1.*_prb.txt")</pre>
```

readBfaToc

Get a list of the sequences in a Maq .bfa file

Description

As coverage needs to know the lengths of the reference sequences, this function is provided which extracts this information from a .bfa file (Maq's "binary FASTA" format).

Usage

```
readBfaToc( bfafile )
```

Arguments

bfafile The file name of the .bfa file.

Value

An integer vector with one element per reference sequence found in the .bfa file, each vector element named with the sequence name and having the sequence length as value.

Author(s)

```
Simon Anders, EMBL-EBI, <sanders@fs.tum.de>
```

(Note: The C code for this function incorporates code from Li Heng's MAQ software, (c) Li Heng and released by him under GPL 2.

40 readFasta

readFasta Read and write FASTA files to or from ShortRead objects

Description

readFasta reads all FASTA-formated files in a directory dirPath whose file name matches pattern pattern, returning a compact internal representation of the sequences and quality scores in the files. Methods read all files into a single R object; a typical use is to restrict input to a single FASTA file.

writeFasta writes an object to a single file, using mode="w" (the default) to create a new file or mode="a" append to an existing file. Attempting to write to an existing file with mode="w" results in an error.

Usage

```
readFasta(dirPath, pattern = character(0), ...,
    nrec=-1L, skip=0L)
## S4 method for signature 'character'
readFasta(dirPath, pattern = character(0), ...,
    nrec=-1L, skip=0L)
writeFasta(object, file, mode="w", ...)
## S4 method for signature 'DNAStringSet'
writeFasta(object, file, mode="w", ...)
```

Arguments

dirPath	A character vector giving the directory path (relative or absolute) or single file name of FASTA files to be read.
pattern	The (grep-style) pattern describing file names to be read. The default (character(0)) results in (attempted) input of all files in the directory.
object	An object to be output in fasta format.
file	A length 1 character vector providing a path to a file to the object is to be written to.
mode	A length 1 character vector equal to either 'w' or 'a' to write to a new file or append to an existing file, respectively.
	$Additional\ arguments\ used\ by\ methods\ or,\ for\ write \texttt{Fasta},\ write\ .\ \texttt{XStringSet}.$
nrec	See ?read.DNAStringSet.
skip	See ?read.DNAStringSet.

Value

readFasta returns a DNAStringSet. containing sequences and qualities contained in all files in dirPath matching pattern. There is no guarantee of order in which files are read.

writeFasta is invoked primarily for its side effect, creating or appending to file file. The function returns, invisibly, the length of object, and hence the number of records written. There is a writeFasta method for any class derived from ShortRead.

readFastq 41

Author(s)

Martin Morgan

Examples

```
showMethods("readFasta")
showMethods("writeFasta")
f1 <- system.file("extdata", "someORF.fa", package="Biostrings")
rfa <- readFasta(f1)
sread(rfa)
id(rfa)

sp <- SolexaPath(system.file('extdata', package='ShortRead'))
rfq <- readFastq(analysisPath(sp), pattern="s_1_sequence.txt")
file <- tempfile()
writeFasta(rfq, file)
readLines(file, 8)
writeFasta(sread(rfq), file) # no 'id's</pre>
```

readFastq

Read and write FASTQ-formatted files

Description

readFastq reads all FASTQ-formated files in a directory dirPath whose file name matches pattern pattern, returning a compact internal representation of the sequences and quality scores in the files. Methods read all files into a single R object; a typical use is to restrict input to a single FASTQ file.

writeFastq writes an object to a single file, using mode="w" (the default) to create a new file or mode="a" append to an existing file. Attempting to write to an existing file with mode="w" results in an error.

Usage

```
readFastq(dirPath, pattern=character(0), ...)
## S4 method for signature 'character'
readFastq(dirPath, pattern=character(0), ..., withIds=TRUE)
writeFastq(object, file, mode="w", full=FALSE, ...)
```

Arguments

dirPath A character vector (or other object; see methods defined on this generic) giving

the directory path (relative or absolute) or single file name of FASTQ files to be

read.

pattern The (grep-style) pattern describing file names to be read. The default (character(0))

results in (attempted) input of all files in the directory.

42 readFastq

object	An object to be output in fastq format. For methods, use showMethods(object, where=getNamesp
file	A length 1 character vector providing a path to a file to the object is to be written to.
mode	A length 1 character vector equal to either 'w' or 'a' to write to a new file or append to an existing file, respectively.
full	A logical(1) indicating whether the identifier line should be repeated full=TRUE or omitted full=FALSE on the third line of the fastq record.
	Additional arguments. In particular, qualityType and filter:
	qualityType: Representation to be used for quality scores, must be one of Auto (choose Phred-like if any character is ASCII-encoded as less than 59) FastqQuality (Phred-like encoding), SFastqQuality (Illumina encoding).
	filter: An object of class srFilter, used to filter objects of class ShortReadQ at input.
withIds	logical(1) indicating whether identifiers should be read from the fastq file.

Details

The fastq format is not quite precisely defined. The basic definition used here parses the following four lines as a single record:

```
@HWI-EAS88_1_1_1_1001_499
GGACTTTGTAGGATACCCTCGCTTTCCTTCTCTGT
+HWI-EAS88_1_1_1_1001_499
]]]]]]]]]Y]Y]]]]]]]VCHVMPLAS
```

The first and third lines are identifiers preceded by a specific character (the identifiers are identical, in the case of Solexa). The second line is an upper-case sequence of nucleotides. The parser recognizes IUPAC-standard alphabet (hence ambiguous nucleotides), coercing . to - to represent missing values. The final line is an ASCII-encoded representation of quality scores, with one ASCII character per nucleotide.

The encoding implicit in Solexa-derived fastq files is that each character code corresponds to a score equal to the ASCII character value minus 64 (e.g., ASCII @ is decimal 64, and corresponds to a Solexa quality score of 0). This is different from BioPerl, for instance, which recovers quality scores by subtracting 33 from the ASCII character value (so that, for instance, !, with decimal value 33, encodes value 0).

The BioPerl description of fastq asserts that the first character of line 4 is a !, but the current parser does not support this convention.

writeFastq creates files following the specification outlined above, using the IUPAC-standard alphabet (hence, sequences containing '.' when read will be represented by '-' when written).

Value

readFastq returns a single R object (e.g., ShortReadQ) containing sequences and qualities contained in all files in dirPath matching pattern. There is no guarantee of order in which files are read.

writeFastq is invoked primarily for its side effect, creating or appending to file file. The function returns, invisibly, the length of object, and hence the number of records written.

readIntensities 43

Author(s)

Martin Morgan

See Also

```
The IUPAC alphabet in Biostrings.
```

```
http://www.bioperl.org/wiki/FASTQ_sequence_format for the BioPerl definition of fastq.
```

Solexa documentation 'Data analysis - documentation : Pipeline output and visualisation'.

Examples

```
showMethods(readFastq)
showMethods(writeFastq)

sp <- SolexaPath(system.file('extdata', package='ShortRead'))
rfq <- readFastq(analysisPath(sp), pattern="s_1_sequence.txt")
sread(rfq)
id(rfq)
quality(rfq)

## SolexaPath method 'knows' where FASTQ files are placed
rfq1 <- readFastq(sp, pattern="s_1_sequence.txt")
rfq1

file <- tempfile()
writeFastq(rfq, file)
readLines(file, 8)</pre>
```

 ${\tt readIntensities}$

Read Illumina image intensity files

Description

 $readIntensities\ reads\ image\ 'intensity'\ files\ (such\ as\ Illumina's\ _int.\ txt\ and\ (optionally)\ _nse.\ txt)$ into a single object.

Usage

```
readIntensities(dirPath, pattern=character(0), ...)
```

Arguments

dirPath	Directory path or other object (e.g., SolexaPath) for which methods are defined.
pattern	A length 1 character vector representing a regular expression to be combined with dirPath, as described below, to match files to be summarized.
	Additional arguments used by methods.

44 readIntensities

Details

Additional methods are defined on specific classes, see, e.g., SolexaPath.

The readIntensities, character-method contains an argument type that determines how intensities are parsed. Use the type argument to readIntensities, character-method, as described below. All readIntensities, character methods accepts the folling arguments:

with Variability: Include estimates of variability (i.e., from parsing _nse files).

verbose: Report on progress when starting to read each file.

The supported types and their signatures are:

```
type="RtaIntensity" Intensities are read from Illumina _cif.txt and _cnf.txt-style files. The signature for this method is
```

```
dirPath, pattern=character(0), ..., type="RtaIntensity", lane=integer(0), cycles=integer(0), cycles=int
```

lane: integer(1) identifying the lane in which cycles and tiles are to be processed.

cycles: integer() enumerating cycles to be processed.

cycleIteration: integer(1) identifying the iteration of the base caller to be summarized

tiles: integer() enumerating tile numbers to be summarized.

laneName, cycleNames, tileNames, posNames: character() vectors identifying the lane and cycle directories, and the 'pos' and tile file names (excluding the '.cif' or '.cnf' extension) to be processed.

The dirPath and pattern arguments are combined as list.files(dirPath, pattern), and must identify a single directory. Most uses of this function will focus on a single tile (specified with, e.g., tiles=1L); the laneName, cycleNames, tileNames, and posNames parameters are designed to work with the default Illumina pipeline and do not normally need to be specified.

```
type="IparIntensity" Intensities are read from Solexa _pos.txt, _int.txt.p, _nse.txt.p-style file triplets. The signature for this method is
```

```
dirPath, pattern=character(0), ...,type="IparIntensity", intExtension="_int.txt.p.gz
Files to be parsed are determined as, e.g., paste(pattern, intExtension, sep="").
```

type="SolexaIntensity" Intensities are read from Solexa _int.txt and _nse.txt-style files.

```
The signature for this method is dirPath, pattern=character(0), ...,type="SolexaIntensity", intExtension="_int.txt",
```

```
Files to be parsed are determined as, e.g., paste(pattern, intExtension, sep="").
```

Value

An object derived from class Intensity.

Author(s)

Martin Morgan mtmorgan@fhcrc.org, Michael Muratet mtmorgan@fhcrc.org, Michael Muratet mtmorgan@hcrc.org, Michael Muratet mtmorgan@hcrc.org), Michael Mura

readPrb 45

Examples

readPrb

Read Solexa prb files as fastq-style quality scores

Description

readPrb reads all _prb.txt files in a directory into a single object. Most methods (see details) do this by identifying the maximum base call quality for each cycle and read, and representing this as an ASCII-encoded character string.

Usage

```
readPrb(dirPath, pattern = character(0), ...)
```

Arguments

dirPath Directory path or other object (e.g., SolexaPath for which methods are defined.

Regular expression matching names of _prb files to be summarized.

Additional arguments, e.g., to srapply, used during evaluation.

Details

The readPrb, character-method contains an argument as that determines the value of the returned object, as follows.

- as="SolexaEncoding" The ASCII encoding of the maximum per cycle and read quality score is encoded using Solexa conventions.
- as="FastqEncoding" The ASCII encoding of the maximum per cycle and read quality score is encoded using Fastq conventions, i.e., ! has value 0.
- as="IntegerEncoding" The maximum per cycle and read quality score is returned as a in integer value. Values are collated into a matrix with number of rows equal to number of reads, and number of columns equal to number of cycles.
- as="array" The quality scores are *not* summarized; the return value is an integer array with dimensions corresponding to reads, nucleotides, and cycles.

46 readQseq

Value

An object of class QualityScore, or an integer matrix.

Author(s)

Martin Morgan mtmorgan@fhcrc.org

Examples

```
fl <- system.file("extdata", package="ShortRead")
sp <- SolexaPath(fl)
readPrb(sp, "s_1.*_prb.txt") # all tiles to a single file</pre>
```

readQseq

Read Solexa qseq files as fastq-style quality scores

Description

readQseq reads all files matching pattern in a directory into a single ShortReadQ-class object. Information on machine, lane, tile, x, and y coordinates, filtering status, and read number are not returned (although filtering status can be used to selectively include reads as described below).

Usage

Arguments

dirPath	Directory path or other object (e.g., SolexaPath) for which methods are defined.
pattern	Regular expression matching names of _qseq files to be summarized.
	Additional argument, passed to I/O functions.
as	character(1) indicating the class of the return type. "XDataFrame" is included
	for backward compatibility, but is no longer supported.
filtered	logical(1) indicating whether to include only those reads passing Solexa filtering?
verbose	logical(1) indicating whether to report on progress during evaluation.

Value

An object of class ShortReadQ.

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

```
fl <- system.file("extdata", package="ShortRead")
sp <- SolexaPath(fl)
readQseq(sp)</pre>
```

readXStringColumns 47

<u> </u>	Read one or more columns into XStringSet (e.g., DNAStringSet) ob- ects
----------	---

Description

This function allows short read data components such as DNA sequence, quality scores, and read names to be read in to XStringSet (e.g., DNAStringSet, BStringSet) objects. One or several files of identical layout can be specified.

Usage

```
readXStringColumns(dirPath, pattern=character(0),
                   colClasses=list(NULL),
                   nrows=-1L, skip=0L,
                   sep = "\t", header = FALSE, comment.char="#")
```

Arg

rguments		
dirPath	A character vector giving the directory path (relative or absolute) of files to be read.	
pattern	The (grep-style) pattern describing file names to be read. The default (character(0)) reads all files in dirPath. All files are expected to have identical numbers of columns.	
colClasses	A list of length equal to the number of columns in a file. Columns with corresponding colClasses equal to NULL are ignored. Other entries in colClasses are expected to be character strings describing the base class for the XStringSet. For instance a column of DNA sequences would be specified as "DNAString". The column would be parsed into a DNAStringSet object.	
nrows	A length 1 integer vector describing the maximum number of XString objects to read into the set. Reads may come from more than one file when dirPath and pattern parse several files and nrow is greater than the number of reads in the first file.	
skip	A length 1 integer vector describing how many lines to skip at the start of each file.	
sep	A length 1 character vector describing the column separator.	
header	A length 1 logical vector indicating whether files include a header line identifying columns. If present, the header of the first file is used to name the returned values.	
comment.char	A length 1 character vector, with a single character that, when appearing at the start of a line, indicates that the entire line should be ignored. Currently there is no way to use comment characters in other than the first position of a line.	

Value

A list, with each element containing an XStringSet object of the type corresponding to the non-NULL elements of colClasses.

48 renewable

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

Examples

```
## valid character strings for colClasses
names(slot(getClass("XString"), "subclasses"))

dirPath <- system.file('extdata', 'maq', package='ShortRead')

colClasses <- rep(list(NULL), 16)
colClasses[c(1, 15, 16)] <- c("BString", "DNAString", "BString")

## read one file
readXStringColumns(dirPath, "out.aln.1.txt", colClasses=colClasses)

## read all files into a single object for each column
res <- readXStringColumns(dirPath, colClasses=colClasses)</pre>
```

renewable

Renew (update) a ShortRead object with new values

Description

Use renew to update an object defined in **ShortRead** with new values. Discover update-able classes and values with renewable.

Usage

```
renewable(x, ...)
renew(x, ...)
```

Arguments

x For renewable: missing, character(1), or a class defined in the ShortRead package. For renew: an instance of a class defined in the ShortRead package.
 ... For renewable, ignored. For renew, named arguments identifying which parts of x are to be renewed.

Details

When invoked with no arguments renewable returns a character vector naming classes that can be renewed.

When invoked with a character(1) or an instance of a **ShortRead** class, a list of the names and values of the elements that can be renewed. When x is a character vector naming a virtual class, then each element of the returned list is a non-virtual descendant of that class that can be used in renewal. This is not fully recursive.

renew is always invoked with the x argument being an instance of a class identified by renewable(). Remaining arguments are name-value pairs identifying the components of x that are to be renewed

report 49

(updated). The name-value pairs must be consistent with renewable(x). The resulting object is checked for validity. Multiple components of the object can be updated in a single call to renew, allowing comparatively efficient complex transformations.

Value

renewable() returns a character vector of renewable classes.

renewable(x) returns a named list. The names correspond to renewable classes, and the elements of the list correspond to renewable components of the class.

renew(x, ...) returns an object of the same class as x, but with components of x replaced by the named values of

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

Examples

```
## discovery
renewable()
renewable("AlignedRead")
renewable("QualityScore") ## instantiable classes
## example data
sp <- SolexaPath(system.file("extdata", package="ShortRead"))</pre>
ap <- analysisPath(sp)</pre>
filt <- chromosomeFilter("chr[[:digit:]+].fa")</pre>
aln <- readAligned(ap, "s_2_export.txt", "SolexaExport",</pre>
                    filter=filt)
## renew chromosomes from 'chr1.fa' to 'chr1', etc
labels <- sub("\\.fa", "", levels(chromosome(aln)))</pre>
renew(aln, chromosome=factor(chromosome(aln), labels=labels))
## multiple changes -- update chromosome, offset position
renew(aln, chromosome=factor(chromosome(aln), labels=labels),
      position=1L+position(aln))
## oops! invalid instances cannot be constructed
try(renew(aln, position=1:10))
```

report

Summarize quality assessment results into a report

Description

This generic function summarizes results from evaluation of qa into a report. Available report formats vary depending on the data analysed.

Usage

```
report(x, ..., dest=tempfile(), type="html")
report_html(x, dest, type, ...)
```

50 report

Arguments

x An object returned by qa, usually derived from class .QA

... Additional arguments used by specific methods.

All methods with type="html" support the argument cssFile, which is a named, length 1 character vector. The value is a path to a CSS file to be incorporated into the report (e.g., system.file("template", "QA.css", package="ShortRead")). The name of cssFile is the name of the CSS file as seen by the html report (e.g.,

"QA.css").

See specific methods for details on additional . . . arguments.

dest The output destination for the final report. For type="html" this is a directory;

for (deprecated) type="pdf" this is a file.

type A text string defining the type of report; available report types depend on the

type of object x; usually this is "html".

Details

report_html is meant for use by package authors wishing to add methods for creating HTML reports; users should always invoke report.

The following methods are defined:

x="BAMQA", ..., dest=tempfile(), type="html" Produce an HTML-based report from an object of class BAMQA.

x="BowtieQA", ..., dest=tempfile(), type="html" Produce an HTML-based report from an object of class BowtieQA.

x="FastqQA", ..., dest=tempfile(), type="html" Produce an HTML-based report from an object of class FastqQA.

x="MAQMapQA", ..., dest=tempfile(), type="html" Produce an HTML-based report from an object of class MAQMapQA.

x="SolexaExportQA", ..., dest=tempfile(), type="html" Produce an HTML-based report from an object of class SolexaExportQA.

x="SolexaExportQA", ..., dest=tempfile(), type="pdf" (Deprecated) Produce an PDF report from an object of class SolexaExportQA.

x="SolexaPath", ..., dest=tempfile(), type="html" Produce an HTML report by first visiting all _export.txt files in the analysisPath directory of x to create a SolexaExportQA instance.

x="SolexaPath", ..., dest=tempfile(), type="pdf" (Deprecated) Produce an PDF report by first visiting all _export.txt files in the analysisPath directory of x to create a SolexaExportQA instance.

x="ANY", ..., dest=tempfile(), type="ANY" This method is used internally

Value

This function is invoked for its side effect; the return value is the name of the directory or file where the report was created.

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

RochePath-class 51

See Also

```
SolexaExportQA
```

Examples

RochePath-class

"RochePath" class representing a Roche (454) experiment location

readPath = experimentPath,

qualPath = readPath,

Description

This class represents the directory location where Roche (454) result files (fasta sequences and qualities) can be found.

Objects from the Class

```
Objects from the class are created with the RochePath constructor:
```

RochePath(experimentPath = NA_character_,

experimentPath character(1) or RochePath pointing to the top-level directory of a Roche experiment.

readPath character() of directories (typically in experimentPath) containing sequence (read) information. The default selects all directories matching list.files(experimentPath, "run").

qualPath character() of directories (typically in experimentPath) containing quality information. The default selects all directories matching list.files(experimentPath, "run").

verbose logical(1) indicating whether invalid paths should be reported interactively.

Slots

RocheSet has the following slots:

```
readPath: Object of class "character", as described in the constructor, above. qualPath: Object of class "character", as described in the constructor, above. basePath: Object of class "character", containing the experimentPath.
```

Extends

```
Class "ExperimentPath", directly. Class ".Roche", directly. Class ".ShortReadBase", by class "ExperimentPath", distance 2. Class ".ShortReadBase", by class ".Roche", distance 2.
```

52 RochePath-class

Methods

```
RochePath has the following methods or functions defined:
```

```
readFasta signature(dirPath = "RochePath",
                                                      pattern=".\.fna$", sample = 1, run = 1, ...):
    Read sequences from files matching list.files(dirPath, pattern) (when dirPath="character")
    or list.files(readPath(dir)[run], pattern)[sample]. The result is a DNAStringSet.
readQual signature(dirPath = "RochePath", reads=NULL,
                                                                   pattern="\.qual$", sample=1, run=1,
    Read quality scores from files matching list.files(qualPath(dirPath)[run])[sample].
    Non-null reads is used as an (optional) template for parsing quality scores.
readFastaQual signature(dirPath = "RochePath",
                                                           fastaPattern = "\.fna$", qualPattern = "\.
    read sequences and quality scores into a ShortReadQ instance.
readFastaQual signature(dirPath = "character",
                                                           fastaPattern = "\.fna$", qualPattern = "\.
    wrapper for method above, coercing dirPath to a RochePath via RochePath(dirPath).
readBaseQuality signature(dirPath = "RochePath", ...): Reads in base and quality infor-
    mation. Currently delegates to readFastaQual, above, but will do more after RochePath
    supports more file types.
read454 signature(dirPath = "RochePath", ...): Pass arguments on to readFastaQual,
    documented above.
readPath signature(object = "RochePath"): return the contents of the readPath slot.
runNames signature(object = "RochePath"): return the basenames of readPath(object).
RocheSet signature(path = "RochePath"): create a RocheSet from path.
Additional methods include:
show signature(object = "RochePath"): Briefly summarize the experiment path locations.
detail signature(x = "RochePath"): Provide additional detail on the Roche path. All file paths
    are presented in full.
```

Author(s)

Michael Lawrence <mflawrence@fhcrc.org>

See Also

ExperimentPath.

```
showClass("RochePath")
```

RocheSet-class 53

RocheSet-class

Roche (454) experiment-wide data container

Description

This class is meant to coordinate all data in a Roche (454) experiment. See SRSet for additional details.

Objects from the Class

Create objects from this class using one of the RocheSet methods documented below

Slots

sourcePath: Object of class "RochePath" The file system location of the data used in this experiment.

readIndex: Object of class "integer" indexing reads included in the experiment; see SRSet for details on data representation in this class.

readCount: Object of class "integer" containing the number of reads associated with each sample; see SRSet for details on data representation in this class.

phenoData: Object of class "AnnotatedDataFrame" with as many rows as there are samples, containing information on experimental design.

readData: Object of class "AnnotatedDataFrame" containing as many rows as there are reads, containing information on each read in the experiment.

Extends

```
Class "SRSet", directly. Class ".Roche", directly. Class ".ShortReadBase", by class "SRSet", distance 2. Class ".ShortReadBase", by class ".Roche", distance 2.
```

Methods

No methods defined with class "RocheSet" in the signature; see SRSet for inherited methods.

Author(s)

Michael Lawrence <mflawrence@fhcrc.org>

See Also

SRSet

```
showClass("RocheSet")
```

S4 RtaIntensity

RtaIntensity

Construct objects of class "RtaIntensity"

Description

RtaIntensity objects contain Illumina image intensity measures created by the RTA pipeline. It will often be more convenient to create this object using readIntensities.

Usage

Arguments

intensity

A matrix of image intensity values. Successive columns correspond to nucleotides A, C, G, T; four successive columns correspond to each cycle. Typically, derived from "_int.txt" files.

measurementError

As intensity, but measuring standard error. Usually derived from "_nse.txt"

readInfo

An object of class AnnotatedDataFrame, containing information described by

 ${\tt RtaIntensityInfo}.$

... Additional arguments, not currently used.

Value

An object of class RtaIntensity.

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

See Also

RtaIntensity, readIntensities.

```
rta <- RtaIntensity(array(runif(60), c(5,4,3)))
intensity(rta)
## subsetting, access, and coercion
as(intensity(rta)[1:2,,], "array")</pre>
```

RtaIntensity-class 55

RtaIntensity-class

Class "RtaIntensity"

Description

Subclass of Intensity for representing image intensity data from the Illumina RTA pipeline.

Objects from the Class

Objects can be created by calls to RtaIntensity or more usually readIntensities.

Slots

Object of RtaIntensity have slots:

readInfo: Object of class "RtaIntensityInfo" representing information about each read.

intensity: Object of class "ArrayIntensity" containing an array of intensities with dimensions read, base, and cycle. Nucleotide are A, C, G, T for each cycle.

measurementError: Object of class "ArrayIntensity" containing measurement errors for each read, cycle, and base, with dimensions like that for intensity.

.hasMeasurementError: Object of class "ScalarLogical" used internally to indicate whether measurement error information is included.

Extends

```
Class "SolexaIntensity", directly.

Class "Intensity", by class "SolexaIntensity", distance 2.

Class ".ShortReadBase", by class "SolexaIntensity", distance 3.
```

Methods

Class "RtaIntensity" inherits accessor, subsetting, and display methods from class SolexaIntensity.

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

See Also

```
SolexaIntensity, readIntensities
```

```
showClass("RtaIntensity")
showMethods(class="RtaIntensity", where=getNamespace("ShortRead"))
```

56 ShortRead-class

ShortRead-class

"ShortRead" class for short reads

Description

This class provides a way to store and manipulate, in a coordinated fashion, uniform-length short reads and their identifiers.

Objects from the Class

Objects from this class are created by readFasta, or by calls to the constructor ShortRead, as outlined below.

Slots

sread: Object of class "DNAStringSet" containing IUPAC-standard, uniform-length DNA strings represent short sequence reads.

id: Object of class "BStringSet" containing identifiers, one for each short read.

Extends

Class ". ShortReadBase", directly.

Methods

Constructors include:

ShortRead signature(sread = "DNAStringSet", id = "BStringSet"): Create a ShortRead object from reads and their identifiers. The length of id must match that of sread.

ShortRead signature(sread = "DNAStringSet", id = "missing"): Create a ShortRead object from reads, creating empty identifiers.

ShortRead signature(sread = "missing", id = "missing", ...): Create an empty ShortRead object.

Methods include:

```
sread signature(object = "AlignedRead"): access the sread slot of object.
```

id signature(object = "AlignedRead"): access the id slot of object.

[signature(x = "ShortRead", i = "ANY", j = "missing"): This method creates a new ShortRead object containing only those reads indexed by i. Additional methods on '[,ShortRead' do not provide additional functionality, but are present to limit inappropriate use.

append signature(x = "ShortRead", values = "ShortRead", length = "missing"): append the sread and id slots of values after the corresponding fields of x.

narrow signature(x = "ShortRead", start = NA, end = NA, width = NA, use.names = TRUE):
 'narrow' sread so that sequences are between start and end bases, according to narrow in
 the IRanges package.

length signature(x = "ShortRead"): returns a integer(1) vector describing the number of
 reads in this object.

ShortRead-class 57

```
width signature(x = "ShortRead"): returns an integer() vector of the widths of each read in
         this object.
    srorder signature(x = "ShortRead"):
    srrank signature(x = "ShortRead"):
    srsort signature(x = "ShortRead"):
    srduplicated signature(x = "ShortRead"): Order, rank, sort, and find duplicates in ShortRead
         objects based on sread(x), analogous to the corresponding functions order, rank, sort, and
         duplicated, ordering nucleotides in the order ACGT.
    srdistance signature(pattern="ShortRead", subject="ANY"): Find the edit distance between
         each read in pattern and the (short) sequences in subject. See srdistance for allowable
         values for subject, and for additional details.
    trimLRPatterns signature(Lpattern = "", Rpattern = "", subject = "ShortRead", max.Lmismatch = 0,
         Remove left and / or right flanking patterns from sread(subject), as described in trimLRPatterns.
         Classes derived from ShortRead (e.g., ShortReadQ, AlignedRead) have corresponding base
         quality scores trimmed, too. The class of the return object is the same as the class of subject,
         except when ranges=TRUE when the return value is the ranges to use to trim 'subject'.
    alphabetByCycle signature(stringSet = "ShortRead"): Apply alphabetByCycle to the sread
         component of stringSet, returning a matrix as described in alphabetByCycle.
    tables signature(x= "ShortRead", n = 50): Apply tables to the sread component of x, re-
         turning a list summarizing frequency of reads in x.
    clean signature(object="ShortRead"): Remove all reads containing non-nucleotide ("N", "-")
         symbols.
    show signature(object = "ShortRead"): provides a brief summary of the object, including its
         class, length and width.
    detail signature(x = "ShortRead"): provides a more extensive summary of this object, dis-
         playing the first and last entries of sread and id.
    writeFasta signature(object, file, ...): write object to file in fasta format. See writeFASTA
         for ... argument values.
Author(s)
    Martin Morgan
See Also
    ShortReadQ
Examples
```

showClass("ShortRead")

showMethods(class="ShortRead", where=getNamespace("ShortRead"))

58 ShortReadQ-class

ShortRead-deprecated Deprecated functions from the ShortRead package

Description

These functions are deprecated, and will become defunct.

Usage

```
uniqueFilter(withSread=TRUE, .name="UniqueFilter")
```

Arguments

withSread A logical(1) indicating whether uniqueness includes the read sequence (withSread=TRUE)

or is based only on chromosome, position, and strand (withSread=FALSE)

.name An optional character(1) object used to over-ride the name applied to default

filters.

Details

See srFilter for details of ShortRead filters.

uniqueFilter selects elements satisfying !srduplicated(x) when withSread=TRUE, and !(duplicated(chromosome duplicated(strand(x))) when withSread=FALSE.

 $The \ behavior \ when \ with Sread=TRUE \ can \ be \ obtained \ with \ occurrence Filter (with Sread=TRUE).$

The behavior when with Sread=FALSE can be obtained using a custom filter

ShortReadQ-class "ShortReadQ" class for short reads and their quality scores

Description

This class provides a way to store and manipulate, in a coordinated fashion, the reads, identifiers, and quality scores of uniform-length short reads.

Objects from the Class

Objects from this class are the result of readFastq, or can be constructed from DNAStringSet, QualityScore, and BStringSet objects, as described below.

Slots

Slots sread and id are inherited from ShortRead. An additional slot defined in this class is:

quality: Object of class "BStringSet" representing a quality score (see readFastq for some discussion of quality score).

ShortReadQ-class 59

Extends

Class "ShortRead", directly. Class ". ShortReadBase", by class "ShortRead", distance 2.

Methods

```
Constructors include:
```

scores falling at or below a.

```
\textbf{ShortReadQ} \text{ signature}(\textbf{sread = "DNAStringSet", quality = "QualityScore", id = "BStringSet"):}
ShortReadQ signature(sread = "DNAStringSet", quality = "BStringSet", id = "BStringSet"):
     Create a ShortReadQ object from reads, their quality scores, and identifiers. When quality
     is of class BStringSet, the type of encoded quality score is inferred from the letters used in
     the scores. The length of id and quality must match that of sread.
ShortReadQ signature(sread = "DNAStringSet", quality = "QualityScore", id = "missing"):
ShortReadQ signature(sread = "DNAStringSet", quality = "BStringSet", id = "missing"):
     Create a ShortReadQ object from reads and their quality scores, creating empty identifiers.
     When quality is of class BStringSet, the type of encoded quality score is inferred from the
     letters used in the scores.
ShortReadQ signature(sread = "missing", quality = "missing", id = "missing", ...):
     Create an empty ShortReadQ object.
See accessors for additional functions to access slot content, and ShortRead for inherited meth-
ods. Additional methods include:
quality inherited from signature(object = "ANY"): access the quality slot of object.
coerce signature(from = "SFastqQuality", to = "QualityScaledDNAStringSet"):
     (Use as(from, "QualityScaledDNAStringSet")) coerce objects of class from to class to,
     using the quality encoding implied by quality(from). See QualityScore for supported
     quality classes and their coerced counterparts.
writeFastq signature(object = "ShortReadQ", file = "character", ...):
writeFastq signature(object = "ShortReadQ", file = "FastqFile", ...): Write object
     to file in fastq format. See ?writeFastq for additional arguments mode and full.
[ signature(x = "ShortReadQ", i = "ANY", j = "missing"): This method creates a new
     ShortReadQ object containing only those reads indexed by i. Additional methods on '[,Short-
     Read' do not provide additional functionality, but are present to limit inappropriate use.
append signature(x = "ShortReadQ", values = "ShortRead", length = "missing"): ap-
     pend the sread, quality and id slots of values after the corresponding fields of x.
narrow signature(x = "ShortReadQ", start = NA, end =
                                                                   NA, width = NA, use.names = TRUE):
     narrow sread and quality so that sequences are between start and end bases, according to
     narrow in the IRanges package.
trimTailw signature(object="ShortReadQ", k="integer",a="character", halfwidth="integer", ..., ra
     trim trailing nucleotides when a window of width 2 * halfwidth + 1 contains k or more quality
```

trimTails signature(object="ShortReadQ", k="integer",a="character", successive=FALSE, ..., ranges trim trailing nucleotides if k nucleotides fall below the quality encoded by a. If successive=FALSE, the k'th failing nucleotide and all subsequent nucleotides are trimmed. If successive=TRUE, failing nucleotides must occur successively; the sequence is trimmed from the first of the successive failing nucleotides.

60 ShortReadQA-class

alphabetByCycle signature(stringSet = "ShortReadQ"): Apply alphabetByCycle to the
 sread component, the quality component, and the combination of these two components
 of stringSet, returning a list of matrices with three elements: "sread", "quality", and
 "both".

alphabetScore signature(object = "ShortReadQ"): See alphabetScore for details.

qa signature(dirPath = "ShortReadQ", lane="character",..., verbose=FALSE): Perform
 quality assessment on the ShortReadQ object using lane to identify the object and returning
 an instance of ShortReadQQA. See qa

detail signature(x = "ShortReadQ"): display the first and last entries of each of sread, id, and quality entries of object.

Author(s)

Martin Morgan

See Also

readFastq for creation of objects of this class from fastq-format files.

Examples

ShortReadQA-class

Quality assessment of fastq files and ShortReadQ objects

Description

These classes contains a list-like structure with summary descriptions derived from visiting one or more fastq files, or from a ShortReadQ object.

Objects from the Class

Objects of the class are usually produced by a qa method.

Slots

.srlist: Object of class "list", containing data frames or lists of data frames summarizing the results of qa.

Extends

```
Class "SRList", directly. Class ".QA", directly. Class ".SRUtil", by class "SRList", distance 2. Class ".ShortReadBase", by class ".QA", distance 2.
```

Snapshot-class 61

Methods

Accessor methods are inherited from the SRList class.

Additional methods defined on this class are:

```
report signature(x="FastqQA", ..., dest=tempfile(), type="html"): produces HTML
    files summarizing QA results. dest should be a directory.
```

report signature(x="ShortReadQA", ..., dest=tempfile(), type="html"): produces HTML files summarizing QA results. dest should be a directory.

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

See Also

qa.

Examples

```
showClass("FastqQA")
```

Snapshot-class

Class "Snapshot"

Description

A Snapshot-class to visualize genomic data from BAM files with zoom and pan functionalities.

Usage

```
Snapshot(files, range, ...)
```

Arguments

files A character() or BamFileList specifying the file(s) to be visualized.

A GRanges object specifing the range to be visualized. range

Additional, optional, arguments to be passed to the Snapshot initialize func-

tion. Arguments include:

functions: A SnapshotFunctionList of functions, in addition to built-in 'fine_coverage', 'coarse coverage', 'multifine coverage', to be used for visualization.

currentFunction: character(1) naming the function, from functinos to be used for data input and visualization. The default chooses a function based on the scale at which the data is being visualized.

annTrack: Annotation track. If built-in visulization functions are to be used, annTrack should be a GRanges instance and the first column of its element-Meatdata would be used to annotate the range.

fac: Character(1) indicating which factor used for grouping the sample files. The factor should be included in the elementMetadata of files, otherwise ignored. Used only to visualize multiple files.

.auto_display: logical(1) indicating whether the visualization is to be updated when show is invoked.

.debug logical(1) indicating whether debug messages are to be printed.

62 Snapshot-class

Methods

```
zoom signature(x = "Snapshot"): Zoom (in or out) the current plot.
```

pan signature(x = "Snapshot"): Pan (right or left) the current plot.

togglefun signature(x = "Snapshot"): Toggle the current functions which imported records are to be immediately evaluated. Note that the active range will be changed to the current active window.

```
togglep signature(x = "Snapshot"): Toggle the panning effects.
```

togglez signature(x = "Snapshot"): Toggle the zooming effects.

Accessors

```
show signature(object = "Snapshot"): Display a Snapshot object.
```

files signature(x = "Snapshot"): Get the files field (object of class BamFileList) of a Snapshot object.

functions signature(x = "Snapshot"): Get the functions field (object of SnapshotFunctionList)
 of a Snapshot object.

vrange signature(x = "Snapshot"): Get the .range field (object of GRanges) of a Snapshot
 object.

getTrellis signature(x = "Snapshot"): Get the trellis object, a field of the SpTrellis object.

Fields

- .debug: Object of class function to display messages while in debug mode
- .auto_display: Object of class logical to automatically disply the coverage plot.
- .range: Object of class GRanges indicating which ranges of records to be imported from BAM fiels.
- .zin: Object of class logical indicating whether the current zooming effect is zoom in.
- .pright: Object of class logical indicating whether the current panning effect is right.
- .data: Object of class data.frame containing coverage a position is represented for each strand and BAM file.
- $. \, {\tt data_dirty:} \,\, {\tt Object} \,\, {\tt of} \,\, {\tt class} \,\, {\tt logical} \,\, {\tt indicating} \,\, {\tt whether} \,\, {\tt to} \,\, {\tt re-evaluate} \,\, {\tt the} \,\, {\tt imported} \,\, {\tt records}.$
- .initial_functions: Object of class SnapshotFunctionList available by the Snapshot object.
- .current_function: Object of class character of the function the imported recored are currently evaluated and visualized.
- annTrack: Default to NULL if not intended to visulize the annotation track. If default visulization function(s) is intended to be used to plot the annotation, annTrack has to be a GRanges instance.
- functions: Object of class SnapshotFunctionList of customized functions to evaluate and visualize the imported records.
- files: Object of class BamFileList to be imported.
- view: Object of class SpTrellis that is essentially a reference class wrapper of Trellis objects.

Snapshot-class 63

Class-Based Methods

```
display(): Display the current Snapshot object.
pan(): Pan (right or left) the current plot.
zoom(): Zoom (in or out) the current plot.
toggle(zoom, pan, currentFunction): Toggle zooming, panning effects or the currentFuction in which the imported records are to be evaluated and visualized.
```

Author(s)

Martin Morgan and Chao-Jen Wong < cwon2@fhcrc.org>

See Also

SpTrellis

```
## example 1: Importing specific ranges of records
file <- system.file("extdata", "SRR002051.chrI-V.bam",</pre>
                     package="yeastNagalakshmi")
which <- GRanges("chrI", IRanges(1, 2e5))</pre>
s <- Snapshot(file, range=which)</pre>
## methods
zoom(s) # zoom in
## zoom in to a specific region
zoom(s, range=GRanges("chrI", IRanges(7e4, 7e4+8000)))
pan(s) # pan right
togglez(s) # change effect of zooming
zoom(s) # zoom out
togglep(s) # chage effect of panning
pan(s)
## accessors
functions(s)
vrange(s)
show(s)
ignore.strand(s)
view(s) ## extract the spTrellis object
getTrellis(s) ## extract the trellis object
## example 2: ignore strand
s <- Snapshot(file, range=which, ignore.strand=TRUE)</pre>
##
## example 3: visulizing annotation track
library(GenomicFeatures)
getAnnGR <- function(txdb, which) {</pre>
    ex <- exonsBy(txdb, by="gene")</pre>
    sn <- keepSeqlevels(ex, seqlevels(which))</pre>
```

```
r <- range(sn)
    gr <- unlist(r)</pre>
    values(gr)[["gene_id"]] <- rep.int(names(r), times=elementLengths(r))</pre>
}
txdbFile <- system.file("extdata", "sacCer2_sgdGene.sqlite",</pre>
                     package="yeastNagalakshmi")
# txdb <- makeTranscriptDbFromUCSC(genome="sacCer2",</pre>
                                     tablename="sgdGene")
txdb <- loadDb(txdbFile)</pre>
which <- GRanges("chrI", IRanges(1, 2e5))</pre>
gr <- getAnnGR(txdb, which)</pre>
## note that the first column of the elementMetadata annotates of the
## range of the elements.
gr
s <- Snapshot(file, range=which, annTrack=gr)</pre>
annTrack(s)
## zoom in to an interesting region
zoom(s, range=GRanges("chrI", IRanges(7e4, 7e4+8000)))
togglez(s) ## zoom out
zoom(s)
pan(s)
## example 4, 5, 6: muliple BAM files with 'multicoarse_covarage'
## and 'mutifine_coverage' view.
## Resolution does not automatically switch for views of multiple
## files. It is important to note if width(which) < 10,000, use
## multifine_coverage. Otherwise use multicoarse_coverage
file <- system.file("extdata", "SRR002051.chrI-V.bam",</pre>
                     package="yeastNagalakshmi")
which <- GRanges("chrI", IRanges(1, 2e5))</pre>
s <- Snapshot(c(file, file), range=which,</pre>
               currentFunction="multicoarse_coverage")
## grouping files and view by 'multicoarse_coverage'
bfiles <- BamFileList(c(a=file, b=file))</pre>
values(bfiles) <- DataFrame(sampleGroup=factor(c("normal", "tumor")))</pre>
values(bfiles)
s <- Snapshot(bfiles, range=which,</pre>
               currentFunction="multicoarse_coverage", fac="sampleGroup")
## grouping files and view by 'multifine_coverage'
which <- GRanges("chrI", IRanges(7e4, 7e4+8000))</pre>
s <- Snapshot(bfiles, range=which,</pre>
               currentFunction="multifine_coverage", fac="sampleGroup")
```

SnapshotFunction-class

Class "SnapshotFunction"

Description

A class to store custom reader and viewer functions for the Snapshot class.

Usage

```
SnapshotFunction(reader, viewer, limits, ...) reader(x, ...) viewer(x, ...) limits(x, ...)
```

Arguments

reader	A function for reading data. The function must take a single argument (a Snapshot instance) and return a data. frame summarizing the file.
viewer	A function for visualizing the data. The function must accept the data.frame created by reader, and return an SpTrellis object representing the view.
limits	An integer(2) indicating the minimum and maximum number of nucleotides the SnapshotFunction is intended to visualize. For instance, a 'fine-scale' viewer displaying a pileup might be appropriate at between 1000 and 50000 nucleotides.
Х	An instance of SnapshotFunction
	Additional arguments, currently unused.

Fields

```
reader: Object of class function for reading data from BAM files and returning a data.frame. viewer: Object of class function for visulazation that returns an SpTrellis object. limits: Object of class integer for the limits of ranges to be visuzlized.
```

Author(s)

Martin Morgan and Chao-Jen Wong

See Also

Snapshot

```
## internally defined function
reader(ShortRead:::.fine_coverage)
viewer(ShortRead:::.fine_coverage)
limits(ShortRead:::.fine_coverage)
```

SolexaExportQA-class Quality assessment summaries from Solexa export and realign files

Description

This class contains a list-like structure with summary descriptions derived from visiting one or more Solexa 'export' or 'realign' files.

Objects from the Class

Objects of the class are usually produced by a qa method.

Slots

.srlist: Object of class "list", containing data frames or lists of data frames summarizing the results of qa.

Extends

```
Class "SRList", directly. Class ".QA", directly. Class ".SRUtil", by class "SRList", distance 2. Class ".ShortReadBase", by class ".QA", distance 2.
```

Methods

Accessor methods are inherited from the SRList class.

Additional methods defined on this class are:

```
report signature(x="SolexaExportQA", ..., dest=tempfile(), type="html"): produces HTML files summarizing QA results. dest should be a directory.
```

```
report signature(x="SolexaExportQA", ..., dest=tempfile(), type="pdf"): (deprecated;
    use type="html" instead) produces a pdf file summarizing QA results. dest should be a file.
```

report signature(x="SolexaRealignQA", ..., dest=tempfile(), type="html"): produces
HTML files summarizing QA results. dest should be a directory.

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

See Also

qa.

```
showClass("SolexaExportQA")
```

SolexaIntensity 67

SolexaIntensity	Construct objects of class "SolexaIntensity" and "SolexaIntensityInfo"

Description

These function constructs objects of SolexaIntensity and SolexaIntensityInfo. It will often be more convenient to create these objects using parsers such as readIntensities.

Usage

Arguments

intensity	A matrix of image intensity values. Successive columns correspond to nu-
	cleotides A, C, G, T; four successive columns correspond to each cycle. Typi-

cally, derived from "_int.txt" files.

measurementError

As intensity, but measuring standard error. Usually derived from "_nse.txt"

files.

readInfo An object of class AnnotatedDataFrame, containing information described by

SolexaIntensityInfo.

lane An integer vector giving the lane from which each read is derived.

tile An integer vector giving the tile from which each read is derived.

x An integer vector giving the tile-local x coordinate of the read from which each

read is derived.

y An integer vector giving the tile-local y coordinate of the read from which each

read is derived.

... Additional arguments, not currently used.

Value

An object of class SolexaIntensity, or SolexaIntensityInfo.

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

See Also

```
{\tt SolexaIntensity}.
```

68 SolexaIntensity-class

SolexaIntensity-class Classes "SolexaIntensity" and "SolexaIntensityInfo"

Description

Instances of Intensity and IntensityInfo for representing image intensity data from Solexa experiments.

Objects from the Class

Objects can be created by calls to SolexaIntensityInfo or SolexaIntensity, or more usually readIntensities.

Slots

Object of SolexaIntensity have slots:

readInfo: Object of class "SolexaIntensityInfo" representing information about each read.

intensity: Object of class "ArrayIntensity" containing an array of intensities with dimensions read, base, and cycle. Nucleotide are A, C, G, T for each cycle.

measurementError: Object of class "ArrayIntensity" containing measurement errors for each read, cycle, and base, with dimensions like that for intensity.

.hasMeasurementError: Object of class "ScalarLogical" used internally to indicate whether measurement error information is included.

Object of SolexaIntensityInfo

data Object of class "data.frame", inherited from AnnotatedDataFrame.

varMetadata Object of class "data.frame", inherited from AnnotatedDataFrame.

 $\label{lem:continuous} \textbf{dimLabels} \ \ \textbf{Object of class "character"}, inherited from \ \textbf{AnnotatedDataFrame}.$

.__classVersion__ Object of class "Versions", inherited from AnnotatedDataFrame.

.init Object of class "ScalarLogical", used internally to indicate whether the user initialized this object.

Extends

Class SolexaIntensity:

Class "Intensity", directly. Class ". ShortReadBase", by class "Intensity", distance 2.

Class SolexaIntensityInfo:

Class "AnnotatedDataFrame", directly Class "IntensityInfo", directly Class "Versioned", by class "AnnotatedDataFrame", distance 2 Class ".ShortReadBase", by class "IntensityInfo", distance 2 Class "IntensityInfo", directly.

SolexaPath-class 69

Methods

Class "SolexaIntensity" inherits accessor and display methods from class Intensity. Additional methods include:

```
[ signature(x = "SolexaIntensity", i="ANY", j="ANY", k="ANY"):
Selects the ith read, jth nucleotide, and kth cycle. Selection is coordinated across intensity, measurement error, and read information.
```

Class "SolexaIntensityInfo" inherits accessor, subsetting, and display methods from class IntensityInfo and AnnotatedDataFrame.

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

See Also

readIntensities

Examples

```
showClass("SolexaIntensity")
sp <- SolexaPath(system.file('extdata', package='ShortRead'))
int <- readIntensities(sp)
int # SolexaIntensity
readInfo(int) # SolexaIntensityInfo
int[1:5,,] # read 1:5</pre>
```

SolexaPath-class

"SolexaPath" class representing a standard output file hierarchy

Description

Solexa produces a hierarchy of output files. The content of the hierarchy varies depending on analysis options. This class represents a standard class hierarchy, constructed by searching a file hierarchy for appropriately named directories.

Objects from the Class

Objects from the class are created by calls to the constructor:

```
SolexaPath(experimentPath,
```

dataPath=.solexaPath(experimentPath, "Data"),

experimentPath character(1) object pointing to the top-level directory of a Solexa run, e.g., /home/solexa/user/080220_HWI-EAS88_0004. This is the only required argument

dataPath (optional) Solexa 'Data' folder .

scanPath (optional) Solexa GoldCrest image scan path.

imageAnalysisPath (optional) Firecrest image analysis path.

baseCallPath (optional) Bustard base call path.

analysisPath (optional) Gerald analysis pipeline path.

70 SolexaPath-class

... Additional arguments, unused by currently implemented methods.

verbose=FALSE (optional) logical vector which, when TRUE results in warnings if paths do not exist.

All paths must be fully-specified.

Slots

SolexaPath has the following slots, containing either a fully specified path to the corresponding directory (described above) or NA if no appropriate directory was discovered.

```
basePath See experimentPath, above.
dataPath See above.
scanPath See above.
imageAnalysisPath See above.
baseCallPath See above.
analysisPath See above.
```

Extends

```
Class ". Solexa", directly. Class ". ShortReadBase", by class ". Solexa", distance 2.
```

Methods

Transforming methods include:

```
readIntensities signature(dirPath = "SolexaPath", pattern=character(0), run, ...):

Use imageAnalysisPath(sp)[run] as the directory path(s) and pattern=character(0) as the pattern for discovering Solexa intensity files. See readIntensities, character-method for additional parameters.
```

```
readPrb signature(dirPath = "SolexaPath", pattern=character(0), run, ...):
```

Use baseCallPath(dirPath)[run] as the directory path(s) and pattern=character(0) as the pattern for discovering Solexa 'prb' files, returning a SFastqQuality object containing the maximum qualities found for each base of each cycle.

The ... argument may include the named argument as. This influences the return value, as explained on the readPrb, character-method page.

```
readFasta signature(dirPath, pattern = character(0), ..., nrec=-1L, skip=0L)
Use analysisPath(dirPath)[run] as the directory path(s) for discovering fasta-formatted files, returning a ShortRead object. The default method reads all files into a single object.
```

```
readFastq signature(dirPath = "SolexaPath", pattern = ".*_sequence.txt",run, ..., qualityType="S
    Use analysisPath(dirPath)[run] as the directory path(s) and pattern=".*_sequence.txt"
    as the pattern for discovering fastq-formatted files, returning a ShortReadQ object. The default
    method reads all sequence files into a single object.
```

```
readBaseQuality signature(dirPath = "SolexaPath", seqPattern = ".*_seq.txt", prbPattern = "s_[1-8]

Use baseCallPath(dirPath)[run] as the directory path(s) and seqPattern=".*_seq.txt"

as the pattern for discovering base calls and prbPattern=".*_prb.txt" as the pattern for discovering quality scores. Note that the default method reads all base call and quality score files into a single object; often one will want to specify a pattern for each lane.
```

SolexaPath-class 71

```
readQseq signature(directory="SolexaPath", pattern=".*_qseq.txt.*", run, ...., filtered=FALSE):
    Use analysisPath(dirPath)[run] as the directory path and pattern=".*_qseq.txt.*"
    as the pattern for discovering read and quality scores in Solexa 'qseq' files. Data from all files
    are read into a single object; often one will want to specify a pattern for each lane. Details are
    as for readQseq,character-method.
```

readAligned signature(dirPath = "SolexaPath", pattern = ".*_export.txt.*", run, ..., filter=srFi
 Use analysisPath(dirPath)[run] as the directory path and pattern=".*_export.txt"
 as the pattern for discovering Eland-aligned reads in the Solexa 'export' file format. Note
 that the default method reads all aligned read files into a single object; often one will want to
 specify a pattern for each lane. Use an object of SRFilter to select specific chromosomes,
 strands, etc.

```
{\bf qa} signature(dirPath="SolexaPath", pattern="character(0)", run, ...):
```

Use analysisPath(dirPath)[run] as the directory path(s) and pattern=".*_export.txt" as the pattern for discovering Solexa export-formatted fileds, returning a SolexaExportQA object summarizing quality assessment. If Rmpi or parallel has been initiated, quality assessment calculations are distributed across available nodes or cores (one node per export file.)

report signature(x, ..., dest=tempfile(), type="pdf"): Use qa(x, ...) to generate quality assessment measures, and use these to generate a quality assessment report at location dest of type type (e.g., 'pdf').

SolexaSet signature(path = "SolexaPath"): create a SolexaSet object based on path.

Additional methods include:

show signature(object = "SolexaPath"): briefly summarize the file paths of object. The experimentPath is given in full; the remaining paths are identified by their leading characters.

detail signature(x = "SolexaPath"): summarize file paths of x. All file paths are presented in full.

Author(s)

Martin Morgan

```
showClass("SolexaPath")
showMethods(class="SolexaPath", where=getNamespace("ShortRead"))
sf <- system.file("extdata", package="ShortRead")</pre>
sp <- SolexaPath(sf)</pre>
sp
readFastq(sp, pattern="s_1_sequence.txt")
nfiles <- length(list.files(analysisPath(sp), "s_[1-8]_export.txt"))</pre>
library(Rmpi)
mpi.spawn.Rslaves(nslaves=nfiles)
report(qa(sp))
## End(Not run)
## Not run:
nfiles <- length(list.files(analysisPath(sp), "s_[1-8]_export.txt"))</pre>
library(parallel)
options(srapply_fapply="parallel")
report(qa(sp))
```

72 SolexaSet-class

```
## End(Not run)

SolexaSet-class "SolexaSet" coordinating Solexa output locations with sample annotations
```

Description

This class coordinates the file hierarchy produced by the Solexa 'pipeline' with annotation data contained in an AnnotatedDataFrame (defined in the **Biobase** package).

Objects from the Class

```
Objects can be created from the constructor:
```

```
SolexaSet(path, ...).
```

- path A character(1) vector giving the fully-qualified path to the root of the directory hierarchy associated with each Solexa flow cell, or an object of class SolexaPath (see SolexaPath for this method).
- ... Additional arguments, especially laneDescription, an AnnotatedDataFrame describing the content of each of the 8 lanes in the Solexa flow cell.

Slots

SolexaSet has the following slots:

```
solexaPath: Object of class "SolexaPath".
```

laneDescription: Object of class "AnnotatedDataFrame", containing information about the samples in each lane of the flow cell.

Extends

```
Class ".Solexa", directly. Class ".ShortReadBase", by class ".Solexa", distance 2.
```

Methods

solexaPath signature(object = "SolexaSet"): Return the directory paths present when this
object was created as a SolexaPath.

laneNames signature(object = "SolexaSet"): Return the names of each lane in the flow cell, currently names are simply 1:8.

show signature(object = "SolexaSet"): Briefly summarize the experiment path and lane description of the Solexa set.

detail signature(x = "SolexaSet"): Provide additional detail on the Solexa set, including the content of solexaPath and the pData and varMetadata of laneDescription.

Methods transforming SolexaSet objects include:

```
readAligned signature(dirPath = "SolexaSet", pattern = ".*_export.txt", run, ..., filter=srFilte
    Use analysisPath(solexaPath(dirPath))[run] as the directory path(s) and pattern=".*_export.txt"
    as the pattern for discovering Eland-aligned reads in the Solexa 'export' file format. Note that
    the default method reads all aligned read files into a single object; often one will want to spec-
    ify a pattern for each lane. Use an object of SRFilter to select specific chromosomes, strands,
    etc.
```

SpTrellis-class 73

Author(s)

Martin Morgan

Examples

SpTrellis-class

Class "SpTrellis"

Description

A reference class to manage the trellis graphics related component of the Snapshot functionality for visualization of genomic data.

Usage

```
SpTrellis(trellis, debug_enabled=FALSE)
```

Arguments

trellis A trellis object for storing the plot of the genome area being visualized.

debug_enabled logical(1) indicating whether class methods should report debugging infor-

mation to the user.

Fields

trellis: Object of class trellis for storing the plot information.

debug_enabled logical(1) indicating whether class methods should report debugging information to the user.

Methods

```
zi signature(x="SpTrellis"): zoom in
zo signature(x="SpTrellis"): zoom out
right signature(x="SpTrellis"): shift to the right
left signature(x="SpTrellis"): shift to the left
restore signature(x="SpTrellis"): restore to the original plot
show signature(x="SpTrellis"): show the current plot
update signature(x="SpTrellis"): update the trellis parameters of the SpTrellis object.
```

74 spViewPerFeature

Author(s)

Chao-Jen cwon2@fhcrc.org

See Also

Snapshot

Examples

```
col <- c("#66C2A5", "#FC8D62")</pre>
x = numeric(1000)
x[sample(1000, 100)] \leftarrow abs(rnorm(100))
df <- data.frame(x = c(x, -x), pos = seq(1, 1e5, length.out=1000),
                 group = rep(c("positive", "negative"), each=1000))
cv <- xyplot(x ~ pos, df, group=group, type="s",</pre>
             col=col, main="yeast chrI:1 - 2e5",
             ylab="Coverage", xlab="Coordinate",
              scales=list(y=list(tck=c(1,0)),
                          x=list(rot=45, tck=c(1,0), tick.number=20)),
              panel=function(...) {
                      panel.xyplot(...)
                      panel.grid(h=-1, v=20)
                      panel.abline(a=0, b=0, col="grey")
              })
s <- SpTrellis(cv)</pre>
zi(s)
zi(s)
left(s)
right(s)
zo(s)
restore(s)
```

spViewPerFeature

Tools to visualize genomic data

Description

Use Snapshot-class to visulatize a specific region of genomic data

Usage

Arguments

GRI

Object GRangeList containing annotation of genomic data. It can be generated by applying exonsBy() or transcriptsBy() to a TranscriptDb instance. See examples below.

name

Character(1) specifying which element in GRL to be visualized.

spViewPerFeature 75

files	Charactor() or BamFileList specifying the file(s) to be visuliazed. If multiple files, local metadata of the files can be hold by setting a DataFrame (values(files) <- DataFrame()). See examples below.
ignore.strand	Logical(1) indicating whether to ignore the strand of the genomic data.
multi.levels	Logical(1) inidicating whether to plot the coverage of multiple files on different panels. If FALSE, the mean coverage of multiple files would be plotted.
fac	Character(1) indicating which column of local metadata (elementMetatdata()) should be used to group the samples. Ignore
	Arguments used for createing a Snapshot object.

Value

A Snapshot instance

Author(s)

Chao-Jen Wong <cwon2@fhcrc.org>

See Also

Snapshot

Examples

```
## Example 1
library(GenomicFeatures)
txdbFile <- system.file("extdata", "sacCer2_sgdGene.sqlite",</pre>
                           package="yeastNagalakshmi")
## either use a txdb file quaried from UCSC or use existing TxDb packages.
txdb <- loadDb(txdbFile)</pre>
grl <- exonsBy(txdb, by="gene")</pre>
file <- system.file("extdata", "SRR002051.chrI-V.bam",</pre>
                          package="yeastNagalakshmi")
s <- spViewPerFeature(GRL=grl, name="YAL001C", files=file)</pre>
## Example 2
## multi-files: using 'BamFileList' and setting up the 'DataFrame'
## holding the phenotype data
bfiles <- BamFileList(c(file, file))</pre>
values(bfiles) <- DataFrame(sampleGroup=factor(c("normal", "tumor")))</pre>
values(bfiles)
s <- spViewPerFeature(GRL=grl, name="YAL001C",</pre>
                       files=bfiles, multi.levels=TRUE, fac="sampleGroup")
```

76 srapply

srapply	Apply-like function for distribution across MPI-based clusters.

Description

This lapply like function evaluates locally or, if **Rmpi** or **parallel** is loaded (and **Rmpi** workers spawned) and options(srapply_fapply) set appropriately (see below), across nodes in a cluster. Errors in evaluation of FUN generate warnings; results are trimmed to exclude results where the error occurs.

Usage

Arguments

X	Tasks to be distributed. X should be an object for which lapply or sapply are defined (more precisely, mpi.parLapply, mpi.parSapply, or mclapply). Performance is best when these objects are relatively small, e.g., file names, compared to the work to be done on each by FUN.
FUN	A function to be applied to each element of X. The function must have or named argument verbose in its signature. It is best if it makes no reference to variables other than those in its argument list. or in loaded packages (the ShortRead package is available on remote nodes).
	Additional arguments, passed to FUN.
fapply	An optional argument defining an lapply-like function to be used in partitioning X. See details, below.
reduce	Optional function accepting a list (the result of fapply and summarizing this. The default reports errors in function evaluation as warnings, returning the remaining values as elements of a list. See details below for additional hints.
USE.NAMES	If TRUE and if X is character, use X as names for the result unless it had names already.
verbose	Report whether evaluation is local or mpi-based; also forwarded to FUN, allowing detailed reports from remote instances.

Details

The default value for fapply is available with ShortRead:::.fapply(). It tests getOption("srapply_fapply") for value "Rmpi" or "parallel".

If Rmpi is indicated, fapply ensures that ShortRead is required on all workers, and then invokes mpi.parLapply with arguments X, FUN, ..., and verbose. The function FUN is wrapped so that errors are returned as objects of class SRError with type RemoteError. If no workers are available, the code evaluates FUN so that errors are reported as with remote evaluation.

If **parallel** is indicated, fapply invokes mclapply with arguments as for mpi.parLapply.

Custom reduce functions might be written as reduce=function(lst) unlist(lst, use.names=TRUE).

srdistance 77

Value

The returned value depends on the value of reduce, but by default is a list with elements containing the results of FUN applied to each of X. Evaluations resulting in an error have been removed, and a warning generated.

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

Examples

```
## ... or 'verbose' required in argument,
srapply(1:10, function(i, ...) i)
## collapse result to vector
srapply(1:10, function(i, ...) i, reduce=unlist)
x <- srapply(1:10, function(i, ...) {
    if (runif(1)<.2) stop("oops") else i
})
length(x) ## trimmed to exclude errors</pre>
```

srdistance

Edit distances between reads and a small number of short references

Description

srdistance calculates the edit distance from each read in pattern to each read in subject. The underlying algorithm pairwiseAlignment is only efficient when both reads are short, and when the number of subject reads is small.

Usage

```
srdistance(pattern, subject, ...)
```

Arguments

An object of class DNAStringSet containing reads whose edit distance is desired.

Subject A short character vector, DNAString or (small) DNAStringSet to serve as reference.

... additional arguments, forward to srapply.

Details

The underlying algorithm performs pairwise alignment from each read in pattern to each sequence in subject. The return value is a list of numeric vectors of distances, one list element for each sequence in subject. The vector in each list element contains for each read in pattern the edit distance from the read to the corresponding subject. The weight matrix and gap penalties used to calculate the distance are structured to weight base substitutions and single base insert/deletions equally. Edit distance between known and ambiguous (e.g., N) nucleotides, or between ambiguous nucleotides, are weighted as though each possible nucleotide in the ambiguity were equally likely.

78 srduplicated

Value

A list of length equal to that of subject. Each element is a numeric vector equal to the length of pattern, with values corresponding to the minimum distance between between the corresponding pattern and subject sequences.

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

See Also

```
pairwiseAlignment
```

Examples

```
sp <- SolexaPath(system.file("extdata", package="ShortRead"))
aln <- readAligned(sp, "s_2_export.txt")
polyA <- polyn("A", 35)
polyT <- polyn("T", 35)

d1 <- srdistance(clean(sread(aln)), polyA)
d2 <- srdistance(sread(aln), polyA)
d3 <- srdistance(sread(aln), c(polyA, polyT))</pre>
```

srduplicated

Order, sort, and find duplicates in XStringSet objects

Description

These generics order, rank, sort, and find duplicates in short read objects, including fastq-encoded qualities. srorder, srrank and srsort differ from the default functions rank, order and sort in that sorting is based on an internally-defined order rather than, e.g., the order implied by LC_COLLATE.

Usage

```
srorder(x, ...)
srrank(x, ...)
srsort(x, ...)
srduplicated(x, ...)
```

Arguments

x The object to be sorted, ranked, ordered, or to have duplicates identified; see the examples below for objects for which methods are defined.

... Additional arguments available for use by methods; usually ignored.

srduplicated 79

Details

Unlike sort and friends, the implementation does not preserve order of duplicated elements. Like duplicated, one element in each set of duplicates is marked as FALSE.

srrank settles ties using the "min" criterion described in rank, i.e., identical elements are ranked equal to the rank of the first occurrence of the sorted element.

The following methods are defined, in addition to methods described in class-specific documenta-

```
srsort signature(x = "XStringSet"):
srorder signature(x = "XStringSet"):
srduplicated signature(x = "XStringSet"):
    Apply srorder, srrank, srsort, srduplicated to XStringSet objects such as those returned by sread.
srsort signature(x = "ShortRead"):
srorder signature(x = "ShortRead"):
srduplicated signature(x = "ShortRead"):
    Apply srorder, srrank, srsort, srduplicated to XStringSet objects to the sread component of ShortRead and derived objects.
```

Value

The functions return the following values:

srorder An integer vector the same length as x, containing the indices that will bring x

into sorted order.

srrank An integer vector the same length as x, containing the rank of each sequence

when sorted.

srsort An instance of x in sorted order.

srduplicated A logical vector the same length as x indicating whether the indexed element

is already present. Note that, like duplicated, subsetting x using the result returned by !srduplicated(x) includes one representative from each set of

duplicates.

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

Examples

```
showMethods("srsort")
showMethods("srorder")
showMethods("srduplicated")

sp <- SolexaPath(system.file('extdata', package='ShortRead'))
rfq <- readFastq(analysisPath(sp), pattern="s_1_sequence.txt")

sum(srduplicated(sread(rfq)))
srsort(sread(rfq))
srsort(quality(rfq))</pre>
```

80 srFilter

srFilter

Functions for user-created and built-in ShortRead filters

Description

These functions create user-defined (srFitler) or built-in instances of SRFilter objects. Filters can be applied to objects from ShortRead, returning a logical vector to be used to subset the objects to include only those components satisfying the filter.

Usage

```
srFilter(fun, name = NA_character_, ...)
## S4 method for signature 'missing'
srFilter(fun, name=NA_character_, ...)
## S4 method for signature 'function'
srFilter(fun, name=NA_character_, ...)
compose(filt, ..., .name)
idFilter(regex=character(0), fixed=FALSE, exclude=FALSE,
         .name="idFilter")
chromosomeFilter(regex=character(0), fixed=FALSE, exclude=FALSE,
                 .name="ChromosomeFilter")
positionFilter(min=-Inf, max=Inf, .name="PositionFilter")
strandFilter(strandLevels=character(0), .name="StrandFilter")
occurrenceFilter(min=1L, max=1L,
                 withSread=c(TRUE, FALSE, NA),
                 duplicates=c("head", "tail", "sample", "none"),
                 .name=.occurrenceName(min, max, withSread,
                                       duplicates))
nFilter(threshold=0L, .name="CleanNFilter")
polynFilter(threshold=0L, nuc=c("A", "C", "T", "G", "other"),
           .name="PolyNFilter")
dustyFilter(threshold=Inf, batchSize=NA, .name="DustyFilter")
srdistanceFilter(subject=character(0), threshold=0L,
                 .name="SRDistanceFilter")
alignQualityFilter(threshold=0L, .name="AlignQualityFilter")
alignDataFilter(expr=expression(), .name="AlignDataFilter")
```

Arguments

fun	An object of class function to be used as a filter. fun must accept a single named argument x , and is expected to return a logical vector such that $x[fun(x)]$ selects only those elements of x satisfying the conditions of fun
name	A character(1) object to be used as the name of the filter. The name is useful for debugging and reference.
filt	A SRFilter object, to be used with additional arguments to create a composite filter.
.name	An optional character(1) object used to over-ride the name applied to default filters.

srFilter 81

to filter based on chromosome. The default (character(0)) performs no filterfixed logical(1) passed to grep, influencing how pattern matching occurs. exclude logical(1) which, when TRUE, uses regex to exclude, rather than include, reads. min numeric(1) numeric(1). For positionFilter, min and max define the closed interval in max which position must be found min <= position <= max. For occurrenceFilter, min and max define the minimum and maximum number of times a read occurs after the filter. strandLevels Either character(0) or character(1) containing strand levels to be selected. ShortRead objects have standard strand levels NA, "+", "-", "*", with NA meaning strand information not available and "*" meaning strand information not relevant. withSread A logical(1) indicating whether uniqueness includes the read sequence (withSread=TRUE), is based only on chromosome, position, and strand (withSread=FALSE), or only the read sequence (withSread=NA), as described for occurrenceFilter below.. duplicates Either character{1}, a function name, or a function taking a single argument. Influence how duplicates are handled, as described for occurrenceFilter below. threshold A numeric(1) value representing a minimum (srdistanceFilter, alignQualityFilter) or maximum (nFilter, polynFilter, dustyFilter) criterion for the filter. The minima and maxima are closed-interval (i.e., x >= threshold, x <= threshold for some property x of the object being filtered). nuc

Either character(0) or a character(1) regular expression used as grep(regex, chromosome(x))

A character vector containing IUPAC symbols for nucleotides or the value

"other" corresponding to all non-nucleotide symbols, e.g., N.

batchSize NA or an integer(1) vector indicating the number of DNA sequences to be

> processed simultaneously by dustyFilter. By default, all reads are processed simultaneously. Smaller values use less memory but are computationally less

efficient.

subject A character() of any length, to be used as the corresponding argument to

srdistance.

A expression to be evaluated with pData(alignData(x)). expr

Additional arguments for subsequent methods; these arguments are not currently

used.

Details

regex

srFilter allows users to construct their own filters. The fun argument to srFilter must be a function accepting a single argument x and returning a logical vector that can be used to select elements of x satisfying the filter with x[fun(x)]

The signature (fun="missing") method creates a default filter that returns a vector of TRUE values with length equal to length(x).

compose constructs a new filter from one or more existing filter. The result is a filter that returns a logical vector with indices corresponding to components of x that pass all filters. If not provided, the name of the filter consists of the names of all component filters, each separated by " o ".

82 srFilter

The remaining functions documented on this page are built-in filters that accept an argument x and return a logical vector of length(x) indicating which components of x satisfy the filter.

```
idFilter selects elements satisfying grep(regex, id(x), fixed=fixed).
```

 $chromosome Filter\ selects\ elements\ satisfying\ grep (regex,\ chromosome (x),\ fixed=fixed).$

positionFilter selects elements satisfying $min \le position(x) \le max$.

strandFilter selects elements satisfying match(strand(x), strand, nomatch=0) > 0.

occurrenceFilter selects elements that occur >=min and <=max times. withSread determines how reads will be treated: TRUE to include the sread, chromosome, strand, and position when determining occurrence, FALSE to include chromosome, strand, and position, and NA to include only sread. The default is withSread=TRUE. duplicates determines how reads with more than max reads are treated. head selects the first max reads of each set of duplicates, tail the last max reads, and sample a random sample of max reads. none removes all reads represented more than max times. The user can also provide a function (as used by tapply) of a single argument to select amongst reads.

nFilter selects elements with fewer than threshold 'N' symbols in each element of sread(x).

polynFilter selects elements with fewer than threshold copies of any nucleotide indicated by

dustyFilter selects elements with high sequence complexity, as characterized by their dustyScore. This emulates the dust command from WindowMaker software. Calculations can be memory intensive; use batchSize to process the argument to dustyFilter in batches of the specified size.

srdistanceFilter selects elements at an edit distance greater than threshold from all sequences in subject.

alignQualityFilter selects elements with alignQuality(x) greater than threshold.

alignDataFilter selects elements with pData(alignData(x)) satisfying expr. expr should be formulated as though it were to be evaluated as eval(expr, pData(alignData(x))).

Value

srFilter returns an object of SRFilter.

Built-in filters return a logical vector of length(x), with TRUE indicating components that pass the filter.

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

See Also

SRFilter.

Examples

```
sp <- SolexaPath(system.file("extdata", package="ShortRead"))
aln <- readAligned(sp, "s_2_export.txt") # Solexa export file, as example
# a 'chromosome 5' filter
filt <- chromosomeFilter("chr5.fa")
aln[filt(aln)]
# filter during input
readAligned(sp, "s_2_export.txt", filter=filt)</pre>
```

SRFilter-class 83

```
# x- and y- coordinates stored in alignData, when source is SolexaExport
xy \leftarrow alignDataFilter(expression(abs(x-500) > 200 \& abs(y-500) > 200))
aln[xy(aln)]
# both filters as a single filter
chr5xy <- compose(filt, xy)</pre>
aln[chr5xy(aln)]
# both filters as a collection
filters <- c(filt, xy)
subsetByFilter(aln, filters)
summary(filters, aln)
# read, chromosome, strand, position tuples occurring exactly once
aln[occurrenceFilter(withSread=TRUE, duplicates="none")(aln)]
# reads occurring exactly once
aln[occurrenceFilter(withSread=NA, duplicates="none")(aln)]
# chromosome, strand, position tuples occurring exactly once
aln[occurrenceFilter(withSread=FALSE, duplicates="none")(aln)]
# custom filter: minimum calibrated base call quality >20
goodq <- srFilter(function(x) {</pre>
    apply(as(quality(x), "matrix"), 1, min) > 20
}, name="GoodQualityBases")
goodq
aln[goodq(aln)]
```

SRFilter-class

"SRFilter" for representing functions operating on ShortRead objects

Description

Objects of this class are functions that, when provided an appropriate object from the ShortRead package, return logical vectors indicating which parts of the object satisfy the filter criterion.

A number of filters are built-in (described below); users are free to create their own filters, using the srFilter function.

Objects from the Class

Objects can be created through srFilter (to create a user-defined filter) or through calls to constructors for predefined filters, as described on the srFilter page.

Slots

.Data: Object of class "function" taking a single named argument x corresponding to the Short-Read object that the filter will be applied to. The return value of the filter function is expected to be a logical vector that can be used to subset x to include those elements of x satisfying the filter.

name: Object of class "ScalarCharacter" representing the name of the filter. The name is useful for suggesting the purpose of the filter, and for debugging failed filters.

84 SRFilterResult-class

Extends

Class "function", from data part. Class ".SRUtil", directly. Class "OptionalFunction", by class "function", distance 2. Class "PossibleMethod", by class "function", distance 2.

Methods

```
srFilter signature(fun = "SRFilter"): Return the function representing the underlying filter;
this is primarily for interactive use to understanding filter function; usually the filter is invoked
as a normal function call, as illustrated below
```

```
name signature(x = "SRFilter"): Return, as a ScalarCharacter, the name of the function.
show signature(object = "SRFilter"): display a brief summary of the filter
coerce signature(from = "SRFilter", to = "FilterRules"): Coerce a filter to a
    FilterRules object of length one.
c signature(x = "SRFilter", ...): Combine filters into a single FilterRules object.
```

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

See Also

srFilter for predefined and user-defined filters.

Examples

```
## see ?srFilter
```

SRFilterResult-class "SRFilterResult" for SRFilter output and statistics

Description

Objects of this class are logical vectors indicating records passing the applied filter, with an associated data frame summarizing the name, input number of records, records passing filter, and logical operation used for all filters in which the result participated.

Usage

```
SRFilterResult(x = logical(), name = NA_character_,
    input = length(x), passing = sum(x), op = NA_character_)
## S4 method for signature 'SRFilterResult, SRFilterResult'
Logic(e1, e2)
## S4 method for signature 'SRFilterResult'
name(x, ...)
stats(x, ...)
## S4 method for signature 'SRFilterResult'
show(object)
```

SRFilterResult-class 85

Arguments

x, object, e1, e2

For SRFilterResult, logical() indicating records that passed filter or, for

others, an instance of SRFilterResult class.

name character() indicating the name by which the filter is to be referred. Internally,

name, input, passing, and op may all be vectors representing columns of a

data. frame summarizing the application of successive filters.

input integer() indicating the length of the original input.

passing integer() indicating the number of records passing the filter.

op character() indicating the logical operation, if any, associated with this filter.

... Additional arguments, unused in methods documented on this page.

Objects from the Class

Objects can be created through SRFilterResult, but these are automatically created by the application of srFilter instances.

Slots

.Data: Object of class "logical" indicating records that passed the filter.

name: Object of class "ScalarCharacter" representing the name of the filter whose results are summarized. The name is either the actual name of the filter, or a combination of filter names and logical operations when the outcome results from application of several filters in a single logical expression.

stats: Object of class "data.frame" summarizing the name, input number of records, records passing filter, and logical operation used for all filters in which the result participated. The data.frame rows correspond either to single filters, or to logical combinations of filters.

Extends

```
Class "logical", from data part. Class ".SRUtil", directly. Class "vector", by class "logical", distance 2. Class "atomic", by class "logical", distance 2. Class "vectorORfactor", by class "logical", distance 3.
```

Methods

```
Logic signature(e1 = "SRFilterResult", e2 = "SRFilterResult"): logic operations on fil-
ters.
```

! signature(x = "SRFilterResult"): Negate the outcome of the current filter results

name signature(x = "SRFilterResult"): The name of the filter that the results are based on.

stats signature(x = "SRFilterResult"): a data.frame as described in the 'Slots' section of this page.

show signature(object = "SRFilterResult"): summary of filter results.

Author(s)

Martin Morgan mailto:mtmorgan@fhcrc.org

See Also

srFilter

86 SRSet-class

Examples

```
fa <- srFilter(function(x) x %% 2 == 0, "Even")
fb <- srFilter(function(x) x %% 2 == 1, "Odd")

x <- 1:10
fa(x) | fb(x)
fa(x) & fb(x)
!(fa(x) & fb(x))</pre>
```

SRSet-class

A base class for Roche experiment-wide data

Description

This class coordinates phenotype (sample) and sequence data, primarily as used on the Roche platform.

Conceptually, this class has reads from a single experiment represented as a long vector, ordered by sample. The readCount slot indicates the number of reads in each sample, so that the sum of readCount is the total number of reads in the experiment. The readIndex field is a light-weight indicator of which reads from all those available that are currently referenced by the SRSet.

Objects from the Class

Objects of this class are not usually created directly, but instead are created by a derived class, e.g., RocheSet.

Slots

sourcePath: Object of class "ExperimentPath", containing the directory path where sequence files can be found.

readIndex: Object of class "integer" indicating specific sequences included in the experiment.

readCount: Object of class "integer" containing the number of reads in each sample included in the experiment. The sum of this vector is the total number of reads.

phenoData: Object of class "AnnotatedDataFrame" describing each sample in the experiment. The number of rows of phenoData equals the number of elements in readCount.

readData: Object of class "AnnotatedDataFrame" containing annotations on all reads.

Extends

```
Class ". ShortReadBase", directly.
```

Methods

SRUtil-class 87

```
readData signature(object="SRSet"):
sourcePath signature(object="SRSet"): Retrieve the corresponding slot from object.
show signature(object = "SRSet"): display the contents of this object.
detail signature(x = "SRSet"): provide more extensive information on the object.
```

Author(s)

Michael Lawrence <mflawrence@fhcrc.org>

Examples

```
showClass("SRSet")
```

SRUtil-class

".SRUtil" and related classes

Description

These classes provide important utility functions in the **ShortRead** package, but may occasionally be seen by the user and are documented here for that reason.

Objects from the Class

Utility classes include:

- .SRUtil-class a virtual base class from which all utility classes are derived.
- SRError-class created when errors occur in **ShortRead** package code.
- SRWarn-class created when warnings occur in ShortRead package code
- SRList-class representing a list (heterogeneous collection) of objects.
- SRVector-class representing a vector (homogeneous collection, i.e., all elements of the same class) of objects.

Objects from these classes are not normally constructed by the user. However, constructors are available, as follows.

```
SRError(type, fmt, ...), SRWarn(type, fmt, ...):
```

type character(1) vector describing the type of the error. type must come from a pre-defined list of types.

fmt a sprintf-style format string for the message to be reported with the error.

... additional arguments to be interpolated into fmt.

```
SRList(...)
```

... elements of any type or length to be placed into the SRList. If the length of ... is 1 and the argument is a list, then the list itself is placed into SRList.

```
SRVector(..., vclass)
```

... elements all satisfying an is relationship with vclass, to be placed in SRVector.

88 SRUtil-class

vclass the class to which all elements in . . . belong. If vclass is missing and length(list(...)) is greater than zero, then vclass is taken to be the class of the first argument of

SRVector errors:

SRVectorClassDisagreement this error occurs when not all arguments . . . satisfy an 'is' relationship with vclass.

Slots

SRError and SRWarn have the following slots defined:

.type: Object of class "character" containing the type of error or warning. .type must come from a pre-defined list of types, see, e.g., ShortRead:::.SRError_types.

.message: Object of class "character" containing a detailed message describing the error or warning.

SRList has the following slot defined:

.srlist: Object of class "list" containing the elements in the list.

SRVector extends SRList, with the following additional slot:

vclass: Object of class "character" naming the type of object all elements of SRVector must be.

Methods

Accessors are available for all slots, and have the same name as the slot, e.g., vclass to access the vclass slot of SRVector. Internal slots (those starting with '.' also have accessors, but these are not exported e.g., ShortRead:::.type.

SRList has the following methods:

length signature(x = "SRList"): return the (integer(1)) length of the SRList.

names signature(x = "SRList"): return a character vector of list element names. The length of the returned vector is the same as the length of x.

names<- signature(x = "SRList", value = "character"): assign value as names for members of x.

[signature(x = "SRList", i = "ANY", j = "missing"): subset the list using standard R list subset paradigms.

[[signature(x = "SRList", i = "ANY", j = "missing"): select element 'i' from the list, using standard R list selection paradigms.

lapply signature(X = "SRList", FUN="ANY"): apply a function to all elements of X, with additional arguments interpreted as with lapply.

sapply signature(X = "SRList"): apply a function to all elements of X, simplifying the result if possible. Additional arguments interpreted as with sapply.

show signature(object = "SRList"): display an informative summary of the object content, including the length of the list represented by object.

detail signature(x = "SRList"): display a more extensive version of the object, as one might expect from printing a standard list in R.

SRVector inherits all methods from SRList, and has the following additional methods:

tables 89

show signature(object = "SRVector"): display an informative summary of the object content, e.g., the vector class (vclass) and length.

detail signature(x = "SRVector"): display a more extensive version of the object, as one might expect from a printing a standard R list.

Author(s)

Martin Morgan

Examples

tables

Summarize XStringSet read frequencies

Description

This generic summarizes the number of times each sequence occurs in an XStringSet instance.

Usage

```
tables(x, n=50, ...)
```

Arguments

x An object for which a tables method is defined.

n An integer (1) value determining how many named sequences will be present in the top portion of the return value.

... Additional arguments available to methods

Details

Methods of this generic summarize the frequency with which each read occurs, There are two components to the summary. The reads are reported from most common to least common; typically a method parameter controls how many reads to report. Methods also return a pair of vectors describing how many reads were represented 1, 2, ... times.

The following methods are defined, in addition to methods described in class-specific documentation:

```
tables signature(x= "XStringSet", n = 50): Apply tables to the XStringSet x.
```

90 trimTails

Value

A list of length two.

top A named integer vector. Names correspond to sequences. Values are the number

of times the corresponding sequence occurs in the XStringSet. The vector is sorted in decreasing order; methods typically include a parameter specifying the

number of sequences to return.

distribution a data.frame with two columns. nOccurrences is the number of times any

particular sequence is represented in the set (1, 2, ...). nReads is the number of

reads with the corresponding occurrence.

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

Examples

trimTails

Trim ends of reads based on nucleotides or qualities

Description

These generic functions remove leading or trailing nucleotides or qualities. trimTails and trimTailw remove low-quality reads from the right end using a sliding window (trimTailw) or a tally of (successive) nucleotides falling at or below a quality threshold (trimTails). trimEnds takes an alphabet of characters to remove from either left or right end.

Usage

trimTails 91

```
..., ranges=FALSE)
## S4 method for signature 'XStringSet'
trimEnds(object, a, left=TRUE, right=TRUE, relation=c("<=", "=="),</pre>
    ..., ranges=FALSE)
## S4 method for signature 'XStringQuality'
trimEnds(object, a, left=TRUE, right=TRUE, relation=c("<=", "=="),</pre>
    ..., ranges=FALSE)
## S4 method for signature 'FastqQuality'
trimEnds(object, a, left=TRUE, right=TRUE, relation=c("<=", "=="),</pre>
    ..., ranges=FALSE)
## S4 method for signature 'ShortRead'
trimEnds(object, a, left=TRUE, right=TRUE, relation=c("<=", "=="),</pre>
    ..., ranges=FALSE)
## S4 method for signature 'ShortReadQ'
trimEnds(object, a, left=TRUE, right=TRUE, relation=c("<=", "=="),</pre>
    ..., ranges=FALSE)
```

Arguments

object	An object for which a methods exist (e.g., ShortReadQ and derived classes); see below to discover these methods.
k	integer(1) describing the number of failing letters required to trigger trimming.
a	For trimTails and trimTailw, a character(1) with nchar(a) == 1L giving the letter at or below which a nucleotide is marked as failing.
	For trimEnds a character() with all nchar() == 1L giving the letter at or below which a nucleotide or quality scores marked for removal.
halfwidth	The half width (cycles before or after the current; e.g., a half-width of 5 would span $5 + 1 + 5$ cycles) in which qualities are assessed.
successive	logical(1) indicating whether failures can occur anywhere in the sequence, or must be successive. If successive=FALSE, then the k'th failed letter and subsequent are removed. If successive=TRUE, the first succession of k failed and subsequent letters are removed.
left, right	logical(1) indicating whether trimming is from the left or right ends.
relation	character(1) selected from the argument values, i.e., "<=" or "==" indicating whether all letters at or below the alphabet(object) are to be removed, or only exact matches.
	Additional arguments, perhaps used by methods.
alphabet	character() (ordered low to high) letters on which quality scale is measured. Usually supplied internally (user does not need to specify). If missing, then set to ASCII characters 0-127.
ranges	logical(1) indicating whether the trimmed object, or only the ranges satisfying the trimming condition, be returned.

Details

trimTailw starts at the left-most nucleotide, tabulating the number of cycles in a window of 2 * halfwidth + 1 surrounding the current nucleotide with quality scores that fall at or below a. The read is trimmed at the first nucleotide for which this number \ge k. The quality of the first or last nucleotide is used to represent portions of the window that extend beyond the sequence.

92 Utilites

trimTails starts at the left-most nucleotide and accumulates cycles for which the quality score is at or below a. The read is trimmed at the first location where this number >= k. With successive=TRUE, failing qualities must occur in strict succession.

trimEnds examines the left, right, or both ends of object, marking for removal letters that correspond to a and relation. The trimEnds, ShortReadQ-method trims based on quality.

Value

An instance of class(object) trimmed to contain only those nucleotides satisfying the trim criterion or, if ranges=TRUE an IRanges instance defining the ranges that would trim object.

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

Examples

```
showMethods(trimTails)

sp <- SolexaPath(system.file('extdata', package='ShortRead'))
rfq <- readFastq(analysisPath(sp), pattern="s_1_sequence.txt")

## remove leading / trailing quality scores <= 'I'
trimEnds(rfq, "I")

## remove leading / trailing 'N's
rng <- trimEnds(sread(rfq), "N", relation="==", ranges=TRUE)
narrow(rfq, start(rng), end(rng))

## remove leading / trailing 'G's or 'C's
trimEnds(rfq, c("G", "C"), relation="==")</pre>
```

Utilites

Utilities for common, simple operations

Description

These functions perform a variety of simple operations.

Usage

```
polyn(nucleotides, n)
```

Arguments

nucleotides A character vector with all elements having exactly 1 character, typically from the IUPAC alphabet.

n An integer(1) vector.

Details

polyn returns a character vector with each element having n characters. Each element contains a single nucleotide. Thus polyn("A", 5) returns AAAAA.

Utilites 93

Value

polyn returns a character vector of length length(nucleotide)

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

Examples

```
polyn(c("A", "N"), 35)
```

Index

!,SRFilterResult-method	qa, 24
(SRFilterResult-class), 84	qa2, <u>27</u>
*Topic IO	QualityScore, 30
readXStringColumns,47	readAligned, 33
*Topic classes	readBaseQuality,38
.QA-class,3	readFasta, 40
AlignedDataFrame-class, 6	readFastq,41
AlignedRead-class, 8	readIntensities, 43
BAMQA-class, 12	readPrb, 45
BowtieQA-class, 13	readQseq, 46
ExperimentPath-class, 18	renewable, 48
GappedReads-class, 21	report, 49
Intensity-class, 22	RtaIntensity, 54
MAQMapQA-class, 23	SolexaIntensity, 67
QA-class, 25	srapply, 76
QualityScore-class, 31	srdistance, 77
RochePath-class, 51	srduplicated,78
RocheSet-class, 53	srFilter, 80
RtaIntensity-class, 55	tables, 89
ShortRead-class, 56	trimTails,90
ShortReadQ-class, 58	Utilites, 92
ShortReadQA-class, 60	*Topic methods
Snapshot-class, 61	GappedReads-class, 21
SnapshotFunction-class, 64	*Topic package
SolexaExportQA-class, 66	ShortReadBase-package, 3
SolexaIntensity-class, 68	.QA, 12, 13, 23, 25, 29, 50, 60, 66
SolexaPath-class, 69	.QA-class, 3
SolexaSet-class, 72	.QA2-class(QA-class), 25
SpTrellis-class, 73	. Roche, <i>51</i> , <i>53</i>
SRFilter-class, 83	.Roche-class (ShortReadBase-package), 3
SRFilterResult-class, 84	. SRUtil, <i>12</i> , <i>13</i> , <i>23</i> , <i>60</i> , <i>66</i> , <i>84</i> , <i>85</i>
SRSet-class, 86	.SRUtil-class(SRUtil-class),87
SRUtil-class, 87	. ShortReadBase, 3, 8, 12, 13, 18, 22, 23, 26,
*Topic manip	31, 51, 53, 55, 56, 59, 60, 66, 68, 70,
accessors, 4	72, 86
AlignedDataFrame, 5	.ShortReadBase-class
AlignedRead, 7	(ShortReadBase-package), 3
alphabetByCycle, 10	. Solexa, 70, 72
alphabetScore, 11	.Solexa-class (ShortReadBase-package), 3
clean, 14	[,AlignedRead,ANY,ANY,ANY-method
countLines, 14	(AlignedRead-class), 8
deprecated, 15	[,AlignedRead,ANY,ANY-method
dustyScore, 17	(AlignedRead-class), 8

[,AlignedRead,ANY,missing,ANY-method (AlignedRead-class), 8	[,ShortReadQ,ANY,ANY,ANY-method (ShortReadQ-class),58
[,AlignedRead,ANY,missing-method	[,ShortReadQ,ANY,ANY-method
(AlignedRead-class), 8 [,AlignedRead,missing,ANY,ANY-method	(ShortReadQ-class), 58 [,ShortReadQ,ANY,missing,ANY-method
<pre>(AlignedRead-class), 8 [,AlignedRead,missing,ANY-method</pre>	<pre>(ShortReadQ-class), 58 [,ShortReadQ,ANY,missing-method</pre>
(AlignedRead-class), 8	(ShortReadQ-class), 58
[,AlignedRead,missing,missing,ANY-method	[,ShortReadQ,missing,ANY,ANY-method
(AlignedRead-class), 8	(ShortReadQ-class), 58
[,AlignedRead,missing,missing-method	[,ShortReadQ,missing,ANY-method
(AlignedRead-class), 8	(ShortReadQ-class), 58
[,GappedReads,ANY,ANY-method	[,ShortReadQ,missing,missing,ANY-method (ShortReadQ-class),58
(GappedReads-class), 21	[,ShortReadQ-class), 36 [,ShortReadQ,missing,missing-method
[,IntensityMeasure,ANY,ANY,ANY-method	(ShortReadQ-class), 58
(Intensity-class), 22	[,SolexaIntensity,ANY,ANY,ANY-method
[,IntensityMeasure,ANY,ANY-method	(SolexaIntensity-class), 68
(Intensity-class), 22	[,SolexaIntensity,ANY,ANY-method
[,IntensityMeasure,ANY,missing,ANY-method	(SolexaIntensity-class), 68
(Intensity-class), 22	[,SolexaIntensity,ANY,missing,ANY-method
[,IntensityMeasure,missing,ANY,ANY-method	(SolexaIntensity-class), 68
(Intensity-class), 22 [,IntensityMeasure,missing,missing,ANY-methology	d[,SolexaIntensity,missing,ANY,ANY-method
(Intensity-class), 22	(SolexaIntensity-class), 68
[,MatrixQuality,ANY,missing,ANY-method	[,SolexaIntensity,missing,missing,ANY-method
(QualityScore-class), 31	(SolexaIntensity-class), 68
[,MatrixQuality,ANY,missing-method	[[,ArrayIntensity,ANY,ANY-method
(QualityScore-class), 31	(Intensity-class), 22
[,QualityScore,ANY,missing,ANY-method	[[,MatrixQuality,ANY,missing-method
(QualityScore-class), 31	(QualityScore-class), 31
[,QualityScore,ANY,missing-method	[[,QualityScore,ANY,missing-method
(QualityScore-class), 31	(QualityScore-class), 31
[,SRList,ANY,missing,ANY-method	[[,SRList,ANY,missing-method
(SRUtil-class), 87	(SRUtil-class), 87
[,SRList,ANY,missing-method	%in%, AlignedRead, RangesList-method
(SRUtil-class), 87	(AlignedRead-class), 8
[,ShortRead,ANY,ANY,ANY-method	accessors, 4, 8, 59
(ShortRead-class), 56	alignData (accessors), 4
[,ShortRead,ANY,ANY-method	alignDataFilter (srFilter), 80
(ShortRead-class), 56	AlignedDataFrame, 5, 5, 6
[,ShortRead,ANY,missing,ANY-method	AlignedDataFrame-class, 6
(ShortRead-class), 56	AlignedRead, 6, 7, 7, 8, 34–37, 57
[,ShortRead,ANY,missing-method	AlignedRead-class, 8
(ShortRead-class), 56	alignQuality (accessors), 4
[,ShortRead,missing,ANY,ANY-method	alignQualityFilter(srFilter), 80
(ShortRead-class), 56	alphabet,FastqQuality-method
[,ShortRead,missing,ANY-method	(QualityScore-class), 31
(ShortRead-class), 56	alphabetByCycle, 10, 32, 57, 60
[,ShortRead,missing,missing,ANY-method	alphabetByCycle,BStringSet-method
(ShortRead-class), 56	(alphabetByCycle), 10
[,ShortRead,missing,missing-method	alphabetByCycle,FastqQuality-method
(ShortRead-class), 56	(QualityScore-class), 31

alphabetByCycle,ShortRead-method	<pre>class:GappedReads (GappedReads-class),</pre>
(ShortRead-class), 56	21
alphabetByCycle,ShortReadQ-method	clean, 14
(ShortReadQ-class), 58	clean, DNAStringSet-method (clean), 14
alphabetFrequency, 32	clean, ShortRead-method
alphabetFrequency,FastqQuality-method	(ShortRead-class), 56
(QualityScore-class), 31	<pre>close.ShortReadFile (FastqFile-class),</pre>
alphabetScore, 11, 32, 60	19
alphabetScore,FastqQuality-method	<pre>coerce,AlignedRead,GappedAlignments-method</pre>
(QualityScore-class), 31	(AlignedRead-class), 8
alphabetScore, SFastqQuality-method	coerce, AlignedRead, GappedReads-method
(QualityScore-class), 31	(AlignedRead-class), 8
alphabetScore, ShortReadQ-method	coerce, AlignedRead, GRanges-method
(ShortReadQ-class), 58	(AlignedRead-class), 8
analysisPath (accessors), 4	coerce, AlignedRead, RangedData-method
AnnotatedDataFrame, 6, 68, 69, 72	(AlignedRead-class), 8
annTrack (Snapshot-class), 61	coerce, AlignedRead, RangesList-method
annTrack, Snapshot-method	(AlignedRead-class), 8
(Snapshot-class), 61	coerce, FastqQuality, matrix-method
append,.ShortReadBase,.ShortReadBase,ANY-met	
(ShortReadBase-package), 3	coerce, FastqQuality, numeric-method
append, AlignedDataFrame, AlignedDataFrame, mis	
(AlignedDataFrame-class), 6	coerce, FastqQuality, PhredQuality-method
append, AlignedRead, AlignedRead, missing-metho	
(AlignedRead-class), 8	coerce, PairwiseAlignedXStringSet, AlignedRead-method
append, MatrixQuality, MatrixQuality, missing-m	
(QualityScore-class), 31	coerce, SFastqQuality, matrix-method
append, QualityScore, QualityScore, missing-met	
(QualityScore-class), 31	coerce, SFastqQuality, SolexaQuality-method
append, ShortRead, ShortRead, missing-method	(QualityScore-class), 31
(ShortRead-class), 56	coerce, ShortReadQ, QualityScaledDNAStringSet-method
append, ShortReadQ, ShortReadQ, missing-method	(ShortReadQ-class), 58
(ShortReadQ-class), 58	coerce, SRFilter, FilterRules-method
	(SRFilter-class), 83
ArrayIntensity (Intensity-class), 22	compose (srFilter), 80
ArrayIntensity-class (Intensity-class), 22	countLines, 14
	coverage, 39
atomic, 85	
BAMQA, 25, 50	coverage, AlignedRead-method
BAMQA-class, 12	(AlignedRead-class), 8
- /	data.frame, 65
baseCallPath (accessors), 4	dataPath (accessors), 4
basePath (deprecated), 15	deprecated, 15
BowtieQA, 50	detail,.ShortReadBase-method
BowtieQA-class, 13	(SRUtil-class), 87
BStringSet, 30	detail, AlignedRead-method
c,GappedReads-method	·
• • • •	(AlignedRead-class), 8
(GappedReads-class), 21	detail, ExperimentPath-method
c, SRFilter-method (SRFilter-class), 83	(ExperimentPath-class), 18
chromosome (accessors), 4	detail,QualityScore-method
chromosome, AlignedRead-method	(QualityScore-class), 31
(AlignedRead-class), 8	detail,RochePath-method
chromosomeFilter(srFilter), 80	(RochePath-class), 51

detail, ShortRead-method	FastqQuality-class
(ShortRead-class), 56	(QualityScore-class), 31
detail,ShortReadQ-method	FastqSampler, 29
(ShortReadQ-class), 58	FastqSampler (FastqFile-class), 19
detail, SolexaPath-method	FastqSampler-class (FastqFile-class), 19
(SolexaPath-class), 69	FastqSamplerList(FastqFile-class), 19
detail, SolexaSet-method	FastqSamplerList,ANY-method
(SolexaSet-class), 72	(FastqFile-class), 19
detail, SRList-method (SRUtil-class), 87	FastqSamplerList,character-method
detail, SRSet-method (SRSet-class), 86	(FastqFile-class), 19
<pre>detail, SRVector-method (SRUtil-class),</pre>	FastqSamplerList-class
87	(FastqFile-class), 19
<pre>dim, Intensity-method (Intensity-class),</pre>	FastqStreamer (FastqFile-class), 19
22	<pre>FastqStreamer-class(FastqFile-class),</pre>
dim,MatrixQuality-method	19
(QualityScore-class), 31	FastqStreamerList(FastqFile-class), 19
DNAStringSet, 21, 40	FastqStreamerList,ANY-method
dustyFilter (srFilter), 80	(FastqFile-class), 19
dustyScore, 17, 82	FastqStreamerList,character-method
dustyScore,DNAStringSet-method	(FastqFile-class), 19
(dustyScore), 17	FastqStreamerList-class
dustyScore,ShortRead-method	(FastqFile-class), 19
(dustyScore), 17	files (Snapshot-class), 61
	<pre>files, Snapshot-method (Snapshot-class),</pre>
ExperimentPath, <i>51</i> , <i>52</i> , <i>86</i>	61
ExperimentPath (ExperimentPath-class),	FilterRules, 84
18	flag (qa2), 27
experimentPath (accessors), 4	flag,.QA2-method(qa2),27
experimentPath, SRSet-method	flag, QAFrequentSequence-method (qa2), 27
(SRSet-class), 86	flag,QAReadQuality-method(qa2), 27
ExperimentPath-class, 18	flag, QASource-method (qa2), 27
•	function, 65, 84
fac (Snapshot-class), 61	functions (Snapshot-class), 61
fac, Snapshot-method (Snapshot-class), 61	functions, Snapshot-method
FastqFile (FastqFile-class), 19	(Snapshot-class), 61
FastqFile-class, 19	CommodAlimmmento 21
FastqFileList (FastqFile-class), 19	GappedAlignments, 21
FastqFileList-class (FastqFile-class),	GappedAlignments-class, 21
19	GappedReads (GappedReads-class), 21
FastqFileReader-class	GappedReads-class, 21
(FastqFile-class), 19	getTrellis (Snapshot-class), 61
FastqQA, 25, 50	getTrellis, Snapshot-method
FastqQA (ShortReadQA-class), 60	(Snapshot-class), 61
FastqQA-class (ShortReadQA-class), 60	GRanges, 9, 61
FastqQuality, 32	grep, 15, 33, 38–41, 47, 81
FastqQuality (QualityScore), 30	id (accessors), 4
FastqQuality,BStringSet-method	id, ShortRead-method (ShortRead-class),
(QualityScore), 30	56
FastqQuality,character-method	idFilter(srFilter),80
(QualityScore), 30	ignore.strand (Snapshot-class), 61
FastqQuality,missing-method	ignore.strand, Snapshot erass), or
(QualityScore), 30	(Snapshot-class), 61

<pre>imageAnalysisPath(accessors), 4</pre>	names, SRList-method (SRUtil-class), 87
IntegerQuality (QualityScore), 30	names<-,SRList,character-method
IntegerQuality-class	(SRUtil-class), 87
(QualityScore-class), 31	narrow, <i>32</i> , <i>56</i> , <i>59</i>
Intensity, 44, 55, 68, 69	narrow,FastqQuality-method
intensity (Intensity-class), 22	(QualityScore-class), 31
Intensity-class, 22	narrow,GappedReads-method
IntensityInfo, 68, 69	(GappedReads-class), 21
<pre>IntensityInfo-class (Intensity-class),</pre>	narrow,MatrixQuality-method
22	(QualityScore-class), 31
IntensityMeasure-class	narrow,ShortRead-method
(Intensity-class), 22	(ShortRead-class), 56
is, 87	narrow,ShortReadQ-method
	(ShortReadQ-class), 58
laneDescription (accessors), 4	nFilter (srFilter), 80
laneNames (accessors), 4	NumericQuality, 31, 33
laneNames, AnnotatedDataFrame-method	NumericQuality(QualityScore), 30
(SolexaSet-class), 72	NumericQuality-class
laneNames,SolexaSet-method	(QualityScore-class), 31
(SolexaSet-class), 72	
lapply, 88	occurrenceFilter(srFilter), 80
lapply,SRList,ANY-method	open.ShortReadFile(FastqFile-class), 19
(SRUtil-class), 87	OptionalFunction, 84
lapply, SRList-method (SRUtil-class), 87	overlap, 9
left,SpTrellis-method	
(SpTrellis-class), 73	pairwiseAlignment, 77, 78
length,MatrixQuality-method	pan (Snapshot-class), 61
(QualityScore-class), 31	pan, Snapshot-method (Snapshot-class), 61
length,QualityScore-method	phenoData, 86
(QualityScore-class), 31	phenoData, SRSet-method (SRSet-class), 86
length,ShortRead-method	pileup (deprecated), 15
(ShortRead-class), 56	polyn (Utilites), 92
length, SRList-method (SRUtil-class), 87	polynFilter(srFilter), 80
limits (SnapshotFunction-class), 64	position (accessors), 4
list.files, <i>15</i> , <i>24</i>	position,AlignedRead-method
Logic, SRFilterResult, SRFilterResult-method	(AlignedRead-class), 8
(SRFilterResult-class), 84	positionFilter (srFilter), 80
logical, 85	PossibleMethod, <i>84</i>
MACHINE A 25 50	OA 20
MAQMapQA, 25, 50	QA, 29
MAQMapQA (MAQMapQA-class), 23	QA (qa2), 27
MAQMapQA-class, 23	qa, 3, 12, 13, 23, 24, 24, 49, 50, 60, 61, 66
matchPattern, 28	qa, character-method (qa), 24
MatrixQuality (QualityScore), 30	qa,list-method(qa),24
MatrixQuality-class	qa, ShortReadQ-method
(QualityScore-class), 31	(ShortReadQ-class), 58
measurementError (Intensity-class), 22	qa, SolexaPath-method
nama (SDEiltar-class) 92	(SolexaPath-class), 69
name (SRFilter-class), 83	QA-class, 25
name, SRFilter-method (SRFilter-class),	qa2, 25, 27
83	qa2,FastqSampler-method (qa2), 27
name, SRFilterResult-method	qa2,QAAdapterContamination-method
(SRFilterResult-class), 84	(qa2), 27

qa2,QACollate-method(qa2),27	qseq (GappedReads-class), 21
qa2,QAFastqSource-method(qa2),27	qseq,GappedReads-method
qa2,QAFrequentSequence-method(qa2),27	(GappedReads-class), 21
qa2,QANucleotideByCycle-method(qa2),27	QualityScore, <i>11</i> , <i>30</i> , 30, <i>46</i> , <i>59</i>
qa2,QANucleotideUse-method(qa2),27	QualityScore-class, 31
qa2,QAQualityByCycle-method(qa2),27	qualPath (RochePath-class), 51
qa2,QAQualityUse-method(qa2),27	qwidth,GappedReads-method
qa2,QAReadQuality-method(qa2),27	(GappedReads-class), 21
qa2,QASequenceUse-method(qa2),27	
QAAdapterContamination, 25	RangedData, 9
QAAdapterContamination (qa2), 27	Ranges, 9
QAAdapterContamination-class	RangesList, 9
(QA-class), 25	rank, 79
QACollate, 25	rbind, 28
QACollate (qa2), 27	rbind,.QA-method(.QA-class),3
QACollate, missing-method (qa2), 27	rbind, QASummary-method (qa2), 27
QACollate, QAFastqSource-method (qa2), 27	read454 (RochePath-class), 51
QACollate-class (QA-class), 25	read454,RochePath-method
QAData (qa2), 27	(RochePath-class), 51
QAData-class (QA-class), 25	readAligned, 7, 8, 10, 25, 30, 33
QAFastqSource, 25	readAligned,character-method
QAFastqSource (qa2), 27	(readAligned), 33
QAFastqSource-class (QA-class), 25	readAligned, SolexaPath-method
QAFiltered (qa2), 27	(SolexaPath-class), 69
QAFiltered-class (QA-class), 25	readAligned, SolexaSet-method
QAFlagged (qa2), 27	(SolexaSet-class), 72
QAFlagged-class (QA-class), 25	readBamGappedReads, 21
QAFrequentSequence, 25	readBaseQuality, 38
QAFrequentSequence (qa2), 27	readBaseQuality,character-method
QAFrequentSequence-class (QA-class), 25	(readBaseQuality), 38
QANucleotideByCycle, 25	readBaseQuality,RochePath-method
QANucleotideByCycle (qa2), 27	(RochePath-class), 51
QANucleotideByCycle-class (QA-class), 25	readBaseQuality,SolexaPath-method
QANucleotideUse, 25	(SolexaPath-class), 69
QANucleotideUse (qa2), 27	readBfaToc, 15, 39
QANucleotideUse-class (QA-class), 25	readCount (SRSet-class), 86
QAQualityByCycle, 25	readData (SRSet-class), 86
QAQualityByCycle (qa2), 27	reader (SnapshotFunction-class), 64
QAQualityByCycle-class (QA-class), 25	readFasta, 40
QAQualityUse, 25	readFasta, character-method (readFasta)
QAQualityUse (qa2), 27	40
QAQualityUse-class (QA-class), 25	readFasta,RochePath-method
QAReadQuality, 25	(RochePath-class), 51
QAReadQuality (qa2), 27	readFasta, SolexaPath-method
QAReadQuality-class (QA-class), 25	(SolexaPath-class), 69
	readFastaQual (RochePath-class), 51
QASequenceUse, 25 QASequenceUse (qa2), 27	readFastaQual, character-method
	(RochePath-class), 51
QASequenceUse-class (QA-class), 25	readFastaQual,RochePath-method
QASource-class (QA-class), 25	(RochePath-class), 51
QASummary-class (QA-class), 25	readFastq, 20, 30, 41, 58, 60
<pre>qnarrow,GappedReads-method (GappedReads-class),21</pre>	readFastq,character-method(readFastq) 41
(Gappeuneaus-Class), 41	41

readFastq,FastqFile-method	27
(FastqFile-class), 19	report,QANucleotideByCycle-method
readFastq,SolexaPath-method	(qa2), 27
(SolexaPath-class), 69	report, QANucleotideUse-method (qa2), 27
readGappedAlignments, 21	report,QAQualityByCycle-method(qa2),27
readGappedReads (GappedReads-class), 21	report,QAQualityUse-method(qa2),27
readIndex (SRSet-class), 86	report,QAReadQuality-method(qa2),27
readInfo(Intensity-class), 22	report, QASequenceUse-method(qa2), 27
readIntensities, 22, 23, 43, 54, 55, 67, 69	report, QASource-method(qa2), 27
readIntensities,character-method	report,SolexaExportQA-method
(readIntensities), 43	(SolexaExportQA-class), 66
readIntensities,SolexaPath-method	report,SolexaPath-method
(SolexaPath-class), 69	(SolexaPath-class), 69
readPath (RochePath-class), 51	report_html (report), 49
readPrb, 39, 45	${\tt report_html,BAMQA-method(BAMQA-class)},$
readPrb, character-method (readPrb), 45	12
readPrb,SolexaPath-method	report_html,BowtieQA-method
(SolexaPath-class), 69	(BowtieQA-class), 13
readQseq, 46	report_html,FastqQA-method
readQseq, character-method (readQseq), 46	(ShortReadQA-class), 60
readQseq,SolexaPath-method	report_html,MAQMapQA-method
(SolexaPath-class), 69	(MAQMapQA-class), 23
readQual (RochePath-class), 51	report_html,ShortReadQQA-method
readQual,character-method	(ShortReadQA-class), 60
(RochePath-class), 51	report_html,SolexaExportQA-method
readQual,RochePath-method	(SolexaExportQA-class), 66
(RochePath-class), 51	report_html,SolexaRealignQA-method
readXStringColumns, 35, 39, 47	(SolexaExportQA-class), 66
renew (renewable), 48	restore (SpTrellis-class), 73
renew, . ShortReadBase-method	restore, SpTrellis-method
(renewable), 48	(SpTrellis-class), 73
renewable, 48	right, SpTrellis-method
renewable,.ShortReadBase-method	(SpTrellis-class), 73
(renewable), 48	RochePath, 51
renewable, character-method (renewable),	RochePath (RochePath-class), 51 RochePath-class, 51
48	RocheSet, 52, 86
renewable, missing-method (renewable), 48	RocheSet (RocheSet-class), 53
report, 24, 25, 29, 49	RocheSet, character-method
report, ANY-method (report), 49	(RochePath-class), 51
report,BAMQA-method(BAMQA-class),12	RocheSet, RochePath-method
report,BowtieQA-method	(RochePath-class), 51
(BowtieQA-class), 13	RocheSet-class, 53
report,FastqQA-method	RtaIntensity, 54, 54
(ShortReadQA-class), 60	RtaIntensity-class, 55
report,MAQMapQA-method	runNames (RochePath-class), 51
(MAQMapQA-class), 23	runNames, RochePath-method
report, QA-method (qa2), 27	(RochePath-class), 51
report,QAAdapterContamination-method	(1.001.01.01.01.00), 51
(qa2), 27	sapply, 88
report, QAFiltered-method (qa2), 27	sapply, SRList-method (SRUtil-class), 87
report, QAFlagged-method (qa2), 27	ScanBamParam, 37
report, QAFrequentSequence-method (qa2),	scanPath (accessors), 4

SFastqQuality, 70	show,IntensityMeasure-method
SFastqQuality (QualityScore), 30	(Intensity-class), 22
SFastqQuality,BStringSet-method	show, NumericQuality-method
(QualityScore), 30	(QualityScore-class), 31
SFastqQuality,character-method	show,QAAdapterContamination-method
(QualityScore), 30	(qa2), 27
SFastqQuality,missing-method	show, QACollate-method (qa2), 27
(QualityScore), 30	show, QAFastqSource-method (qa2), 27
SFastqQuality-class	show, QAFrequentSequence-method (qa2), 27
(QualityScore-class), 31	show, QAReadQuality-method (qa2), 27
ShortRead, 8, 40, 58, 59, 70, 79	show, QASummary-method (qa2), 27
ShortRead (ShortRead-class), 56	show,RochePath-method
ShortRead, DNAStringSet, BStringSet-method	(RochePath-class), 51
(ShortRead-class), 56	show, ShortRead-method
ShortRead, DNAStringSet, missing-method	(ShortRead-class), 56
(ShortRead-class), 56	<pre>show, Snapshot-method (Snapshot-class),</pre>
ShortRead, missing, missing-method	61
(ShortRead-class), 56	show, SnapshotFunction-method
ShortRead-class, 56	(SnapshotFunction-class), 64
ShortRead-deprecated, 58	show, SolexaExportQA-method
ShortReadBase-package, 3	(SolexaExportQA-class), 66
ShortReadFile-class (FastqFile-class),	show, SolexaPath-method
19	(SolexaPath-class), 69
ShortReadQ, 4, 8, 28, 39, 42, 46, 52, 57, 60,	show, SolexaSet-method
70, 91	(SolexaSet-class), 72
ShortReadQ (ShortReadQ-class), 58	show,SpTrellis-method
${\tt ShortReadQ,DNAStringSet,BStringSet,BStringSet}$	
(ShortReadQ-class), 58	<pre>show,SRFilter-method(SRFilter-class),</pre>
ShortReadQ, DNAStringSet, BStringSet, missing-model and the string Set in the stri	
(ShortReadQ-class), 58	show, SRFilterResult-method
ShortReadQ, DNAStringSet, QualityScore, BStringSet, QualityScore, QualityScore, BStringSet, QualityScore, Qualit	
(ShortReadQ-class), 58	show, SRList-method (SRUtil-class), 87
ShortRead Q, DNA String Set, Quality Score, missing	-ภษณ์พอฟิRSet-method(SRSet-class), 86
(ShortReadQ-class), 58	show, SRVector-method (SRUtil-class), 87
ShortReadQ, missing, missing, missing-method	Snapshot, <i>61</i> , <i>65</i> , <i>73</i> – <i>75</i>
(ShortReadQ-class), 58	Snapshot (Snapshot-class), 61
ShortReadQ-class, 58	Snapshot,BamFileList,GRanges-method
ShortReadQA-class, 60	(Snapshot-class), 61
ShortReadQQA, 60	Snapshot, character, GRanges-method
ShortReadQQA-class (ShortReadQA-class),	(Snapshot-class), 61
60	Snapshot, character, missing-method
show, .QA-method(.QA-class), 3	(Snapshot-class), 61
show,.ShortReadBase-method	Snapshot-class, 61
(ShortReadBase-package), 3	SnapshotFunction
show, AlignedRead-method	(SnapshotFunction-class), 64
(AlignedRead-class), 8	SnapshotFunction-class, 64
show,ExperimentPath-method	SnapshotFunctionList, 61
(ExperimentPath-class), 18	SnapshotFunctionList
show,FastqQuality-method	(SnapshotFunction-class), 64
(QualityScore-class), 31	SnapshotFunctionList, ANY-method
show,Intensity-method	(SnapshotFunction-class), 64
(Intensity-class), 22	${\tt SnapshotFunctionList,SnapshotFunction-method}$

(SnapshotFunction-class), 64	(srduplicated), 78
SnapshotFunctionList-class	sread, 79
(SnapshotFunction-class), 64	sread (accessors), 4
SolexaExportQA, 25, 50, 51, 71	SRError (SRUtil-class), 87
SolexaExportQA (SolexaExportQA-class),	SRError-class (SRUtil-class), 87
66	SRFilter, 33, 71, 72, 80, 82
SolexaExportQA-class, 66	srFilter, 42, 58, 80, 83–85
SolexaIntensity, 22, 55, 67, 67	srFilter, function-method (srFilter), 80
SolexaIntensity-class, 68	srFilter, missing-method (srFilter), 80
SolexaIntensityInfo, 67	srFilter, SRFilter-method
SolexaIntensityInfo(SolexaIntensity),	(SRFilter-class), 83
67	SRFilter-class, 83
SolexaIntensityInfo-class	SRFilterResult, 85
(SolexaIntensity-class), 68	SRFilterResult (SRFilterResult-class),
SolexaPath, 18, 24, 34, 43–46, 72	84
SolexaPath (SolexaPath-class), 69	SRFilterResult-class, 84
solexaPath (accessors), 4	SRList, 12, 13, 23, 60, 61, 66
SolexaPath-class, 69	SRList (SRUtil-class), 87
SolexaRealignQA-class	SRList-class (SRUtil-class), 87
(SolexaExportQA-class), 66	srorder (srduplicated), 78
SolexaSet, <i>34</i> , <i>71</i>	srorder, AlignedRead-method
SolexaSet (SolexaSet-class), 72	(AlignedRead-class), 8
SolexaSet, character-method	srorder, FastqQuality-method
(SolexaSet-class), 72	(QualityScore-class), 31
SolexaSet, SolexaPath-method	srorder, ShortRead-method
(SolexaPath-class), 69	(ShortRead-class), 56
SolexaSet-class, 72	srorder, XStringSet-method
sourcePath (SRSet-class), 86	(srduplicated), 78
sprintf, 87	srrank (srduplicated), 78
SpTrellis (SpTrellis - sless) 73	<pre>srrank,AlignedRead-method</pre>
SpTrellis (SpTrellis-class), 73	· -
SpTrellis-class, 73	srrank, FastqQuality-method
spViewPerFeature, 74	(QualityScore-class), 31
srapply, 24, 45, 76	srrank, ShortRead-method
srdistance, 57, 77, 81	(ShortRead-class), 56
srdistance, DNAStringSet, character-method	srrank, XStringSet-method
(srdistance), 77	(srduplicated), 78
srdistance, DNAStringSet, DNAString-method	SRSet, 53
(srdistance), 77	SRSet-class, 86
srdistance, DNAStringSet, DNAStringSet-method	srsort, 33
(srdistance), 77	srsort (srduplicated), 78
srdistance, ShortRead, ANY-method	srsort, FastqQuality-method
(ShortRead-class), 56	(QualityScore-class), 31
srdistanceFilter(srFilter), 80	srsort, ShortRead-method
srduplicated, 78	(ShortRead-class), 56
srduplicated, AlignedRead-method	srsort,XStringSet-method
(AlignedRead-class), 8	(srduplicated), 78
srduplicated, FastqQuality-method	SRUtil-class, 87
(QualityScore-class), 31	SRVector (SRUtil-class), 87
srduplicated, ShortRead-method	SRVector-class (SRUtil-class), 87
(ShortRead-class), 56	SRWarn (SRUtil-class), 87
<pre>srduplicated,XStringSet-method</pre>	SRWarn-class(SRUtil-class),87

stats (SRFilterResult-class), 84	trimTailw,XStringQuality-method
stats, SRFilterResult-method	(trimTails), 90
(SRFilterResult-class), 84	
strand, AlignedRead-method	uniqueFilter (ShortRead-deprecated), 58
(AlignedRead-class), 8	Utilites, 92
strandFilter (srFilter), 80	
Strainar litter (Strillter), 60	vclass (accessors), 4
	vector, 85
tables, <i>57</i> , 89	Versioned, <i>6</i> , <i>68</i>
tables, ShortRead-method	view (Snapshot-class), 61
(ShortRead-class), 56	view, Snapshot-method (Snapshot-class),
tables, XStringSet-method (tables), 89	61
tapply, 82	viewer (SnapshotFunction-class), 64
togglefun (Snapshot-class), 61	• •
togglefun, Snapshot-method	vrange (Snapshot-class), 61
(Snapshot-class), 61	vrange, Snapshot-method
togglep (Snapshot-class), 61	(Snapshot-class), 61
togglep, Snapshot-method	
(Snapshot-class), 61	width,FastqQuality-method
togglez (Snapshot-class), 61	(QualityScore-class), 31
togglez, Snapshot-method	width,MatrixQuality-method
	(QualityScore-class), 31
(Snapshot-class), 61	width,NumericQuality-method
trellis-class (Snapshot-class), 61	(QualityScore-class), 31
trimEnds (trimTails), 90	width,QualityScore-method
trimEnds,FastqQuality-method	(QualityScore-class), 31
(trimTails), 90	width, ShortRead-method
<pre>trimEnds,ShortRead-method(trimTails),</pre>	(ShortRead-class), 56
90	write.XStringSet, 40
<pre>trimEnds, ShortReadQ-method(trimTails),</pre>	writeFASTA, 57
90	writeFasta (readFasta), 40
trimEnds,XStringQuality-method	writeFasta,DNAStringSet-method
(trimTails), 90	
<pre>trimEnds, XStringSet-method(trimTails),</pre>	(readFasta), 40
90	writeFasta, ShortRead-method
trimLRPatterns, 57	(ShortRead-class), 56
trimLRPatterns, ShortRead-method	writeFastq, 20, 59
(ShortRead-class), 56	writeFastq(readFastq),41
trimTails, 90	writeFastq,ShortReadQ,character-method
	(ShortReadQ-class), 58
trimTails,BStringSet-method	writeFastq,ShortReadQ,FastqFile-method
(trimTails), 90	(ShortReadQ-class), 58
trimTails,FastqQuality-method	
(QualityScore-class), 31	XStringSet, <i>10</i> , <i>79</i> , <i>89</i>
trimTails,ShortReadQ-method	
(ShortReadQ-class), 58	yield, <i>20</i>
trimTails,XStringQuality-method	yield (FastqFile-class), 19
(trimTails),90	yield,FastqFileReader-method
trimTailw(trimTails),90	(FastqFile-class), 19
trimTailw,BStringSet-method	yield,FastqSampler-method
(trimTails), 90	(FastqFile-class), 19
trimTailw,FastqQuality-method	yield,FastqStreamer-method
(QualityScore-class), 31	(FastqFile-class), 19
trimTailw,ShortReadQ-method	(1 45 64, 110 61455), 17
(ShortReadQ-class), 58	zi(SpTrellis-class),73
(Siloi cheddy cidss), so	(opin cirio crass), 15

```
zi, SpTrellis-method (SpTrellis-class),
73
zo (SpTrellis-class), 73
zo, SpTrellis-method (SpTrellis-class),
73
zoom (Snapshot-class), 61
zoom, Snapshot-method (Snapshot-class),
61
```