# Analyse RT-PCR data with ddCt

Jitao David Zhang, Rudolf Biczok and Markus Ruschhaupt March 30, 2012

#### Abstract

Quantitative real–time PCR (qRT–PCR or RT–PCR for short) is a laboratory technique based on the polymerase chain reaction, and is commonly used to amplify and quantify a targeted nucleotide molecule simultaneously. The data analysis of qRT–PCR experiments can be divided in subsequent steps: importing the data, setting reference sample(s) and housekeeping gene(s), (optional) filtering, applying algorithm for the relative expression, and finally reporting the results in the form of text and figure. ddCt package implements the  $2^{-\Delta \Delta C_T}$  algorithm in a pipeline as described above and performs end–to–end analysis of qRT-PCR experiments in R and Bioconductor. In this vignette we introduce the usage of the pipeline<sup>1</sup>.

This vignette introduces the analysis of RT–PCR data with ddCt by a step–to–step guide through the pipeline. The idea is that the user could read the vignette, modify the parameters to her/his need, and then use the Sweave function in R to perform the analysis. The Sweave function converts the '.Rnw' version into a LATEX file while processing the RT–PCR data in the background. It gives the same result as if the user invokes the pipeline in the way described in the section 3.1.

The vignette begins with a short description of the values to be calculated by the ddCt. Then it explains the parameters to be specified by the user in depth.

## 1 Introduction

Several values are calculated during pipeline by Sweaving this document. The most important ones are explained here in the order of the execution. To simply the discussion, the values discussed are calculated without efficiencies (see below).

- 1. The mean or median of the technical replicates for one gene–sample pair combination is the  $C_T$  value.
- 2. For a sample A the  $C_T$  value of the housekeeping gene (or the median of the  $C_T$  values of all housekeeping genes) is subtracted from the corresponding  $C_T$  value of a gene Gene 1. The result is the  $dC_T$  value for Gene 1 and sample A. This value is used for the t-test and the  $Wilcoxon\ test$ .

<sup>&</sup>lt;sup>1</sup>For more general information about qRT–PCR and the  $2^{-\Delta\Delta C_T}$  algorithm, please follow the another vignette rtPCR released along the package ddCt.

- 3. For a gene Gene1 the  $dC_T$  value of the reference sample (or the mean of the  $dC_T$  values of all reference samples) is subtracted from the corresponding  $dC_T$  value of Gene1 and a sample B. This is called the  $ddC_T$  value for Gene1 and sample B.
- 4. The transformation  $x \to 2^{-x}$  is applied to each  $ddC_T$  value. The resulting value is called 'exprs' (used to be named as 'level', these two words are exchangable in this document).

Additionally, for each value an error is calculated. This is based on the standard error of the mean (S.E.M) if the mean is used to summarize the replicates. If the median is used for summarization of the individual Ct values, the MAD<sup>2</sup> divided by the square root of replicate number is used as an error estimate.

# 2 Prerequisites

The package ddCt requires the latest version of R to perform properly. Please check the DESCRIPTION file along the package for the requirement of packages.

To run this document and perform the analysis, you need a tab-delimited plain text file containing the individual Ct values of your experiment. This file is usually exported from the software used to measure the Ct values, for example from SDM $^{\odot}$ Software of Roche Coop.. It is important that you export the individual Ct values, since normally the software combines the replicates of the measurements to one value. See the file extdata/Experiment1.txt in the directory of the ddCt package for an example.

Once the file is ready you can go on to set parameters.

# 3 Using ddCt

#### 3.1 Invoke the pipeline

To use the functionality of the ddCt package, one can

- Writting R scripts with the functions implemented in the ddCt. It requires basic programming skills but provides the maximum flexibility. See the vignette rtPCR for a short example and the manual pages of the functions in the package.
- Calling the ddCt script in the scripts sub-directory in the ddCt installation package. This script can be invoked by either typing in the command line, or by modifying the parameters in this vigenette and then run the Sweave function. Actually this vignette is just a wrapper of the ddCt script which describes the parameters. In the vigette we discuss this approach.

When using the ddCt script, One can call the ddCt script in the scripts sub-directory in the ddCt installation package ( $PKG\_DIR$  in the following samples) via the Rscript command in the command line<sup>3</sup>:

<sup>&</sup>lt;sup>2</sup>median absolute deviation, defined as  $MAD = median|x_i - \tilde{x}|$ , where  $\tilde{x}$  is the median of x. Its relationship with standard deviation  $\sigma$  can be expressed as  $\sigma \approx 1.4826 MAD$ 

<sup>&</sup>lt;sup>3</sup>The examples shown here are for Debian Linux system, on the windows system the path separator has to be modified, however the functionality is the same

This command-line approach does not need to invoke a R session and is easy to be built into automatic analysis pipelines. It is however hard to know all the possible parameters and not easy to be used by users who are not familiar with the command line.

Alternatively, one could load the script via the source function in the R prompt. If you do so, you have to pass the parameters as a list to the ddCtExec function, which is an end-to-end function combining data import, analysis and report:

```
> scFile <- system.file("scripts/ddCt.R", package="ddCt")
> inputFile <- system.file("extdata/Experiment1.txt", package="ddCt")
> source(scFile)
> ddCtExec(list(inputFile=inputFile, ...))
```

The parameters for the ddCtExec function are all equal with parameters in the Rscript call. In this example, we use the second way to invoke the ddCt package.

```
> source(system.file("scripts", "ddCt.R", package="ddCt"))
```

#### 3.2 Parameters

Table 1 on the page 3 illustrates the usable parameters in the ddCt.R script. Not all parameters have to be set or they already have default values. Only inputFile, referenceGene and referenceSample are required parameters.

Parameter	Description	Required	Default		
inputFile	a valid SDS input file	Yes	NULL		
referenceGene	houskeeping gene to use	Yes	NULL		
referenceSample	calibration sample to use	Yes	NULL		
loadPath	directory, from which the	No	Working directory		
	input files are loaded				
savePath	output directory	No	Working directory		
confFile	load parameters from this	No	NULL		
	R script				
sampleAnnotationFile	optional annotation file	No	NULL		
columnForGrouping	columns used to group	No	NULL		
	output				
onlyFromAnnotation	only use annotated sam-	No	FALSE		
	ples				
geneAlias	replace specified gene	No	NULL		
	names				
sampleAlias	replace specified sample	No	NULL		
	names				
threshold	threshold value	No	40		
mode	the used algorithm mode	No	median		
Continued on next page					

Table 1: List of parameters

Table 1 – continued from previous page

plotMode  the kind of values which schould be ploted (Ct, level,) algorithm  efficiencies  efficiencies  efficienciesError  genesRemainInOutput  samplesRemainInOutput  genesNotInOutput  genesNotInOutput  groupingBySamples plotPerObject  groupingColor  cutoff brewerColor  cutoff brewerColor  legend ttestPerform ttestCol pairsCol samplesRemainInTTest samplesRemainInTTest samplesRemainInCor samplesRemainInCor samplesRemainInCor samplesRemainInTTest samplesRemainInCor samplesRemainIn	Table 1 – continued from previous page							
schould be ploted (Ct, level,)  algorithm  efficiencies  efficiencies  efficience vaue for each gene  error value  samplesRemainInOutput  samplesRemainInOutput  genesNotInOutput  genesNotInOutput  groupingBySamples  glotPerObject  groupingForPlot  groupingColor  cutoff  brewerColor  cutoff  brewerColor  cutoff  cutoff  brewerColor  cutoff  cuto	Parameter	Description	Required	Default				
level,)   use ddCt or ddCtWihtE   algorithm   use ddCt or ddCtWihtE   algorithm   efficiencies   efficience vaue for each   gene   error value   show only these specified   genes in output   show only these specified   genes in output   dont show these specified   genes in output   dont show these specified   genes in output   groupingBySamples   plot output by   plot output by   plot output by   generate one plot per kind   of value (Ct, level,)   groupingColor   The color for each bar   No   rull   c("#E41A1C", "#377EB8", "#4DAF4A", "#9984EA3", "#FFF00")   value   color for the brewer   No   c("Set3", "Set1", "Accent", "Dark2", "Spectral", "PuOr", "BrBG")   Trest   samplesRemainInCor   use this samples for correlation plot   samples NotlnCor   dont use this samples for No   NULL   samples NotlnCor   No   NULL	plotMode		No	c("level","Ct")				
algorithm  efficiencies  efficience vaue for each gene efficiencesError genesRemainInOutput  samplesRemainInOutput  samplesRemainInOutput  genesNotInOutput  genesNotInOutput  genesNotInOutput  samplesNotInOutput  genesi in output  dont show these specified genessamples in output  dont show these specified genessamples in output  groupingBySamples plot output by generate one plot per kind of value (Ct, level,) groupingColor  cutoff cutoff value  color for the brewer  show legend  ttestPerform  ttestCol pairsCol samplesRemainInTTest  samplesRemainInCor  samplesRemainInCor  samplesRemainInCor  samplesRemainInCor  samplesRotInCor  show legend  tool the pairs samples for correlation plot lagorithm  non non non non non non non non non n		_ ` ` `						
efficiencies efficience vaue for each gene efficienciesError error value show only these specified genessamplesRemainInOutput samplesRemainInOutput show only these specified genessamples in output dont show these specified genessamples in output genesNotInOutput dont show these specified genessamples in output dont show these specified genessamples in output groupingBySamples plot PerObject group also in plot of value (Ct, level,) groupingColor The color for each bar  cutoff brewerColor  show legend ttestPerform ttestCol pairsCol samplesRemainInTest samplesRemainInCor samplesRemainInCor  algorithm efficience vaue for each gene not efficience vaue for each No NULL NO NULL NO NULL TRUE  NO NULL NO NULL NO NULL C("#4511A1C","#377EB8","#4DAF4A", "#984EA3","#FF7F00") NULL C("Set3","Set1","Accent", "Dark2","Spectral","PuOr","BrBG") TRUE  TRUE  TRUE  TRUE  NO NULL NO N								
efficiencies efficience vaue for each gene error value show only these specified genes in output show only these specified genes in output show only these specified genes in output dont show these specified genes in output dont show these specified genes in output dont show these specified genes in output genes in output dont show these specified genes in output dont show these specified genes in output genessamples in output dont show these specified genessamples in output plot output by generate one plot per kind of value (Ct, level,) groupingForPlot group also in plot groupingColor The color for each bar No c("#E41A1C","#377EB8","#4DAF4A", "#984EA3","#FFFF00")  cutoff cutoff value No NULL ("Seat3","Set1","Accent", "Dark2","Spectral","PuOr","BrBG")  legend show legend No TRUE  ttestCol color of the TTest No FALSE  ttestCol pairsCol samplesRemainInTTest use this samples for TTest samplesNotInTTest dont use this samples for correlation plot samplesRemainInCor dont use this samples for No NULL  samplesNotInCor dont use this samples for No NULL  samplesNotInCor	algorithm		No	ddCt				
efficienciesError genesRemainInOutput samplesRemainInOutput samplesRemainInOutput genesNotInOutput genesNotInOutput samplesNotInOutput samplesNotInOutput samplesNotInOutput genesNotInOutput genessamples in output dont show these specified genessamples in output genessamples in output groupingBySamples plot output by plotPerObject generate one plot per kind of value (Ct, level,) groupingForPlot group also in plot groupingColor The color for each bar  cutoff cutoff value brewerColor color for the brewer legend ttestPerform ttestCol pairsCol samplesRemainInTTest samplesRemainInCor samplesRemainInCor samplesRemainInCor show only these specified yenessamples procified genessamples in output No NULL No FALSE No NULL No								
efficienciesError genesRemainInOutput samplesRemainInOutput samplesRemainInOutput genesNotInOutput genesNotInOutput samplesNotInOutput samplesNotInOutput groupingBySamples plotPerObject groupingForPlot groupingColor cutoff brewerColor cutoff brewerColor cutoff brewerColor cutoff testepform ttestCol pairsCol samplesRemainInTTest samplesRemainInTTest samplesRemainInCor samplesRemainInCor genes in output dont show these specified genes in output No NULL  TRUE  NO NULL  c("#E41A1C","#377EB8","#4DAF4A", "#984EA3","#FFF700") NULL  c("\$ext3","Set1","Accent", "Dark?","Spectral","PuOr","BrBG")  TRUE  TRUE  No NULL  volor for the brewer No NULL  volor for the Trest No NULL  volor of the Trest No NULL  volor of the pairs volor of the	efficiencies		No	NULL				
genesRemainInOutput samplesRemainInOutput samplesRemainInOutput genesNotInOutput samplesNotInOutput samplesNotInOutput samplesNotInOutput samplesNotInOutput genessamples in output dont show these specified genessamples in output dont show these specified genes in output dont show these specified genessamples in output plot output by plot output by generate one plot per kind of value (Ct, level,) groupingForPlot group also in plot The color for each bar  cutoff cutoff cutoff value color for the brewer  cutoff color for the brewer  show legend ttestPerform ttestCol pairsCol samplesRemainInTTest samplesRemainInCor  samplesRemainInCor  samplesRemainInCor  samplesRemainInCor  samplesRemainInCor  samplesRemainInCor  samplesNotInCor  show only these specified genes in output show obly these specified show these specified show lublu struct subject output by No TRUE subject output by No NULL subject output by No NULL No NULL No NULL No NULL No No NULL No NULL No NULL No NULL No NULL No NULL No NO NULL NO NULL NO NULL NO NULL NO NO NULL NO NULL NO NULL NO NO NULL NO NULL NO NO NULL NO NULL NO NULL NO NULL NO NO NULL NO NULL NO NO NO NULL NO NO NULL NO NO NO NULL NO NO NO NULL N								
genes in output show only these specified genessamples in output dont show these specified genessamples in output dont show these specified genessin output dont show these specified genessamples in output groupingBySamples plot output by plotPerObject generate one plot per kind of value (Ct, level,) groupingForPlot group also in plot groupingColor The color for each bar  cutoff brewerColor  cutoff color for the brewer  show legend ttestPerform ttestCol pairsCol samplesRemainInTTest samplesRemainInTTest samplesRemainInCor  genes in output show only these specified genessamples in output No NULL No TRUE No NULL ("#E41A1C","#377EB8","#4DAF4A", "#984EA3","#FF7F00") NULL No C'("Set3","Set1","Accent", "Dark2","Spectral","PuOr","BrBG") TRUE  TRUE No TRUE No TRUE No No NULL No FALSE No NULL No NO NULL No NULL No NULL No NO NULL								
samplesRemainInOutput show only these specified genessamples in output dont show these specified genes in output dont show these specified genes in output dont show these specified genes in output groupingBySamples plot output by plot output by generate one plot per kind of value (Ct, level,) groupingForPlot groupingColor The color for each bar No NULL groupingColor The color for each bar No c'"#E41A1C","#377EB8","#4DAF4A", "#984EA3","#FF7F00")  cutoff cutoff value No NULL c'"\$#377EB8","#4DAF4A", "#984EA3","#FF7F00")  cutoff brewerColor color for the brewer No c'"\$#543","Spectral","PuOr","BrBG")  legend testPerform perform a TTest color of the TTest No NULL samplesRemainInTTest samplesRemainInTTest use this samples for TTest dont use this samples for CTTest use this samples for correlation plot dont use this samples for No NULL  samplesNotInCor dont use this samples for No NULL  samplesNotInCor dont use this samples for No NULL  samplesNotInCor dont use this samples for No NULL	genesRemainInOutput		No	NULL				
genesNotInOutput dont show these specified genes in output dont show these specified genes in output dont show these specified genesamples in output groupingBySamples plot output by plot output by generate one plot per kind of value (Ct, level,) groupingForPlot group also in plot No NULL groupingColor The color for each bar No c("#E41A1C","#377EB8","#4DAF4A", "#984EA3","#FF7F00")  cutoff cutoff value No NULL c("Set3","Set1","Accent", "Dark2","Spectral","PuOr","BrBG")  legend show legend testPerform perform a TTest color of the TTest No NULL samplesRemainInTTest and suse this samples for TTest amplesNotInTTest use this samples for correlation plot samplesNotInCor dont use this samples for No NULL  samplesNotInCor No NULL		_						
genesNotInOutput dont show these specified genes in output dont show these specified genes in output dont show these specified genessamples in output groupingBySamples plot output by plot output by generate one plot per kind of value (Ct, level,) groupingForPlot group also in plot No NULL groupingColor The color for each bar No c("#E41A1C","#377EB8","#4DAF4A", "#984EA3","#FF7F00") NULL color for the brewer No NULL color for the brewer No c("Se43","Set1","Accent", "Dark2","Spectral","PuOr","BrBG") TRUE testPerform perform a TTest No FALSE testCol pairsCol color of the pairs use this samples for TTest No NULL samplesRemainInTTest use this samples for TTest use this samples for Tolor use this samples for correlation plot dont use this samples for No NULL samplesNotInCor dont use this samples for No NULL No NULL samplesRemainInCor use this samples for No NULL NO NUL	samplesRemainInOutput	show only these specified	No	NULL				
genes in output dont show these specified genessamples in output plot output by plot PerObject generate one plot per kind of value (Ct, level,) groupingForPlot groupingColor The color for each bar  cutoff brewerColor  legend ttestPerform ttestCol pairsCol samplesRemainInTTest samplesRemainInCor samplesRemainInCor  generate one plot per kind of value (No raule (Ct, level,) group also in plot The color for each bar No cutoff value No No RULL c'("#E41A1C","#377EB8","#4DAF4A", "#984EA3","#FFF700") NULL c'("Set3","Set1","Accent", "Dark2","Spectral","PuOr","BrBG") TRUE  No RALSE No NULL NO NO NULL NO NULL NO NULL NO NO NULL NO NULL NO NULL NO NULL NO NULL NO NO NULL N		genessamples in output						
samplesNotInOutput dont show these specified genessamples in output plot output by plot output by generate one plot per kind of value (Ct, level,) groupingForPlot group also in plot No NULL groupingColor The color for each bar No c'#E41A1C", "#377EB8", "#4DAF4A", "#984EA3", "#FF7F00")  cutoff cutoff value No NULL (c'Set3", "Set1", "Accent", "Dark2", "Spectral", "PuOr", "BrBG")  legend show legend No TRUE  ttestPerform perform a TTest No NULL  ttestCol color of the TTest No NULL  samplesRemainInTTest use this samples for TTest Aont use this samples for TTest  samplesRemainInCor use this samples for correlation plot dont use this samples for No NULL  samplesNotInCor dont use this samples for No NULL  samplesNotInCor dont use this samples for No NULL	genesNotInOutput	dont show these specified	No	NULL				
genessamples in output plot output by plot Object  generate one plot per kind of value (Ct, level,) groupingForPlot group also in plot groupingColor  The color for each bar  No  Cutoff brewerColor  color for the brewer  legend ttestPerform ttestCol pairsCol samplesRemainInTTest samplesRemainInTTest samplesRemainInCor  genessamples in output plot output by generate one plot per kind of value (Ct, level,) No  No  NULL No  C("#E41A1C","#377EB8","#4DAF4A", "#984EA3","#FF7F00")  NULL C("Set3","Set1","Accent", "Dark2","Spectral","PuOr","BrBG")  TRUE  No  NO  TRUE  NO  NULL NO  FALSE  NO  NULL NO  NO  NO  NULL NO  NO  NO  NO  NO  NO  NO  NO  NO  N		genes in output						
groupingBySamples plot output by plotPerObject  generate one plot per kind of value (Ct, level,) groupingForPlot group also in plot The color for each bar  No  cutoff brewerColor  color for the brewer  legend ttestPerform ttestCol pairsCol samplesRemainInTTest samplesRemainInTTest samplesRemainInCor  samplesRemainInCor  samplesRemainInCor  samplesRemainInCor  samplesRoul samplesRemainInCor  samplesNotInCor  plot output by generate one plot per kind No TRUE  No NULL No NULL No NULL No NULL No NO TRUE No NULL N	samplesNotInOutput	dont show these specified	No	NULL				
plotPerObject generate one plot per kind of value (Ct, level,) groupingForPlot group also in plot No NULL groupingColor The color for each bar No c("#E41A1C","#377EB8","#4DAF4A", "#984EA3","#FF7F00")  cutoff cutoff value No NULL brewerColor color for the brewer No c("Set3","Set1","Accent", "Dark2","Spectral","PuOr","BrBG")  legend show legend No TRUE ttestPerform perform a TTest No FALSE ttestCol color of the TTest No NULL pairsCol color of the pairs No NULL samplesRemainInTTest use this samples for TTest No NULL samplesNotInTTest dont use this samples for No NULL samplesRemainInCor use this samples for correlation plot samplesNotInCor dont use this samples for No NULL		genessamples in output						
of value (Ct, level,) groupingForPlot group also in plot The color for each bar  No c("#E41A1C","#377EB8","#4DAF4A", "#984EA3","#FF7F00")  cutoff brewerColor  color for the brewer  color for the brewer  No c("Set3","Set1","Accent", "Dark2","Spectral","PuOr","BrBG")  legend ttestPerform perform a TTest perform a TTest pairsCol color of the TTest color of the pairs samplesRemainInTTest samplesNotInTTest dont use this samples for TTest samplesNotInCor  volume No NULL NO NU	groupingBySamples	plot output by	No	TRUE				
groupingForPlot groupingColor  The color for each bar  No  c("#E41A1C","#377EB8","#4DAF4A",	plotPerObject	generate one plot per kind	No	TRUE				
groupingColor  The color for each bar  No  c("#E41A1C","#377EB8","#4DAF4A",		of value (Ct, level,)						
cutoff brewerColor  cutoff value brewerColor  color for the brewer  color for the brewer  No  c("Set3","Set1","Accent", "Dark2","Spectral","PuOr","BrBG")  REPRETATION OF TRUE  ttestPerform  perform a TTest  ttestCol  pairsCol  samplesRemainInTTest  samplesRemainInTTest  dont use this samples for  TTest  samplesRemainInCor  samplesRemainInCor  samplesRemainInCor  cutoff value  No  NULL	groupingForPlot	group also in plot	No	NULL				
cutoff brewerColor  cutoff value brewerColor  color for the brewer  No  c("Set3","\$\frac{1}{2}\$ ("PuOr","BrBG")  legend ttestPerform perform a TTest testCol pairsCol samplesRemainInTTest samplesNotInTTest  samplesRemainInCor  samplesNotInCor  cutoff value No NULL No NULL No TRUE Test No NULL	groupingColor	The color for each bar	No	c("#E41A1C","#377EB8","#4DAF4A",				
brewerColor  color for the brewer  No  c("Set3", "Set1", "Accent", "Dark2", "Spectral", "PuOr", "BrBG")  No  ttestPerform  perform a TTest  color of the TTest  pairsCol  samplesRemainInTTest  samplesNotInTTest  samplesRemainInCor  samplesNotInCor  color of the pairs  samples for correlation plot  samplesNotInCor  No  c("Set3", "Set1", "Accent", "Dark2", "Spectral", "PuOr", "BrBG")  No  TRUE  No  NULL								
legend show legend perform a TTest No FALSE ttestCol color of the TTest No NULL pairsCol color of the pairs No NULL samplesRemainInTTest use this samples for TTest No NULL samplesNotInTTest dont use this samples for TTest SamplesRemainInCor use this samples for correlation plot samplesNotInCor dont use this samples for No NULL  No NULL  No NULL  No NULL  No NULL  No NULL	cutoff	cutoff value	No	NULL				
legend show legend perform a TTest No FALSE ttestCol color of the TTest No NULL pairsCol color of the pairs No NULL samplesRemainInTTest use this samples for TTest No NULL samplesNotInTTest dont use this samples for TTest samplesRemainInCor use this samples for correlation plot samplesNotInCor dont use this samples for No NULL	brewerColor	color for the brewer	No	c("Set3","Set1","Accent",				
ttestPerform ttestCol ttestCol pairsCol samplesRemainInTTest samplesNotInTTest samplesRemainInCor samplesNotInCor  perform a TTest color of the TTest volume No volume TTest volume TTest volume No volume TTest volume No volume TTest volume No volume TTest volume TTest volume No volume TTest vo				"Dark2","Spectral","PuOr","BrBG")				
ttestCol color of the TTest No NULL No NULL samplesRemainInTTest use this samples for TTest samplesRemainInCor use this samples for correlation plot samplesNotInCor dont use this samples for No NULL NULL	legend	show legend	No	TRUE				
pairsCol color of the pairs use this samples for TTest dont use this samples for TTest samplesRemainInCor use this samples for correlation plot samplesNotInCor dont use this samples for No NULL	ttestPerform	perform a TTest	No	FALSE				
samplesRemainInTTest use this samples for TTest dont use this samples for TTest samplesRemainInCor use this samples for correlation plot samplesNotInCor dont use this samples for No NULL	ttestCol	color of the TTest	No	NULL				
samplesNotInTTest dont use this samples for TTest samplesRemainInCor use this samples for correlation plot samplesNotInCor dont use this samples for No NULL	pairsCol	color of the pairs	No	NULL				
samplesRemainInCor use this samples for correlation plot dont use this samples for No NULL samplesNotInCor dont use this samples for No NULL	samplesRemainInTTest	use this samples for TTest	No	NULL				
samplesRemainInCor use this samples for correlation plot samplesNotInCor dont use this samples for No NULL	samplesNotInTTest	dont use this samples for	No	NULL				
lation plot samplesNotInCor dont use this samples for No NULL		TTest						
lation plot samplesNotInCor dont use this samples for No NULL	samplesRemainInCor	use this samples for corre-	No	NULL				
samplesNotInCor dont use this samples for No NULL		_						
	samplesNotInCor	1 -	No	NULL				
		correlation plot						

# 4 Setting parameters

## 4.1 Input and output directories

First you have to specify the directory where your data is stored (loadPath) and the directory where the results are supposed to be stored (savePath). You can specify the path relatively to the directory you are in when starting the script in R.

```
> params <- list(loadPath = system.file("extdata", package="ddCt"),
+ savePath = getwd())</pre>
```

#### 4.2 Input files

Then you have to specify the name of the exported file from the RT-PCR device containing the individual Ct values (inputFile). Optionally, you may specify a file containing additional annotation data related to the samples (sampleAnnotationFile). Refer to the file 'SampleAnnotation trelated along the package as an example. The file should be a tab-delimited with several columns and one row for each sample. The first column must have the column name of 'Sample' and include the sample names. Also a valid annotation file is in the form of Table 2:

Sample	Variable1	Variable2	Variable3
Sample2	1	0	1
Sample1	1	1	1

Table 2: Simple structure of an annotation file

This file with additional information is important if you want to perform the t-test, if you use special colors for groups of samples, or if you want the samples to be grouped according to one column of the sample annotation file that can be specified through the parameter columnForGrouping. If none of these is important, just set all parameters to NULL. This also holds true for other parameters described below: setting NULL will make the ddCt neglect the parameter. If you have a sample annotation file and want to include only samples from this file into your final object, you can set the parameter 'useOnlySamplesFromAnno'.

```
> params$confFile <- system.file("scripts", "ddCt.conf", package="ddCt")
> params$inputFile <- c("Experiment1.txt")
> params$sampleAnnotationFile = NULL
> params$columnForGrouping = NULL
> params$onlyFromAnnotation = FALSE
```

### 4.3 Change gene names and sample names (optional)

You may change gene names and sample names, for instance to make the figure report of the data pretty. These can be done in the ddCt package.

In the following example, 'Gene1' will become 'Gene A', 'Gene4' will become 'Gene B', and so on. In the final plot, 'Gene A' will be plotted first, then 'Gene B' followed by 'Gene C', 'HK1' and finally 'HK2', which is the result of default ordering or factor. If you want to rename genes, all genes have to be included into this list, otherwise an error will be raised.

The rename of the samples follow a similar way by setting the variable sampleAlias in the params list.

```
> # "Gene5"="Gene C",
> # "Gene2"="HK1",
> # "Gene3"="HK2")
> #params$sampleAlias = NULL
```

**Attention:** If you rename genes or samples, the new names must be used for the parameters in the settings below.

In the commald-line you can set sutch a parameter in this way:

```
Rscript ddCt.R -geneAlias=Gene1="Gene A",Gene2="Gene B" ...
```

#### 4.4 Housekeeping genes and reference samples

To calculate relative expression, one has to specify the housekeeping gene(s) and the reference sample(s). In this sample, one reference sample and two housekeeping genes are used. If more than one object is specified, the names have to be given as shown in this example.

If you specify more than one housekeeping gene, the software will use the mean of the  $C_t$  values of the housekeeping genes for the normalization. If you use more than one sample, the algorithm uses the mean of the chosen samples as the reference.

```
> params$referenceGene <- c("Gene1", "Gene2")
> params$referenceSample <- c("Sample1")</pre>
```

The comand-line version would look like this:

```
Rscript ddCt.R -referenceGene=Gene1,Gene2 ...
```

You may set a threshold to filter the  $C_T$  values, which is used to set an upper boundary of the  $C_T$  value to be considered. Any  $C_T$  value above this threshold will be treated as 'undetermined'.

```
> params$threshold <- 40
```

And next step you have to specify if you want to use the 'mean' or the 'median' to summarize the individual  $C_T$  values for a gene/sample combination. The median is often more appropriate when replicate number is large since it is more robust.

```
> params$mode = "median"
```

#### 4.5 Efficiencies (optional)

You may include amplification efficiency for each gene (see the vignette rtPCR for more background information). There is also the possibility to include error estimates for the efficiencies (for example the standard deviation). These estimates will be used for the error calculation. These efficiencies can be specified in the following way:

```
> params$efficiencies = c("Gene1"=1.9, "Gene2"=2, ...)
```

where "Gene1" represents the gene and 1.9 the efficience which is used for this gene If you use efficiencies, only the raw Ct value and the final 'level' is calculated. There are no 'dCt' or 'ddCt' values. Hence no t-test can be performed if efficiencies are used. In this example we do not use efficiencies.

```
> #params$efficiencies = c("Gene\ A"=1.9,"Gene\ B"=1.8,"HK1"=2,"Gene\ C"=2,"HK2"=2)
> #params$efficienciesError = c("Gene\ A"=0.01,"Gene\ B"=0.1,"HK1"=0.05,"Gene\ C"=0.01,"HK2"=0.
```

# 4.6 Plot parameter (optional)

The following parameters are used to change the final plot. First you have to specify what you want to plot. Here you can specify either or both of the following two choices:

- level For each gene and sample the relative expression to the reference line is plotted
- Ct the raw, unnormalized but merged Ct values are plotted

```
> params$plotMode = c("level","Ct")
```

There may be cases where experiments want to exlcude certain gene and/or sample in the plot. Or sometimes one merely wants to plot a small fraction of genes and/or samples. One could set options in the parameter. For example the following example excludes the "NTC" sample in the plot:

```
> #REMAIN
> params$genesRemainInOutput = NULL
> params$samplesRemainInOutput = NULL
> #NOT
> params$genesNotInOutput = NULL
> params$samplesNotInOutput = NULL
```

Do you want the final plot to be drawn in a way such that the samples are plotted next to each other?

```
> params$groupingBySamples = FALSE
```

set the variable plotPerObject to FALSE if one does not need one plot for each pair.

```
> params$plotPerObject = TRUE
```

This is often useful if you have many genes or samples. Depending on the parameter groupingBySamples either each gene (groupingBySamples=TRUE) will have its own plot, or each sample (groupingBySamples=FALSE) will have its own plot.

If you have a single plot for each individual gene, then you may color the sample bars according to one parameter of your sample annotation file (if you have specified such a file in the beginning of this script). You may also specify the colors (maybe with the help of RColorBrewer).

- > #params\$groupingForPlot = NULL
- > #params\$groupingColor = c("#E41A1C","#377EB8","#4DAF4A","#984EA3","#FF7F00")

You may specify a cutoff for the y axis for all plots. Then all plots have the same scale.

#### > #params\$cutoff = NULL

With the parameter brewerColor you can specify color sets you want to use in order to color the individual bars. For additional information have a look at the description of the RColorBrewer package.

```
> #params$brewerColor = c("Set3", "Set1", "Accent", "Dark2", "Spectral", "PuOr", "BrBG")
```

Set legende as TRUE in case you want the legend along the plot. Sometimes the plot will be messed up if the legend is plotted, so please have a try.

> params\$legende = TRUE

## 4.7 t-test and Wilcoxon test specification (optional)

You may perform either t-test or Wilcoxon test if you have specified a sample annotation file above. **Attention:** The t-test and Wilcoxon test will not work if efficiencies are used for the calculation.

#### > params\$ttestPerform = FALSE

If you want to do a t-test/Wilcoxon test, you have to specify the name of the column of your sample information file that includes the group information needed for the tests. If there are more than two groups in this column, tests for each possible pair of groups are performed.

#### > #params\$ttestCol = NULL

If you want to perform a paired t-test/Wilcoxon test, you have to specify a column of your sample information file that contains information describing the pairing of the samples. A pair of samples must have the same number/parameter in that column. Please have a look at the example.

#### > #params\$pairsCol = NULL

You can specify whether you want to exclude some samples from the t-test. Here we again want to exclude the 'NTC' sample.

- > #params\$samplesRemainInTTest = NULL
  > #params\$samplesNotInTTest = NULL

### 4.8 Housekeeping gene correlation (optional)

If more than one housekeeping gene is used, the correlation (*Pearson* or *Spearman*) for each pair is calculated, together with a correlation plot. You may specify some samples that are not supposed to be used for the correlation calculation, for example negative control. In the default setting those samples are excluded that are also excluded from the plot.

```
> #params$samplesRemainInCor = NULL
> #params$samplesNotInCor = NULL
```

#### 4.9 Execution

Now all parameters have been set and you can pass them to the ddCtExec function

> ddCtExec(params)

At last all you have to do is to start R, change the directory if you like, and then type the following:

> Sweave("rtPCR-usage.Rnw")

The output is stored in a directory called Result\_YOU\_INPUT\_FILENAMES

## 5 Results

Various files are created during the calculation process. The most important files are the following:

- A tab-delimited text file containing all calculated values for each gene/sample combination, e.g. Ct+error, dCt+error, ddCt+error, level+error. Name of the file: 'allValues' followed by the name of the file containing the individual Ct values and extension.
- A .pdf file including the plots, either of the 'levels' or the raw Ct values (depending on your choice) **Name of the file:** 'Level' or 'Ct' followed by 'Result', the name of the file containing the individual Ct values and the extension. In this example the name is 'LevelResultTest.pdf'.
- A tab-delimited text file containing the level and the error of the level. This file can be used in Excel to create own plots with error bars. Name of the file: 'LevelPlusError' followed by the name of the file containing the individual Ct values and extension.
- A tab-delimited text file containing the results from the t-test and the Wilcoxon test (if these tests have been performed). Name of the file: 'ttest' followed by the name of the groups used in the test, followed by the name of the file containing the individual Ct values and extension. In this example we have: 'ttestG1G2Test.txt'.
- An R object containing all calculated values for each gene/sample combination, e.g. Ct+error, dCt+error, ddCt+error, level+error. Furthermore the object includes additional sample information. Name of the file: 'Result' followed by the name of the file containing the individual Ct values and extension.

There are additional tab-delimited text files and .html files containing 'Ct','dCt', 'ddCt' and 'level' information. All this information is also included in the 'allValues' file, but in a different format.

# 6 Session Info

The script has been running in the following session:

- R version 2.15.0 (2012-03-30), x86\_64-unknown-linux-gnu
- Locale: LC\_CTYPE=en\_US.UTF-8, LC\_NUMERIC=C, LC\_TIME=en\_US.UTF-8, LC\_COLLATE=C, LC\_MONETARY=en\_US.UTF-8, LC\_MESSAGES=en\_US.UTF-8, LC\_PAPER=C, LC\_NAME=C, LC\_ADDRESS=C, LC\_TELEPHONE=C, LC\_MEASUREMENT=en\_US.UTF-8, LC\_IDENTIFICATION=C
- Base packages: base, datasets, grDevices, graphics, methods, stats, utils
- Other packages: Biobase 2.16.0, BiocGenerics 0.2.0, RColorBrewer 1.0-5, ddCt 1.10.0, lattice 0.20-6, xtable 1.7-0
- Loaded via a namespace (and not attached): grid 2.15.0, tools 2.15.0