Package 'methylumi'

March 26, 2013

Type Package
Title Handle Illumina methylation data
Version 2.4.0
Date 2012-06-13
Author Sean Davis, Pan Du, Sven Bilke, Tim Triche, Jr., Moiz Bootwalla
Depends Biobase, methods, R (>= 2.13), scales, reshape2, ggplot2
Imports Biobase, graphics, lattice, annotate, genefilter, Annotation Dbi, stats 4, Bioc Generics
Suggests lumi, lattice, limma, xtable, IlluminaHumanMethylation27k.db (>= 1.4.4), IlluminaHumanMethylation450k.db, SQN, MASS,matrixStats, parallel
biocViews DNAMethylation, TwoChannel, Preprocessing, QualityControl,CpGIsland
Maintainer Sean Davis <sdavis2@mail.nih.gov></sdavis2@mail.nih.gov>
Description This package provides classes for holding and manipulating Illumina methylation data. Based on eSet, it can contain MIAME information, sample information, feature information, and multiple matrices of data. An ''intelligent'' import function,methylumiR can read the Illumina text files and create a MethyLumiSet. methylumIDAT can directly read raw IDAT files from HumanMethylation27 and HumanMethylation450 microarrays. Normalization,background correction, and qual ity control features for GoldenGate,Infinium, and Infinium HD arrays are also included.
License GPL-2
R topics documented:
methylumi-package CpGs estimateM extractBarcodeAndPosition featureFilter getAssayDataNameSubstitutions methylData-class MethyLumi-accessors MethyLumi-class

2 methylumi-package

	MethyLumi-strippers	10
	methylumIDAT	11
	methylumiGenerics	12
	MethyLumiM-class	
	MethyLumiQC-class	
	methylumiR	
	MethyLumiSet-class	
	mldat	
	normalizeMethyLumiSet	
	plotSampleIntensities	
	qcplot	
	varFilter	
Index		24
meth	aylumi-package Handle Illumina methylation data	

Description

This package contains a class structure for handling methylation data from Illumina as well as utility functions for loading the data from files generated by Illumina. Normalization that attempts to correct for dye bias is also included.

Important data classes include: MethyLumiSet and MethyLumiQC, both of which are subsets of the MethyLumi class, which is a subset of the eSet class.

A worked example of the use of the package can be found by typing: openVignette().

A full listing of the available documentation can be obtained by typing help.start() and selecting methylumi from the Packages link or by typing library(help="methylumi").

Details

Package: methylumi Type: Package License: GPL

Author(s)

Sean Davis <sdavis2@mail.nih.gov>

References

http://watson.nci.nih.gov/~sdavis/software/R

See Also

Biobase

CpGs 3

CpGs	Data frame describing loci on the 27 and 450k arrays.	

Description

Data frame describing loci on the 27 and 450k arrays.

Usage

```
data(CpGs)
```

Examples

```
data(CpGs)
head(CpGs)
```

estimateM

Estimate methylation M-value matrix

Description

Estimate methylation M-value matrix from MethyLumiM-class object or eSet-class object, which include methylated and unmethylated probe intensities

Usage

```
estimateM(methyLumiM, returnType=c("ExpressionSet", "matrix"), offset=100)
```

Arguments

methyLumiM MethyLumiM-class object or eSet-class object, which include methylated and

unmethylated probe intensities

returnType determine whether return an ExpressionSet (MethyLumiM in this case) or ma-

trix object

offset offset added to the methylated and unmethylated probe intensities when estimat-

ing the M-value

Details

M-value is the log2 ratio between Illumina methylated and unmethylated probe intensities. As variations of small intensities can cause big changes in the ratio estimation, so an offset is added to methylated and unmethylated probe intensities when estimating the M-value.

Please check the lumi package for more details of estimateM function.

Value

A MethyLumiM or matrix object of methylation M-value

Author(s)

Pan DU

References

Du, P., Zhang, X, Huang, C.C., Jafari, N., Kibbe, W.A., Hou, L., and Lin, S.M., (2010) 'Comparison of Beta-value and M-value methods for quantifying methylation levels by microarray analysis', (under review)

extractBarcodeAndPosition

Extract the Barcode and Position Information from Sentrix ID

Description

The sentrix IDs from an illumina sentrix array contain positional information that might be useful. This function simply extracts that information from the ID itself.

Usage

extractBarcodeAndPosition(sentrixids)

Arguments

sentrixids A character vector of sentrix IDs that look like: 1632405013_R001_C001

Value

A data.frame with three columns:

sentrix numeric, the sentrix ID
row numeric, the sentrix row
column numeric, the sentrix column

Author(s)

Sean Davis <sdavis2@mail.nih.gov>

See Also

methylumiR

Examples

 $extractBarcodeAndPosition(c('12341234_R001_C001'))$

featureFilter 5

featureFilter

Annotation-based Filtering of Features (CpG sites) in a MethyLumiSet or MethyLumiM object

Description

Features with insufficient annotation carry little value for the subsequent data analysis. The function featureFilter provides options of filtering features (CpG sites) from a MethyLumiSet (or MethyLumiM) object based on available annotation data.

Usage

```
featureFilter(eset, require.entrez=FALSE, require.GOBP=FALSE, require.GOCC=FALSE, require.GOMF=FALSE, exclude.ChrX=FALSE, require.closeToTSS=FALSE, range.DistToTSS=c(-500, 300), require.CpGisland=FALSE, ...)
```

Arguments

eset A MethyLumiSet or MethyLumiM object.

require.GOBP, require.GOCC, require.GOMF

If TRUE, filter out features whose target genes are not annotated to at least one GO term in BP, CC and MF ontology, respectively.

exclude.ChrX If TRUE, filter out features in chromosome X to avoid gender effect.

require.closeToTSS

If TRUE, filter out features that are not close to transcription start site (TSS). Features without annotation of distance to TSS will also be removed. Can only used for GoldenGate platform.

range.DistToTSS

Ignored if require.colseToTSS is FALSE. A vector of numeric values of length 2, indicating the range of tolerable distance from transcription start site (TSS) in basepair (bp). If require.clostToTSS is TRUE, features whose distance to TSS falls outside this designated range will be removed. The default value is c(-500,300), where -500 represents the distance to TSS from the left and 300 the distance from the right.

require.CpGisland

If TRUE, filter out features that are not in CpG islands.

Unused, but available for specializing methods.

Value

The function featureFilter returns a list consisting of:

eset The filtered MethyLumiSet or MethyLumiM object.

filter.log A list giving details of how many probe sets where removed for each annotation-

based filtering step performed.

Author(s)

Chao-Jen Wong < cwon2@fhcrc.org>

References

R. Bourgon, R. Gentleman, W. Huber, *Independent filtering increases power for detecting differentially expressed genes*, PNAS, vol. 107, no. 21, pp:9546-9551.

See Also

nsFilter

getAssayDataNameSubstitutions

Return a data.frame of AssayData name substitutions.

Description

The Illumina methylation platforms use two distinct platforms, the "goldengate" platform and the "infinium" platform. Each of these uses different file formats as well as different assay techologies. To make the downstream data handling more straightforward and uniform between the two different systems, a simple mapping from the column names in the output files from the Illumina software is used to convert things from Red/Green or Cy5/Cy3 to unmethylated/methylated. This function simply returns that mapping.

Usage

getAssayDataNameSubstitutions()

Details

A file in the extdata directory called "substitutions.txt" contains two columns. The function loads this file and uses the first column as a match against column names in the data file (with the "sample part" removed). If matched, the second column gives the replacement.

Value

A data.frame with two columns, regex and replacement.

Author(s)

Sean Davis <seandavi@gmail.com>

Examples

getAssayDataNameSubstitutions()

methylData-class 7

methylData-class

Class "methylData", superclass for MethyLumiSet and MethyLumiM

Description

A superclass (virtual) for MethyLumiSet and MethyLumiM.

Objects from the Class

A virtual Class: No objects may be created from it.

Methods

```
diagnostics signature(x = "methylData"): diagnostic plots of data
```

methylated.N signature(object = "methylData"): accessor for assayData element of the same name

 $\label{eq:methylated.N} \textbf{methylated.N} \textbf{-} \ \operatorname{signature}(\operatorname{object} = "methylData", \ value = "matrix") \textbf{:} \ replace \ method \ for \ assayData \ element \ of \ the \ same \ name$

```
plotNAs \ \operatorname{signature}(\operatorname{object} = "methylData") \colon ...
```

 $\label{eq:pval.detect} \textbf{pval.detect} \ \ \text{signature} (\textbf{object} = "methylData") \textbf{:} \ \ \textbf{accessor} \ \ \textbf{for} \ \ \textbf{assayData} \ \ \textbf{element} \ \ \textbf{of} \ \ \textbf{the} \ \ \textbf{same} \\ \textbf{name}$

pval.detect<- signature(object = "methylData", value = "numeric"): replace method for assayData element of the same name

unmethylated.N signature(object = "methylData"): accessor for assayData element of the
same name

unmethylated.N<- signature(object = "methylData", value = "matrix"): replace method for assayData element of the same name

Author(s)

Tim Triche, Jr.

See Also

MethyLumiSet,MethyLumiM

Examples

```
showClass("methylData")
```

8 MethyLumi-class

MethyLumi-accessors

methylumi accessors

Description

These functions serve as getters and setters for information in methylumi classes.

Usage

```
betas(object)
pvals(object)
methylated(object)
unmethylated(object)
getHistory(object)
QCdata(object)
```

Arguments

object

an object of class MethyLumi or a subclass

Details

See the methods definitions in MethyLumiSet and MethyLumiQC for details.

Author(s)

Sean Davis <sdavis2@mail.nih.gov>

See Also

normalize Methy Lumi Set, Methy Lumi QC, e Set

MethyLumi-class

The base class for storing Illumina Methylation data

Description

This class inherits from eSet from the Biobase package and is used as a base class for the other two methylumi classes, MethyLumiSet and MethyLumiQC.

Objects from the Class

The MethyLumi class is a virtual class and is not meant to be instantiated. Instead, one should instantiate a MethyLumiSet or a MethyLumiQC object.

MethyLumi-class 9

Slots

```
assayData: Object of class "AssayData"

phenoData: Object of class "AnnotatedDataFrame"

featureData: Object of class "AnnotatedDataFrame" that will hold the annotation columns from
the Beadstudio output, if they are available.

experimentData: Object of class "MIAME"

apportation: Object of class "keyperters"; note that this slot is not currently used but may be used
```

annotation: Object of class "character"; note that this slot is not currently used, but may be used in the future to store the character name of the annotation package, if available.

```
.__classVersion__: Object of class "Versions"
```

Extends

Class "eSet", directly. Class "VersionedBiobase", by class "eSet", distance 2. Class "Versioned", by class "eSet", distance 3.

Methods

 $\label{eq:pvals} \textbf{pvals} \mbox{--} \ \ \mbox{signature} (\mbox{object} = \mbox{"MethyLumi"}, \ \mbox{value} = \mbox{"matrix"}) \mbox{: Set the assayData slot of the same name and stores the P-values from BeadStudio}$

pvals signature(object = "MethyLumi"): Get the assayData slot of the same name

betas<- signature(object = "MethyLumi", value = "matrix"): Set the assayData slot of the same name and represents the methylation values for the samples, analogous to exprs() in gene expression data.

betas signature(object = "MethyLumi"): Get the assayData slot of the same name

methylated<- signature(object = "MethyLumi", value = "matrix"): Set the assayData slot that represents the Methylated single-channel signal

 $\label{eq:methylated} \begin{tabular}{ll} methylated & signature (object = "MethyLumi"): Get the assay Data slot that represents the Methylated single-channel signal \\ \end{tabular}$

unmethylated<- signature(object = "MethyLumi", value = "matrix"): Set the assayData slot
that represents the Unmethylated single-channel signal</pre>

unmethylated signature(object = "MethyLumi"): Get the assayData slot that represents the Unmethylated single-channel signal

controlTypes signature(object = "MethyLumi": Find the unique control type beeds in the QC-data slot.

qcplot signature(object = "MethyLumi",what,...): Plot of QC data. This plot can be useful for diagnosing the problems associated with specific samples or arrays. The value for "what" is one of the control types (which can be found by using controlTypes() on the object.

summary signature(object = "MethyLumi",...): summary method for MethyLumi objects.

Author(s)

Sean Davis <sdavis2@mail.nih.gov>

See Also

```
methylumiR, MethyLumiSet, MethyLumiQC, eSet
```

Examples

```
\#\# The class structure
showClass("MethyLumi")
\#\# read in some data
\#\# Read in sample information
samps <- read.table(system.file("extdata/samples.txt",
                       package = "methylumi"),sep="\t",header=TRUE)
\#\# Perform the actual data reading
\#\# This is an example of reading data from a
## Sentrix Array format file (actually two files,
## one for data and one for QC probes)
mldat <- methylumiR(system.file('extdata/exampledata.samples.txt',
     package='methylumi'),
    qcfile=system.file('extdata/exampledata.controls.txt',
      package="methylumi"),
    sampleDescriptions=samps)
mldat
## Get history information
getHistory(mldat)
## Get QC data, which is another eSet-derived object
QCdata(mldat)
```

MethyLumi-strippers

Strip excessive probe-level data from MethyLumiSets

Description

450k datasets with probe-level stderrs, out-of-band intensities, and bead numbers can become huge. These functions help to manage their growth in memory, at least until preprocessing and QC is completed, whereupon the summary data can be exported to a RangedData-based object of some sort for integration.

Usage

```
stripMethyLumiSet(object)
stripBeadNs(object)
stripBeadSDs(object)
stripOOB(object)
```

Arguments

object

an object of class MethyLumi or a subclass

Author(s)

Tim Triche, Jr. <tim.triche@gmail.com>

methylumIDAT 11

methylumIDAT	methylumIDAT
--------------	--------------

Description

Read a directory of methylumi idat files and return a MethylumiSet.

Usage

```
methylumIDAT(barcodes = NULL, pdat = NULL, parallel = F, n = F, n.sd = F, oob = T, idatPath=getwd(), ...)
```

Arguments

barcodes A vector of barcodes to read. Either this argument or pdat must be specified.

pdat A data.frame describing the samples. A special column named "barcodes" can be used to specify the barcodes to be read.

parallel If TRUE, an attempt will be made to process using multiple cores on a multicore machine.

n

n.sd
oob
idatPath The path to the directory containing the idat files.
...

Details

Read a set of .idat files and return a MethylumiSet object.

Value

A MethylumiSet object.

Author(s)

Tim Triche, Jr.

Examples

12 MethyLumiM-class

methylumiGenerics

Generics defined in methylumi

Description

See the individual classes for details of methods.

Author(s)

Sean Davis, Pan Du, and Tim Triche, Jr.

MethyLumiM-class

Class "MethyLumiM": for Illumina Methylation microarray data using logRatios

Description

MethyLumiM is a class inherited from ExpressionSet-class. It is designed for Illumina Methylation microarray data. The exprs dataMatrix included in the assayData slot of MethyLumiM object includes a matrix of M-values, which is the log2 ratio of methylated and unmethylated probe intensities. The MethyLumiM class include a boxplot function uniquely designed for two-mode histogram data. It also include a coerce function to map from MethyLumi-class, MethyLumiSet-class or other eSet-class inherited object to MethyLumiM class object.

Objects from the Class

Objects can be created by calls of the form new("MethyLumiM", exprs, methylated, unmethylated, detection, meth The "exprs" is a matrix of M-values, which is the log2 ratio of methylated and unmethylated probe intensities; "methylated" and "unmethylated" are intensity matrix measured by methylated and unmethylated probes of Illumina Infinium methylation microarray; "detection" is the detection p-value outputted by Illumina GenomeStudio software; "methylated.N" and "unmethylated.N" are bead numbers for methylated and unmethylated probes. "exprs", "methylated" and "unmethylated" information are required for MethyLumiM class. When creating a new MethyLumiM object, the information of "exprs", "methylated", "unmethylated" and "detection" can also be provided directly through "assayData".

Slots

MethyLumiM-class 13

Extends

Class "ExpressionSet", directly. Class "eSet", by class "ExpressionSet", distance 2. Class "VersionedBiobase", by class "ExpressionSet", distance 3. Class "Versioned", by class "ExpressionSet", distance 4.

Methods

boxplot signature(x = "MethyLumiM"): plot distribution of M-value

coerce signature(from = "eSet", to = "MethyLumiM"): map from MethyLumi-class, MethyLumiSet-class or other eSet-class inherited object to MethyLumiM class object. MethyLumiM object will only keep "exprs", "methylated", "unmethylated" and "detection" data matrix in the assayData.

 $\label{eq:getHistory} \textbf{getHistory} \ \ \textbf{signature} (\textbf{object} = \texttt{"MethyLumiM"}) \textbf{:} \ \textbf{access the operation history of MethyLumiM object}.$

initialize signature(.Object = "MethyLumiM"): class initialization

methylated signature(object = "MethyLumiM"): retrieve the data matrix measured by methylated probes

methylated<- signature(object = "MethyLumiM"): set the data matrix measured by methylated probes

 $\begin{tabular}{ll} \textbf{unmethylated} & signature (object = "MethyLumiM"): retrieve the data matrix measured by unmethylated probes \\ \end{tabular}$

unmethylated<- signature(object = "MethyLumiM"): set the data matrix measured by unmethylated probes

methylated.N signature(object = "MethyLumiM"): retrieve the data matrix keeping the number of beads of methylated probes

methylated.N<- signature(object = "MethyLumiM"): set the data matrix keeping the number of beads of methylated probes

unmethylated.N signature(object = "MethyLumiM"): retrieve the data matrix keeping the number of beads of unmethylated probes

unmethylated.N<- signature(object = "MethyLumiM"): set the data matrix keeping the number of beads of unmethylated probes

detection signature(object = "MethyLumiM"): retrieve detection data matrix in AssayData-class

detection<- signature(object = "MethyLumiM"): set detection data matrix in AssayData-class

 $\textbf{controlData} \ \ \text{signature} (object = "MethyLumiM") : \textbf{retrieve the } controlData \ \textbf{in } MethyLumiQC\text{-}class$

controlData <- signature(object = "MethyLumiM"): set controlData in MethyLumiQC-class

 $\label{eq:dataType} {\bf dataType} \ {\bf signature} ({\bf object} = "MethyLumiM") {\bf :} \ {\bf retrieve} \ {\bf the} \ {\bf dataType}, \ {\bf by} \ {\bf default} \ {\bf it} \ {\bf is} \ "M", \ {\bf it} \ {\bf can} \ {\bf also} \ {\bf be} \ "Beta" \ {\bf or} \ "Intensity"$

dataType<- signature(object = "MethyLumiM"): set dataType in MethyLumiM-class, the value can be "M", "Beta" or "Intensity"

Author(s)

Pan DU

References

1. Du, P., Zhang, X, Huang, C.C., Jafari, N., Kibbe, W.A., Hou, L., and Lin, S.M., (2010) 'Comparison of Beta-value and M-value methods for quantifying methylation levels by microarray analysis'

See Also

MethyLumi-class and MethyLumiSet-class

Examples

```
showClass("MethyLumiM")
```

MethyLumiQC-class

Class "MethyLumiQC" for holding Illumina methylation QC data

Description

This class inherits from the MethyLumi class (and therefore, from eSet in Biobase) and is designed to hold QC data from Illumina Beadstudio output. These data can be potentially useful when determining the cause for quality problems.

Objects from the Class

 $Objects\ can\ be\ created\ by\ calls\ of\ the\ form\ new ("MethyLumiQC", assayData, phenoData, featureData, experiment) and the property of the form\ new ("MethyLumiQC", assayData, phenoData, featureData, experiment) and the property of the form\ new ("MethyLumiQC", assayData, phenoData, featureData, experiment) and the property of the form\ new ("MethyLumiQC", assayData, phenoData, featureData, experiment) and the property of the form\ new ("MethyLumiQC", assayData, phenoData, featureData, experiment). The property of the form\ new ("MethyLumiQC", assayData, phenoData, featureData, experiment) and the property of the form\ new ("MethyLumiQC", assayData, phenoData, featureData, experiment). The property of the property of$

Slots

```
assayData: Object of class "AssayData"

phenoData: Object of class "AnnotatedDataFrame"

featureData: Object of class "AnnotatedDataFrame" containing the annotation columns from
        the Illumina Beadstudio output. In particular, the names of the probes describe the types of
        control probes.

experimentData: Object of class "MIAME"

annotation: Object of class "character", not currently used

. classVersion : Object of class "Versions"
```

Extends

```
Class "MethyLumi", directly. Class "eSet", by class "MethyLumi", distance 2. Class "VersionedBiobase", by class "MethyLumi", distance 3. Class "Versioned", by class "MethyLumi", distance 4.
```

Methods

```
initialize signature(.Object = "MethyLumiQC")
Cy3.N signature(object = "MethyLumiQC"): ...
Cy3<- signature(object = "MethyLumiQC", value = "matrix"): ...
Cy5.N signature(object = "MethyLumiQC"): ...
Cy5<- signature(object = "MethyLumiQC", value = "matrix"): ...
QCdata<- signature(object = "MethyLumiSet", value = "MethyLumiQC"): ...
combine signature(x = "MethyLumiQC", y = "MethyLumiQC"): ...
controlData<- signature(object = "MethyLumiSet", value = "MethyLumiQC"): ...
controlTypes signature(object = "MethyLumiQC"): determine the character vector of control types from the QCdata information</pre>
```

methylumiR 15

```
hist signature(x = "MethyLumiQC"): ...
intensitiesByChannel signature(object = "MethyLumiQC"): ...
methylated signature(object = "MethyLumiQC"): ...
negctls.stderr signature(object = "MethyLumiQC", channel = "character"): ...
negctls.stderr signature(object = "MethyLumiQC", channel = "missing"): ...
negctls signature(object = "MethyLumiQC", channel = "character"): ...
negctls signature(object = "MethyLumiQC", channel = "missing"): ...
negnorm signature(object = "MethyLumiQC", channel = "character"): ...
negnorm signature(object = "MethyLumiQC", channel = "missing"): ...
normctls signature(object = "MethyLumiQC"): ...
qcplot signature(object = "MethyLumiQC", what, ...): QC plots of various controltypes
unmethylated signature(object = "MethyLumiQC"): ...
```

Author(s)

Sean Davis <sdavis2@mail.nih.gov>

See Also

methylumiR, MethyLumiSet, MethyLumi, eSet

Examples

showClass("MethyLumiQC")

methylumiR

Load data from Illumina methylation platform

Description

This function is useful for loading Illumina methylation data into a MethyLumiSet object. Sample information can be supplied and will then be incorporated into the resulting phenoData slot.

Usage

```
methylumiR(filename, qcfile=NULL, sampleDescriptions = NULL, sep = NULL, ...)
```

Arguments

filename A filename of the excel-like file from BeadStudio acfile A filename of the excel-like file from BeadStudio

sampleDescriptions

A data frame that contains at least one column, SampleID (case insensitive). This column MUST match the part of the column headers before the .Avg Beta, etc. Also, if a column called SampleLabel (case insensitive), it is used for sam-

ple labels, IF the sampleLabel column contains unique identifiers

seperator used in the BeadStudio (or GenomeStudio) output file. If it is NULL, sep

the function will automatically estimate it.

Passed into read.delim()

16 MethyLumiSet-class

Details

Used to construct a MethyLumiSet object....

Value

A MethyLumiSet object

Author(s)

Sean Davis <sdavis2@mail.nih.gov>

See Also

MethyLumiSet-class, MethyLumiQC-class

Examples

MethyLumiSet-class

Class "MethyLumiSet" for containing Illumina methylation data

Description

This class inherits from the MethyLumi class (and therefore, from eSet in Biobase) and is designed to hold both the intensities and the calculated betas, as well as pvalues if present.

Objects from the Class

Objects can be created by calls of the form new("MethyLumiSet", assayData, phenoData, featureData, experiment! An object of this type is the main storage class for methylation data from Illumina. Subsetting, etc., works as normal (rows represent genes, columns represent samples). There is also a rudimentary history tracking system that is modeled after that from the lumi package.

Slots

```
QC: Object of class "QCDataOrNULL", containing either the MethyLumiQC object or NULL history: Object of class "data.frame", containing a running history of transforms to the data contained herein assayData: Object of class AssayData phenoData: Object of class AnnotatedDataFrame
```

MethyLumiSet-class 17

```
featureData: Object of class AnnotatedDataFrame, containing the annotation columns from the Illumina Beadstudio output experimentData: Object of class MIAME annotation: Object of class "character", not currently used .__classVersion__: Object of class "Versions" protocolData: Object of class "AnnotatedDataFrame" that contains protocol information, including scan date if available
```

Extends

Class "MethyLumi", directly. Class "methylData", directly. Class "eSet", by class "MethyLumi", distance 2. Class "VersionedBiobase", by class "MethyLumi", distance 3. Class "Versioned", by class "MethyLumi", distance 4.

Methods

```
[ signature(x = "MethyLumiSet"): subsetting, genes as rows, samples as columns
betas<- signature(object = "MethyLumiSet", value = "matrix"): Set the assayData slot of
the same name
betas signature(object = "MethyLumiSet"): Get the assayData slot of the same name
boxplot signature(x = "MethyLumiSet"): boxplot of all sample betas
combine signature(x = "MethyLumiSet", y = "MethyLumiSet")
corplot signature(x = "MethyLumiSet"): returns m-values
getHistory signature(object = "MethyLumiSet"): returns a data.frame containing the history
for this object
hist signature(x = "MethyLumiSet"): histogram of the betas for the data
initialize signature(.Object = "MethyLumiSet")
pairs signature(x = "MethyLumiSet"): pairs plot of the betas for the object. Note that pairs
plots of more than a few samples are not very useful.</pre>
```

plotSampleIntensities signature(x = "MethyLumiSet"): The intensities as output by the Beadstudio software often show a considerable amount of dye bias. This method shows a graphical example of this dye bias. In short, for each of the Cy3 and Cy5 channels, a cutoff in beta is used to calculate which Cy3 and Cy5 values should be plotted at high-methylation and low-methylation status. Any offset between Cy3 and Cy5 when plotted in this way likely represents dye bias and will lead to biases in the estimate of beta.

 $\label{eq:QCdata} \mbox{\bf QCdata$<-} \ \ \mbox{signature} (\mbox{object} = \mbox{\tt "MethyLumiSet"}, \ \mbox{value} = \mbox{\tt "MethyLumiQC"}) : \mbox{assign QC data} \ \mbox{to the QC slot}$

```
\label{eq:QCdata} QCdata \ \ \mathrm{signature} (\mathrm{object} = "MethyLumiSet") \text{: retrieve the QC data}.
```

show signature(object = "MethyLumiSet")

methylated<- signature(object = "MethyLumiSet", value = "matrix"): Set the assayData slot associated with methylated intensity

 $\label{eq:methylated} \textbf{methylated} \ \ \textbf{signature} (object = "MethyLumiSet") : \textbf{Get the assayData slot associated with methylated intensity}$

unmethylated<- signature(object = "MethyLumiSet", value = "matrix"): Set the assayData
slot associated with unmethylated intensity</pre>

18 MethyLumiSet-class

```
unmethylated signature(object = "MethyLumiSet"): Get the assayData slot associated with
        unmethylated intensity
    qcplot signature(object = "MethyLumiSet", what, ...): QC plots of various controltypes
    controlTypes signature(object = "MethyLumiSet"): determine the character vector of control
        types from the QCdata information
    Cy3.N signature(object = "MethyLumiSet"): ...
    Cy5.N signature(object = "MethyLumiSet"): ...
    combine27k450k signature(x = "MethyLumiSet", y = "MethyLumiSet"): ...
    controlData signature(object = "MethyLumiSet"): ...
    controlData<- signature(object = "MethyLumiSet", value = "MethyLumiQC"): ...
    featureFilter signature(eset = "MethyLumiSet"): ...
    intensities.IB signature(x = "MethyLumiSet", channel = "character"): ...
   intensities.IB signature(x = "MethyLumiSet", channel = "missing"): ...
    intensities.M signature(x = "MethyLumiSet", channel = "character"): ...
   intensities.M signature(x = "MethyLumiSet", channel = "missing"): ...
    intensities.OOB.allelic signature(x = "MethyLumiSet", channel = "character", allele = "character"):
   intensities.OOB.allelic signature(x = "MethyLumiSet", channel = "missing", allele = "missing"):
   intensities.OOB signature(x = "MethyLumiSet", channel = "character"): ...
   intensities.OOB signature(x = "MethyLumiSet", channel = "missing"): ...
   intensities.U signature(x = "MethyLumiSet", channel = "character"): ...
   intensities.U signature(x = "MethyLumiSet", channel = "missing"): ...
    intensitiesByChannel signature(object = "MethyLumiSet"): ...
    negctls.stderr signature(object = "MethyLumiSet", channel = "character"): ...
    negctls.stderr signature(object = "MethyLumiSet", channel = "missing"): ...
    negctls signature(object = "MethyLumiSet", channel = "character"): ...
    negctls signature(object = "MethyLumiSet", channel = "missing"): ...
    negnorm signature(object = "MethyLumiSet", channel = "character"): ...
    negnorm signature(object = "MethyLumiSet", channel = "missing"): ...
    normctls signature(object = "MethyLumiSet"): ...
    plotSampleIntensities signature(x = "MethyLumiSet"): ...
    probeNAs signature(object = "MethyLumiSet"): ...
    sampleNAs signature(object = "MethyLumiSet"): ...
    total.intensity signature(object = "MethyLumiSet"): ...
    varFilter signature(eset = "MethyLumiSet"): ...
Author(s)
    Sean Davis & Tim Triche, Jr.
See Also
    methylumiR, normalizeMethyLumiSet, methylumIDAT, MethyLumiQC, eSet
```

Examples

```
showClass("MethyLumiSet")\\
```

mldat 19

mldat

Example SAM format Illumina methylation dataset

Description

This is an example MethyLumiSet object.

Usage

data(mldat)

Examples

data(mldat)

normalizeMethyLumiSet Normalize a MethyLumiSet, accounting for dye bias

Description

The Illumina GoldenGate methylation platform uses two colors, one to represent the unmethylated state and the other to represent the methylated state. This function corrects that dye bias and recalculates the betas based on the corrected intensities.

For HumanMethylation27 data, the function does nothing.

For HumanMethylation450 data, the function delegates to normalizeViaControls() the task of scaling red and green intensities against a reference array (chip) which defaults to the first chip in a set. The code to do this is based on code from the 'minfi' package and uses the built-in normalization controls to scale the channels of the samples, so that a consistent degree of dye bias is maintained for Infinium II probes across an experiment or set of experiments.

Usage

```
normalizeMethyLumiSet(x, beta.cuts = c(0.2, 0.8), mapfun = c("atan", "ratio"))
```

Arguments

X	A MethyLumiSet	object
---	----------------	--------

beta.cuts Two numeric values with the first less than the second and between 0 and 1,

representing the beta cutoffs that will be used when determining the median

intensities to which to correct. See details below.

mapfun Either "atan" or "ratio". See details below.

20 plotSampleIntensities

Details

For HumanMethylation450 data, the function delegates to normalizeViaControls() the task of scaling red and green intensities against a reference array (chip) which defaults to the first chip in a set. The code to do this is based on code from the 'minfi' package and uses the built-in normalization controls to scale the channels of the samples, so that a consistent degree of dye bias is maintained for Infinium II probes across an experiment or set of experiments. The remainder of the documentation below is specific to GoldenGate data.

The Illumina GoldenGate methylation platform uses two colors, one to represent the unmethylated state and the other to represent the methylated state. This function corrects that dye bias and recalculates the betas based on the corrected intensities.

As a first step, the medians for each of Cy3 and Cy5 are calculated at high and low betas, representing the (nearly) fully methylated state and the (nearly) fully unmethylated states. Values of Cy3 and Cy5 that are negative are set to zero for this process. Then, the Cy5 medians are adjusted to match those of the Cy3 channel, thereby correcting the dye bias.

To map the new intensities back to betas, one of two map functions can be used. The default is the atan(Cy3/Cy5). The ratio maps using the function (Cy3/Cy3+Cy5). The differences should be very small, but we feel that the atan map function is probably the mathematically appropriate way of doing this.

Value

A new "MethyLumiSet" that contains the corrected betas and the adjusted intensities.

Author(s)

Sean Davis <sdavis2@mail.nih.gov>

Examples

plotSampleIntensities

Plot the sample intensities.

Description

The Illumina methylation platforms all show a significant dye bias. The plotSampleIntensities method shows the density plots for the two channels allowing direct visualization of the effect.

Usage

```
plotSampleIntensities(x,beta.cuts,s)
```

qcplot 21

Arguments

x an object of class MethyLumi or a subclass

beta.cuts cutoffs for low and high beta values

s sample number to plot

Examples

```
data(mldat)
plotSampleIntensities(mldat,s=1)
```

qcplot

Methods for dealing with control data for Illumina methylation data.

Description

The qcplot function simply generates a plot of the control probe information for a given control Type.

Usage

```
\begin{array}{l} \operatorname{qcplot}(\operatorname{object}, \operatorname{controltype}, \ldots) \\ \operatorname{controlTypes}(\operatorname{object}, \ldots) \end{array}
```

Arguments

object An object of class MethyLumiSet or MethyLumiQC

controltype A single character value representing the bead type to plot from the quality con-

trol data. The available types are accessible via the controlTypes method.

... passed to plot function

Details

The descriptions of the various control types can be obtained from the Illumina methylation user's guides.

Author(s)

Sean Davis <sdavis2@mail.nih.gov>

See Also

MethyLumiSet, MethyLumiQC

Examples

```
data(mldat)
controlTypes(mldat)
qcplot(mldat,controlTypes(mldat)[3])
```

22 varFilter

varFilter	Variation-based Filtering of Features (CpG sites) in a MethyLumiSet
	or MethyLumiM object

Description

The function varFilter removes features exhibiting little variation across samples. Such non-specific filtering can be advantageous for downstream data analysis.

Usage

```
varFilter(eset, var.func=IQR, var.cutoff=0.5, filterByQuantile=TRUE, ...)
```

Arguments

eset An MethyLumiSet or MethyLumiM object.

var.func The function used as the per-feature filtering statistics.

var.cutoff A numeric value indicating the cutoff value for variation. If filterByQuantile is

TRUE, features whose value of var.func is less than var.cutoff-quantile of all var.func value will be removed. It FALSE, features whose values are less than

var.cutoff will be removed.

filterByQuantile A logical indicating whether var.cutoff is to be interprested as a quantile of all

var.func (the default), or as an absolute value.

... Unused, but available for specializing methods.

Details

This function is a counterpart of functions nsFilter and varFilter available from the genefilter package. See R. Bourgon et. al. (2010) and nsFilter for detail.

It is proven that non-specific filtering, for which the criteria does not depend on sample class, can increase the number of discoverie. Inappropriate choice of test statistics, however, might have adverse effect. limma's moderated t-statistics, for example, is based on empirical Bayes approach which models the conjugate prior of gene-level variance with an inverse of χ^2 distribution scaled by observed global variance. As the variance-based filtering removes the set of genes with low variance, the scaled inverse χ^2 no longer provides a good fit to the data passing the filter, causing the limma algorithm to produce a posterior degree-of-freedom of infinity (Bourgon 2010). This leads to two consequences: (i) gene-level variance estimate will be ignore, and (ii) the p-value will be overly optimistic (Bourgon 2010).

Value

The function featureFilter returns a list consisting of:

eset The filtered MethyLumiSet or MethyLumiM object.

filter.log Shows many low-variant features are removed.

Author(s)

Chao-Jen Wong < cwon2@fhcrc.org >

varFilter 23

References

R. Bourgon, R. Gentleman, W. Huber, *Independent filtering increases power for detecting differentially expressed genes*, PNAS, vol. 107, no. 21, pp:9546-9551, 2010.

See Also

nsFilter

Examples

```
data(mldat)
## keep top 75 percent
filt <- varFilter(mldat, var.cutoff=0.25)
filt$filter.log
dim(filt$eset)</pre>
```

Index

*Topic IO	combine, MethyLumiQC, MethyLumiQC-method
getAssayDataNameSubstitutions, 6	(MethyLumiQC-class), 14
MethyLumi-accessors, 8	combine, MethyLumiSet, MethyLumiSet-method
methylumi-package, 2	(MethyLumiSet-class), 16
methylumIDAT, 11	combine27k450k (methylumiGenerics), 12
methylumiR, 15	combine27k450k,MethyLumiSet,MethyLumiSet-method
*Topic classes	(MethyLumiSet-class), 16
methylData-class, 7	compare.chips (methylumiGenerics), 12
MethyLumi-class, 8	controlData (methylumiGenerics), 12
MethyLumiM-class, 12	controlData,MethyLumiM-method
MethyLumiQC-class, 14	(MethyLumiM-class), 12
MethyLumiSet-class, 16	controlData,MethyLumiSet-method
*Topic datasets	(MethyLumiSet-class), 16
CpGs, 3	controlData<- (methylumiGenerics), 12
mldat, 19	controlData<-,MethyLumiM,ANY-method
*Topic hplot	(MethyLumiM-class), 12
plotSampleIntensities, 20	controlData<-,MethyLumiSet
qcplot, 21	(MethyLumiSet-class), 16
*Topic manip	controlData<-,MethyLumiSet,MethyLumiQC-method
extractBarcodeAndPosition, 4	(MethyLumiQC-class), 14
normalizeMethyLumiSet, 19	controlTypes, 21
*Topic methods	controlTypes (qcplot), 21
estimateM, 3	controlTypes,MethyLumiM-method
[,MethyLumiM-method	(MethyLumiM-class), 12
(MethyLumiM-class), 12	control Types, Methy Lumi QC-method
[,MethyLumiSet-method	(MethyLumiQC-class), 14
(MethyLumiSet-class), 16	controlTypes,MethyLumiSet-method
AnnotatedDataFrame, 16, 17	(MethyLumiSet-class), 16
AssayData, 16	corplot (MethyLumi-accessors), 8
·	corplot,MethyLumiSet-method
betas (MethyLumi-accessors), 8	(MethyLumiSet-class), 16
betas,MethyLumiM-method	CpGs, 3
(MethyLumiM-class), 12	cy3 (methylumiGenerics), 12
betas,MethyLumiSet-method	Cy3.N (methylumiGenerics), 12
(MethyLumiSet-class), 16	Cy3.N,MethyLumiQC-method
betas<- (MethyLumi-accessors), 8	(MethyLumiQC-class), 14
betas<-,MethyLumiSet,matrix-method	Cy3.N,MethyLumiSet-method
(MethyLumiSet-class), 16	(MethyLumiSet-class), 16
Biobase, 2	Cy3<- (methylumiGenerics), 12
boxplot, Methy Lumi Set-method	Cy3<-,MethyLumiQC,matrix-method
(MethyLumiSet-class), 16	(MethyLumiQC-class), 14
combine,MethyLumiM,MethyLumiM-method	cy5 (methylumiGenerics), 12
(MethyLumiM-class) 12	Cv5 N (methylumiGenerics), 12

INDEX 25

Cy5.N,MethyLumiQC-method	IDATsToDFs (methylumiGenerics), 12
(MethyLumiQC-class), 14	IDATtoDF (methylumiGenerics), 12
Cy5.N,MethyLumiSet-method	intensities.IB (methylumiGenerics), 12
(MethyLumiSet-class), 16	intensities. IB, Methy Lumi Set, character-method
Cy5<- (methylumiGenerics), 12	(MethyLumiSet-class), 16
Cy5<-,MethyLumiQC,matrix-method	intensities.IB,MethyLumiSet,missing-method
(MethyLumiQC-class), 14	(MethyLumiSet-class), 16
	intensities.M (methylumiGenerics), 12
dataType (methylumiGenerics), 12	intensities.M,MethyLumiSet,character-method
dataType,MethyLumiM-method	(MethyLumiSet-class), 16
(MethyLumiM-class), 12	intensities.M,MethyLumiSet,missing-method
dataType<- (methylumiGenerics), 12	(MethyLumiSet-class), 16
data Type <-, Methy Lumi M, ANY-method	intensities.OOB (methylumiGenerics), 12
(MethyLumiM-class), 12	intensities.OOB,MethyLumiSet,character-method
dataType<-,MethyLumiM-method	(MethyLumiSet-class), 16
(MethyLumiM-class), 12	intensities.OOB,MethyLumiSet,missing-method
designIItoMandU (methylumiGenerics), 12	(MethyLumiSet-class), 16
designItoMandU (methylumiGenerics), 12	intensities.OOB.allelic,MethyLumiSet,character,character-met
detection (methylumiGenerics), 12	(MethyLumiSet-class), 16
detection, MethyLumiM-method	7.
(MethyLumiM-class), 12	intensities.OOB.allelic,MethyLumiSet,missing,missing-method
detection<- (methylumiGenerics), 12	(MethyLumiSet-class), 16
detection<-,MethyLumiM-method	intensities.U (methylumiGenerics), 12
(MethyLumiM-class), 12	intensities.U,MethyLumiSet,character-method
diagnostics (methylumiGenerics), 12	(MethyLumiSet-class), 16
diagnostics, methylData-method	intensities.U,MethyLumiSet,missing-method
(methylData-class), 7	(MethyLumiSet-class), 16
· · · · · · · · · · · · · · · · · · ·	intensitiesByChannel (methylumiGenerics),
eSet, 2, 8, 9, 13–15, 17, 18	12
estimateM, 3	intensitiesByChannel,MethyLumiM-method
ExpressionSet, 13	(MethyLumiM-class), 12
exprs,MethyLumiSet-method	intensitiesByChannel,MethyLumiQC-method
(MethyLumiSet-class), 16	(MethyLumiQC-class), 14
extractBarcodeAndPosition, 4	intensitiesByChannel,MethyLumiSet-method
C	(MethyLumiSet-class), 16
featureFilter, 5	
featureFilter,MethyLumiM-method	lumiCSV (methylumiGenerics), 12
(featureFilter), 5	lumIDAT (methylumiGenerics), 12
featureFilter,MethyLumiSet-method	
(featureFilter), 5	mergeProbeDesigns (methylumiGenerics),
getAssayDataNameSubstitutions, 6	12
getBarcodes (methylumiGenerics), 12	methylated (MethyLumi-accessors), 8
getHistory (MethyLumi-accessors), 8	methylated,MethyLumi-method
	(MethyLumi-class), 8
getHistory,MethyLumiM-method (MethyLumiM close) 12	methylated,MethyLumiM-method
(MethyLumiM-class), 12	(MethyLumiM-class), 12
getHistory,MethyLumiSet-method	methylated,MethyLumiQC-method
(MethyLumiSet-class), 16	(MethyLumiQC-class), 14
getProtocolData (methylumiGenerics), 12	methylated,MethyLumiSet-method
hist,MethyLumiQC-method	(MethyLumiSet-class), 16
(MethyLumiQC-class), 14	methylated.N (methylumiGenerics), 12
hist, Methy Lumi Set-method	methylated.N (methylumGenerics), 12 methylated.N,methylData-method
(MethyLumiSet-class) 16	(methylData-class) 7
LIVICILIY DIHIHOCE-CIASSI III	LIUGHIYI ZAJAECIASSI. Z

26 INDEX

metnylated.N,MetnyLumiM-metnod	negctis,MetnyLumiSet,missing-method
(MethyLumiM-class), 12	(MethyLumiSet-class), 16
methylated.N<- (methylumiGenerics), 12	negctls.stderr, MethyLumiQC, character-method
methylated.N<-,methylData,matrix-method	(MethyLumiQC-class), 14
(methylData-class), 7	negctls.stderr,MethyLumiQC,missing-method
methylated.N<-,MethyLumiM,ANY-method	(MethyLumiQC-class), 14
(MethyLumiM-class), 12	negctls.stderr,MethyLumiSet,character-method
methylated<- (MethyLumi-accessors), 8	(MethyLumiSet-class), 16
methylated<-,MethyLumi,matrix-method	negctls.stderr,MethyLumiSet,missing-method
(MethyLumi-class), 8	(MethyLumiSet-class), 16
methylated<-,MethyLumiM,ANY-method	negnorm (methylumiGenerics), 12
(MethyLumiM-class), 12	negnorm,MethyLumiM,character-method
methylated<-,MethyLumiSet,matrix-method	(MethyLumiM-class), 12
(MethyLumiSet-class), 16	negnorm, Methy Lumi M, missing-method
methylData, 17	(MethyLumiM-class), 12
methylData-class, 7	negnorm, Methy Lumi QC, character-method
MethyLumi, 2, 14, 15, 17	(MethyLumiQC-class), 14
methylumi (methylumi-package), 2	negnorm, Methy Lumi QC, missing-method
MethyLumi-accessors, 8	(MethyLumiQC-class), 14
MethyLumi-class, 8	negnorm, Methy Lumi Set, character-method
methylumi-package, 2	(MethyLumiSet-class), 16
MethyLumi-strippers, 10	negnorm, MethyLumiSet, missing-method
methylumi.bgcorr (methylumiGenerics), 12	(MethyLumiSet-class), 16
methylumi.diagnostics	normalizeMethyLumiSet, 8, 18, 19
(methylumiGenerics), 12	normctls (methylumiGenerics), 12
methylumiCSV (methylumiGenerics), 12	normetis, Methy Lumi M-method
The state of the s	(MethyLumiM-class), 12
methylumIDAT, 11, 18	normctls,MethyLumiQC-method
methylumiGenerics, 12	(MethyLumiQC-class), 14
MethyLumiM, 7	normctls, Methy Lumi Set-method
MethyLumiM (MethyLumiM-class), 12	· · · · · · · · · · · · · · · · · · ·
MethyLumiM-class, 12	(MethyLumiSet-class), 16
MethyLumiQC, 2, 8, 9, 16, 18, 21	nsFilter, 6, 22, 23
MethyLumiQC-class, 14	
methylumiR, 9, 15, 15, 18	pairs,MethyLumiSet-method
MethyLumiSet, 2, 7–9, 15, 19, 21	(MethyLumiSet-class), 16
MethyLumiSet-class, 16	plotNAs (methylumiGenerics), 12
MIAME, 17	plotNAs,methylData-method
mldat, 19	(methylData-class), 7
	plotNegOob (methylumiGenerics), 12
NChannelSetToMethyLumiSet	plotSampleIntensities, 20
(methylumiGenerics), 12	plot Sample Intensities, Methy Lumi Set-method
negctls (methylumiGenerics), 12	(MethyLumiSet-class), 16
negctls, Methy Lumi M, character-method	probeNAs (methylumiGenerics), 12
(MethyLumiM-class), 12	probeNAs, MethyLumiM-method
negctls,MethyLumiM,missing-method	(MethyLumiM-class), 12
(MethyLumiM-class), 12	probeNAs, MethyLumiSet-method
negctls, Methy LumiQC, character-method	(MethyLumiSet-class), 16
(MethyLumiQC-class), 14	psummary (methylumiGenerics), 12
negctls, Methy Lumi QC, missing-method	pval.detect (methylumiGenerics), 12
(MethyLumiQC-class), 14	pval.detect,methylData-method
negctls,MethyLumiSet,character-method	(methylData-class), 7
(MethyLumiSet-class), 16	pval.detect<- (methylumiGenerics), 12

INDEX 27

pval.detect < -, methyl Data, numeric-method	total. intensity, Methy Lumi M-method
(methylData-class), 7	(MethyLumiM-class), 12
pvals (MethyLumi-accessors), 8	total. intensity, Methy Lumi Set-method
pvals,MethyLumi-method	(MethyLumiSet-class), 16
(MethyLumi-class), 8	
pvals,MethyLumiM-method	unmethylated (MethyLumi-accessors), 8
(MethyLumiM-class), 12	unmethylated, MethyLumi-method
pvals<- (MethyLumi-accessors), 8	(MethyLumi-class), 8
pvals<-,MethyLumi,matrix-method	unmethylated, MethyLumiM-method
(MethyLumi-class), 8	(MethyLumiM-class), 12
	unmethylated, MethyLumiQC-method
qc.probe.plot (methylumiGenerics), 12	(MethyLumiQC-class), 14
QCdata (MethyLumi-accessors), 8	unmethylated, MethyLumiSet-method
QCdata,MethyLumiM-method	(MethyLumiSet-class), 16
(MethyLumiM-class), 12	unmethylated.N (methylumiGenerics), 12
QCdata,MethyLumiSet-method	unmethylated.N,methylData-method
(MethyLumiSet-class), 16	(methylData-class), 7
QCdata<- (MethyLumi-accessors), 8	unmethylated.N,MethyLumiM-method
QCdata<-,MethyLumiSet	(MethyLumiM-class), 12
(MethyLumiSet-class), 16	unmethylated.N<- (methylumiGenerics),
QCdata<-,MethyLumiSet,MethyLumiQC-method	12
(MethyLumiQC-class), 14	unmethylated.N<-,methylData,matrix-method
qcplot, 21	(methylData-class), 7
qcplot,Methylumi-method	unmethylated.N<-,MethyLumiM,ANY-method
(MethyLumi-class), 8	(MethyLumiM-class), 12
	unmethylated<- (MethyLumi-accessors), 8
qcplot,MethyLumiQC-method	unmethylated<-,MethyLumi,matrix-method
(MethyLumiQC-class), 14	(MethyLumi-class), 8
qcplot,MethyLumiSet-method	unmethylated<-,MethyLumiM,ANY-method
(MethyLumiSet-class), 16	(MethyLumiM-class), 12
IM (1 I IDATE) (1 1 'C' ')	unmethylated<-,MethyLumiSet,matrix-method
readMethyLumIDAT (methylumiGenerics),	(MethyLumiSet-class), 16
12	(MethyLumiset-class), 10
sampleNAs (methylumiGenerics), 12	varFilter, 22
sampleNAs,MethyLumiM-method	varFilter, MethyLumiM-method (varFilter),
(MethyLumiM-class), 12	22
sampleNAs,MethyLumiSet-method	varFilter,MethyLumiSet-method
(MethyLumiSet-class), 16	(varFilter), 22
show, Methy Lumi Set-method	Versioned, 9, 13, 14, 17
(MethyLumiSet-class), 16	VersionedBiobase, 9, 13, 14, 17
stripBeadNs (MethyLumi-strippers), 10	, , , , , , , , , , , , , , , , , , , ,
stripBeadSDs (MethyLumi-strippers), 10	
stripBeadSDs (MethyLumi-strippers), 10 stripMethyLumiSet (MethyLumi-strippers),	
10	
stripOOB (MethyLumi-strippers), 10	
subset27kas450k (methylumiGenerics), 12	
subset450kas27k (methylumiGenerics), 12	
subsetCommonProbes	
(methylumiGenerics), 12	
summary, Methy Lumi-method	
(MethyLumi-class), 8	
total.intensity (methylumiGenerics), 12	
obai.moensity (methyrumidenerics), 12	