Package 'beer'

October 18, 2022

```
Title Bayesian Enrichment Estimation in R
Version 1.0.0
Description BEER implements a Bayesian model for analyzing
     phage-immunoprecipitation sequencing (PhIP-seq) data. Given a PhIPData
     object, BEER returns posterior probabilities of enriched antibody responses,
     point estimates for the relative fold-change in comparison to negative
     control samples, and more. Additionally, BEER provides a convenient
     implementation for using edgeR to identify enriched antibody responses.
License MIT + file LICENSE
Encoding UTF-8
LazyData false
Depends R (>= 4.2.0), PhIPData (>= 1.1.1), rjags
Imports cli, edgeR, BiocParallel, methods, progressr, stats,
     SummarizedExperiment, utils
Suggests testthat (>= 3.0.0), BiocStyle, covr, knitr, rmarkdown,
     dplyr, ggplot2, spelling
SystemRequirements JAGS (4.3.0)
biocViews Software, StatisticalMethod, Bayesian, Sequencing, Coverage
URL https://github.com/athchen/beer/
BugReports https://github.com/athchen/beer/issues
Config/testthat/edition 3
VignetteBuilder knitr
RoxygenNote 7.1.2
Language en-US
git_url https://git.bioconductor.org/packages/beer
git_branch RELEASE_3_15
git_last_commit 4487873
git_last_commit_date 2022-04-26
```

Type Package

2 .beadsRRBeer

Date/Publication 2022-10-18

Author Athena Chen [aut, cre] (https://orcid.org/0000-0001-6900-2264),
Rob Scharpf [aut],
Ingo Ruczinski [aut]

Maintainer Athena Chen <achen70@jhu.edu>

R topics documented:

	ISRRBeer Function to run the heads-only round robin using BEER	_
Index		26
	summarizeRunOne	25
	e	23 24
		23
	C	22
	8	21
	<u> </u>	19
	edgeROne	18
	brewOne	17
		15
		14
	• •	13
	.tidyInputsPrior	13
		12
	· ·	12
	e e	11
		11
	e	10
	getABMOM	9
	getABMLEProp	8
	getABMLE	8
	.getABEdgeR	7
	.edgeRBeads	6
	.checkOverwrite	6
	.brewSamples	4 5
	.beadsRREdgeR	3
	.beadsRRBeer	2

Description

Each sample is run in comparison to all other beads-only samples to approximate the false positive rate of detecting enrichments.

.beadsRREdgeR 3

Usage

```
.beadsRRBeer(
  object,
  prior.params = list(method = "edgeR", a_pi = 2, b_pi = 300, a_phi = 1.25, b_phi =
      0.1, a_c = 80, b_c = 20, fc = 1),
  beads.args = list(lower = 1),
  jags.params = list(n.chains = 1, n.adapt = 1000, n.iter = 10000, thin = 1, na.rm =
  TRUE, burn.in = 0, post.thin = 1, seed = as.numeric(format(Sys.Date(), "%Y%m%d"))),
  sample.dir = NULL,
  assay.names = c(phi = NULL, phi_Z = "logfc", Z = "prob", c = "sampleInfo", pi =
      "sampleInfo"),
  summarize = TRUE,
  BPPARAM = bpparam()
)
```

Arguments

object	PhIPData object
prior.params	named list of prior parameters
beads.args	named list of parameters supplied to estimating beads-only prior parameters (a_0,b_0)
jags.params	named list of parameters for running MCMC using JAGS
sample.dir	path to temporarily store RDS files for each sample run, if NULL then $[base::tempdir]$ is used to temporarily store MCMC output and cleaned afterwards.
assay.names	named vector indicating where MCMC results should be stored in the PhIPData object
summarize	logical indicating whether to return a PhIPData object.
BPPARAM	[BiocParallel::BiocParallelParam] passed to BiocParallel functions.

Value

vector of process IDs or a PhIPData object

Deadsincugen Function to run the bedas-only round room using eagen	. beadsRREdgeR	Function to run the beads-only round robin using edgeR
--	----------------	--

Description

Each sample is run in comparison to all other beads-only samples to approximate the false positive rate of detecting enrichments.

4 .brewSamples

Usage

```
.beadsRREdgeR(
  object,
  threshold.cpm = 0,
  threshold.prevalence = 0,
  assay.names = c(logfc = "logfc", prob = "prob"),
  BPPARAM = BiocParallel::bpparam()
)
```

Arguments

object A PhIPData object of only beads-only samples.

threshold.cpm CPM threshold to be considered present in a sample threshold.prevalence

proportion of beads-only samples that surpass threshold.cpm.

assay.names named vector specifying the assay names for the log2(fold-change) and exact

test p-values. If the vector is not names, the first and second entries are used as

defaults

BPPARAM [BiocParallel::BiocParallelParam] passed to BiocParallel functions.

Value

vector of process IDs

.brewSamples

Run BEER for all samples

Description

Encapsulated function to run each sample against all beads-only samples. The code is wrapped in this smaller function to (1) modularize the code and (2) make sure the cli output colors don't change.

Usage

```
.brewSamples(
  object,
  sample.id,
  beads.id,
  se.matrix,
  prior.params,
  beads.prior,
  beads.args,
  jags.params,
  tmp.dir,
  BPPARAM
)
```

.checkCounts 5

Arguments

object	PhIPData object
sample.id	vector of sample IDs to iterate over
beads.id	vector of IDs corresponding to beads-only samples
se.matrix	matrix indicating which peptides are clearly enriched
prior.params	list of prior parameters
beads.prior	data frame of beads-only prior parameters
beads.args	named list of parameters supplied to estimating beads-only prior parameters (a_0,b_0)
jags.params	list of JAGS parameters
tmp.dir	directory to store JAGS samples

[BiocParallel::BiocParallelParam] passed to BiocParallel functions.

Value

BPPARAM

vector of process id's for internal checking of whether functions were parallelized correctly.

. checkCounts Function to check that the counts matrix only contains integers	
---	--

Description

Function to check that the counts matrix only contains integers

Usage

```
.checkCounts(object)
```

Arguments

object PhIPData object

Value

nothing if all counts are integers, and error otherwise

6 .edgeRBeads

.checkOverwrite

Function to check whether an assay will be overwritten

Description

If the an assay is not specified (e.g. with NA), then .checkOverwrite() will return FALSE (rather than NA).

Usage

```
.checkOverwrite(object, assay.names)
```

Arguments

object PhIPData object

assay.names character vector of assay names

Value

logical vector indicating whether data in an assay will be overwritten

.edgeRBeads

Estimate edgeR dispersion parameters from the beads-only data

Description

Wrapper function to estimate edgeR dispersion parameters from beads-only samples. Peptides can be pre-filtered based on a minimum read count per million (cpm) and the proportion of beads-only samples that surpass the cpm threshold.

Usage

```
.edgeRBeads(object, threshold.cpm = 0, threshold.prevalence = 0)
```

Arguments

object PhIPData object (can have actual serum samples)
threshold.cpm CPM threshold to be considered present in a sample
threshold.prevalence

proportion of beads-only samples that surpass threshold.cpm.

Value

a DGEList object with common, trended, and tagwise dispersion estimates

.getABEdgeR 7

.getABEdgeR

Derive beta shape parameters using edgeR dispersion estimates

Description

Given a PhIPData object, beads-only shape parameters are estimated by first deriving the peptide-specific edgeR dispersion estimate ϕ^{edgeR} . ϕ^{edgeR} corresponds to the squared coefficient of variation for the proportion of reads pulled for a given peptide. Using ϕ^{edgeR} to derive an estimate of the variance for the proportion of reads pulled by a single peptide, the mean and variance are converted to shape parameters of a beta distribution.

Usage

```
.getABEdgeR(
  object,
  threshold.cpm = 0,
  threshold.prevalence = 0,
  lower = 1,
  upper = Inf
)
```

Arguments

Value

dataframe with rows corresponding to peptides and columns corresponding to estimated shape parameters of the beta distribution.

See Also

```
[.edgeRBeads()] for estimating \phi^{edgeR}
```

8 .getABMLEProp

.getABMLE	Wrapper function to derive MLE estimates of a, b from beads-only samples

Description

Wrapper function to derive MLE estimates of a, b from beads-only samples

Usage

```
.getABMLE(
  object,
  prop.offset = 1e-08,
  optim.method = "default",
  lower = 1,
  upper = Inf
)
```

Arguments

```
object a PhIPData object

prop.offset offset to use when the proportion of reads is 0.

optim.method optimization method passed to [stats::optim].

lower lower bound for the shape parameters.

upper upper bound for the shape parameters.
```

Value

```
a data frame of MLE estimates of a, b
```

See Also

[stats::optim] for available optimization methods

.getABMLEProp	Helper function to derive MLE estimates of a, b from a vector of pro-
	portions

Description

Helper function to derive MLE estimates of a, b from a vector of proportions

.getABMOM 9

Usage

```
.getABMLEProp(
  prop,
  prop.offset = 1e-08,
  optim.method = "default",
  lower = 1,
  upper = Inf
)
```

Arguments

prop vector of proportions.

prop.offset offset to use when the proportion of reads is 0.

optim.method optimization method passed to [stats::optim].

lower lower bound for the shape parameters.

upper upper bound for the shape parameters.

Value

a data frame of MLE estimates of a, b

See Also

[stats::optim] for available optimization methods

 $\begin{tabular}{lll} \it Wrapper function to derive MOM estimates of a, b from beads-only \\ \it samples \end{tabular}$

Description

Wrapper function to derive MOM estimates of a, b from beads-only samples

Usage

```
.getABMOM(
  object,
  offsets = c(mean = 1e-08, var = 1e-08),
  lower = 1,
  upper = Inf,
  ...
)
```

10 .getABMOMProp

Arguments

```
object a PhIPData object.

offsets vector defining the offset to use when the mean and/or variance are zero.

lower lower bound for the shape parameters.

upper upper bound for the shape parameters.

... parameters passed to [base::mean] and [stats::var].
```

Value

a data frame with MOM estimates of a, b

.getABMOMProp Helper function to derive MOM estimates of a, b from a vector of proportions

Description

Helper function to derive MOM estimates of a, b from a vector of proportions

Usage

```
.getABMOMProp(
  prop,
  offsets = c(mean = 1e-08, var = 1e-08),
  lower = 1,
  upper = Inf,
  ...
)
```

Arguments

```
prop vector of proportions.

offsets vector defining the offset to use when the mean and/or variance are zero.

lower lower bound for the shape parameters.

upper upper bound for the shape parameters.

... parameters passed to [base::mean] and [stats::var].
```

Value

```
a data frame with MOM estimates of a, b
```

.guessEnrichedEdgeR 11

.guessEnrichedEdgeR

Guess super-enriched peptides based on edgeR fold-change estimates

Description

Guess super-enriched peptides based on edgeR fold-change estimates

Usage

```
.guessEnrichedEdgeR(object, threshold = 15, fc.name = "logfc")
```

Arguments

object PhIPData object.

threshold minimum estimated fc for a peptide to be considered super-enriched.

fc.name assay name corresponding to the assay that stores the edgeR estimated log2

fold-changes.

Value

logical matrix of the with the same dimensions as object indicating which peptides are considered super-enriched.

.guessEnrichedMLE

Guess enriched peptides based on MLE estimates of the true fold-

change

Description

Guess enriched peptides based on MLE estimates of the true fold-change

Usage

```
.guessEnrichedMLE(object, beads.prior, threshold = 15)
```

Arguments

object PhIPData object.

beads.prior data.frame of prior parameters for beads-only samples.

threshold minimum estimated fc for a peptide to be considered super-enriched.

Value

logical matrix of the with the same dimensions as object indicating which peptides are considered super-enriched.

12 .tidyInputsJAGS

.tidyAssayNames

Clean-up specified assay names

Description

Tidy inputs related to 'assay.names'. Supplies default values for missing parameters and ensures that all required parameters are present.

Usage

```
.tidyAssayNames(assay.names)
```

Arguments

assay.names

named list specifying where to store each assay.

Value

tidied list of assay.names

 $. \verb|tidyInputsJAGS||$

Clean inputs for JAGS parameters

Description

Tidy inputs related to 'jags.params'. Supplies default values for missing parameters and ensures that all required parameters are present.

Usage

```
.tidyInputsJAGS(jags.params)
```

Arguments

jags.params

named list of JAGS parameters

Value

tidied list of JAGS parameters.

.tidyInputsPrior

.tidyInputsPrior	Clean up inputs for prior estimation
------------------	--------------------------------------

Description

Tidy inputs related to 'prior.parameters'. Supplies default values for missing parameters and ensures that all required parameters are present.

Usage

```
.tidyInputsPrior(prior.params, object, beads.args)
```

Arguments

prior.params named list of prior parameters

object PhIPData object

beads.args parameters used to estimate a_0, b_0

Value

tidied list of prior parameters.

.tidyInputsSE Clean up inputs for identifying super-enriched peptides

Description

Tidy inputs related to 'se.params'. Supplies default values for missing parameters and ensures that all required parameters are present.

Usage

```
.tidyInputsSE(se.params, beads.prior)
```

Arguments

se.params named list of parameters for super-enriched estimation

beads.prior data.frame with beads-only parameters

Value

tidied list of parameters for identifying super-enriched peptides.

14 beadsRR

beadsRR	Beads-only round robin

Description

To approximate the false positive rate of each approach, each beads-only sample is run individually against all other samples. For BEER, this means that the sample to be compared is encoded as an actual sample, and prior parameters for beads-only samples are re-estimated. Thus, the beads-only round robin also serves to assess how similar the beads-only samples are to one another.

Usage

```
beadsRR(object, method, BPPARAM = BiocParallel::bpparam(), ...)
```

Arguments

object PhIPData object method one of 'beer' or 'edgeR' specifying which method to use. [BiocParallel::BiocParallelParam] passed to BiocParallel functions. **BPPARAM** parameters passed to the method specific functions. See the *Details* section

below for additional information.

Details

If method == 'beer', then valid parameters include prior.params, beads.args, jags.params, sample.dir, assay.names, and summarize. A description of the first four parameters can be found in brew. summarize is a logical value indicating whether a PhIPData object with the summarized results should be returned. When running beadsRR, summarize typically does not need to be changed.

When method == 'edgeR', threshold.cpm, threshold.prevalence, and assay.names are valid additional parameters that can be supplied to beadsRR. See edgeR for additional details on each of these parameters.

Value

a PhIPData object

See Also

brew for BEER parameters, edgeR for edgeR parameters, and [BiocParallel::BiocParallelParam] for parallelization.

Examples

```
sim_data <- readRDS(system.file("extdata", "sim_data.rds", package = "beer"))</pre>
beadsRR(sim_data, method = "beer")
beadsRR(sim_data, method = "edgeR")
```

brew 15

brew

Bayesian Enrichment Estimation in R (BEER)

Description

Run BEER to estimate posterior probabilities of enrichment, sample-specific attenuation constants, relative fold-changes in comparison to beads-only samples, and proportion of peptides enriched per sample as described in Chen et. al. See *Details* for more information on input parameters.

Usage

Arguments

object	PhIPData object
prior.params	named list of prior parameters
beads.args	named list of parameters supplied to estimating beads-only prior parameters (a_0,b_0)
se.params	named list of parameters specific to identifying clearly enriched peptides
jags.params	named list of parameters for running MCMC using JAGS
sample.dir	path to temporarily store RDS files for each sample run, if NULL then $[base::tempdir]$ is used to temporarily store MCMC output and cleaned afterwards.
assay.names	named vector indicating where MCMC results should be stored in the PhIPData object
beadsRR	logical value specifying whether each beads-only sample should be compared to all other beads-only samples.
BPPARAM	[BiocParallel::BiocParallelParam] passed to BiocParallel functions.

16 brew

Details

prior.params. List of prior parameters. Parameters include,

• method: method used to estimate beads-only prior parameters a_0, b_0. Valid methods include 'custom' or any of the methods specified in getAB. If method = 'custom' is specified, a_0 and b_0 must be included in the list of prior parameters. 'edgeR' is used as the default method for estimating a_0, b_0.

- a_pi, b_pi: prior shape parameters for the proportion of peptides enriched in a sample. Defaults to 2 and 300, respectively.
- a_phi, b_phi: prior shape parameters of the gamma distribution that describe the valid range of enriched-fold changes. The shift is specified by fc. The default values of a_phi and b_phi are 1.25 and 0.1, respectively.
- a_c, b_c: prior shape parameters for the attenuation constant. Default values for a_c and b_c are 80 and 20.
- fc: minimum fold change for an enriched-peptide. fc describes the shift in the gamma distribution.

beads.args. Named list of parameters supplied to getAB. The estimation method used is specified in prior.params, but other valid parameters include lower and upper bounds for elicited parameters. As JAGS recommends that a,b>1 for the beta distribution, beads.args defaults to list(lower = 1).

se.params. Named list of parameters supplied to guessEnriched. By default list(method = 'mle') is used to identify clearly enriched peptides.

jags.params. Named list of parameters for MCMC sampling. By default, BEER only runs one chain with 1,000 adaptation iteration and 10,000 sampling iterations. If unspecified, BEER uses the current date as the seed.

sample.dir. Path specifying where to store the intermediate results. If NULL, then results are stored in the default temporary directory. Otherwise, the MCMC samples for running BEER on each sample is stored as a single RDS file in the specified directory.

assay.names. Named list specifying where to store the point estimates. If NULL, estimates are not added to the PhIPData object. Valid exported estimates include,

- phi: fold-change estimate after marginalizing over the posterior probability of enrichment. By default point estimates are not exported.
- phi_Z: fold-change estimate presuming the peptide is enriched. By default phi_Z estimates are stored in 'logfc' assay.
- Z: posterior probability of enrichment. Estimates are stored in the 'prob' assay by default.
- c: attenuation constant estimates. Stored in 'sampleInfo' by default.
- pi: point estimates for the proportion of peptides enriched in a sample. Stored in 'sampleInfo' by default.

Value

A PhIPData object with BEER results stored in the locations specified by assay. names.

brewOne 17

See Also

[BiocParallel::BiocParallelParam] for subclasses, beadsRR for running each beads-only sample against all remaining samples, getAB for more information about valid parameters for estimating beads-only prior parameters, guessEnriched for more information about how clearly enriched peptides are identified, and [rjags::jags.model] for MCMC sampling parameters.

Examples

```
sim_data <- readRDS(system.file("extdata", "sim_data.rds", package = "beer"))
## Default back-end evaluation
brew(sim_data)
## Serial
brew(sim_data, BPPARAM = BiocParallel::SerialParam())
## Snow
brew(sim_data, BPPARAM = BiocParallel::SnowParam())</pre>
```

brew0ne

Run BEER for one sample

Description

This function is not really for external use. It's exported for parallelization purposes. For more detailed descriptions see brew.

Usage

```
brewOne(
  object,
  sample,
  prior.params,
  n.chains = 1,
  n.adapt = 1000,
  n.iter = 10000,
  thin = 1,
  na.rm = TRUE,
  ...,
  seed = as.numeric(format(Sys.Date(), "%Y%m%d"))
)
```

Arguments

```
object PhIPData object
sample sample name
prior.params vector of prior parameters
```

18 edgeROne

n.chains	number of chains to run
n.adapt	number of iterations to use as burn-in.
n.iter	number of iterations for the MCMC chain to run (after n.adapt)
thin	thinning parameter
na.rm	what to do with NA values (for JAGS)
	extra parameters for JAGS
seed	number/string for reproducibility purposes.

Value

nothing, saves the the results to an RDS in either a temp directory or the specified directory.

Examples

```
sim_data <- readRDS(system.file("extdata", "sim_data.rds", package = "beer"))
beads_prior <- getAB(subsetBeads(sim_data), "edgeR")
brewOne(sim_data, "9", list(
    a_0 = beads_prior[["a_0"]],
    b_0 = beads_prior[["b_0"]],
    a_pi = 2, b_pi = 300,
    a_phi = 1.25, b_phi = 0.1,
    a_c = 80, b_c = 20,
    fc = 1
))</pre>
```

edgeR0ne

Run edgeR for one sample against all the beads-only samples.

Description

This function is not really for external use. It's exported for parallelization purposes. For more detailed descriptions see runEdgeR.

Usage

```
edgeROne(object, sample, beads, common.disp, tagwise.disp, trended.disp)
```

Arguments

object	PhIPData object
sample	sample name of the sample to compare against beads-only samples
beads	sample names for beads-only samples
common.disp	edgeR estimated common dispersion parameter
tagwise.disp	edgeR estimated tagwise dispersion parameter
trended.disp	edgeR estimated trended dispersion parameter

getAB 19

Value

list with sample name, log2 fc estimate, and log10 p-value

Examples

```
sim_data <- readRDS(system.file("extdata", "sim_data.rds", package = "beer"))
beads_disp <- beer:::.edgeRBeads(sim_data)
edgeROne(
    sim_data, "9", colnames(sim_data)[sim_data$group == "beads"],
    beads_disp$common.dispersion, beads_disp$tagwise.disp,
    beads_disp$trended.disp
)</pre>
```

getAB

Estimate beads-only shape parameters

Description

Beta shape parameters are estimated using the proportion of reads-pulled per peptide across the beads-only samples. Currently, only three estimation methods are supported: edgeR, method of moments (MOM), maximum likelihood (MLE). Note that edgeR can only be used on PhIPData objects while MOM and MLE methods can also be applied to vectors of values between 0 and 1. Parameters that can be passed to each method are listed in the details.

Usage

```
getAB(object, method = "mom", ...)
```

Arguments

object a PhIPData object or a vector

method one of c("edgeR", "mle", "mom") designating which method to use to estimate

beads-only prior parameters. MOM is the default method used to estimate shape

parameters.

... parameters passed to specific estimating functions. See details for more infor-

mation

Details

edgeR derived estimates rely on edgeR's peptide-specific dispersion estimates, denoted ϕ^{edgeR} . ϕ^{edgeR} corresponds to the squared coefficient of variation for the proportion of reads pulled for a given peptide. Using ϕ^{edgeR} to derive an estimate of the variance for the proportion of reads pulled by a single peptide, the mean and variance are transformed into shape parameters satisfying the lower and upper bounds. When method = "edgeR", the following additional parameters can be specified.

20 getAB

- threshold.cpm: CPM threshold to be considered present in a sample.
- threshold.prevalence: proportion of beads-only samples that surpass threshold.cpm.
- lower: minimum value of the beta shape parameters.
- upper: maximum value of the beta shape parameters.

Method of Moments (MOM) estimates are derived by transforming the sample mean and variance to shape parameters of the beta distribution. For method = "mom", the following parameters can be adjusted:

- offsets: vector defining the offset to use when the mean and/or variance are zero.
- lower: lower bound for the shape parameters.
- upper: upper bound for the shape parameters.
- ...: parameters passed to [base::mean] and [stats::var].

Maximum Likelihood (MLE) estimates rely on [stats::optim] to derive shape parameters that maximize the likelihood of observed data. By default the L-BFGS-B optimization method is used. Parameters for MLE estimates include:

- prop. offset: offset to use when the proportion of reads is 0.
- optim.method: optimization method passed to [stats::optim].
- lower: lower bound for the shape parameters.
- upper: upper bound for the shape parameters.

Value

a data frame of beta shape parameters where each row corresponds to a peptide.

Examples

```
## PhIPData object
sim_data <- readRDS(system.file("extdata", "sim_data.rds", package = "beer"))
getAB(sim_data, method = "edgeR")
getAB(sim_data, method = "mle")
getAB(sim_data, method = "mom")

## Vector of proportions
prop <- rbeta(100, 2, 8)
getAB(prop, method = "mle")
getAB(prop, method = "mom")</pre>
```

guessEnriched 21

guessEnriched Identifying clearly enriched peptides

Description

As clearly enriched peptides will always have a 100% posterior probability of enrichment, BEER removes these peptides a priori to running the model. Clearly enriched peptides can be identified using edgeR estimated fold-changes or maximum likelihood estimates based on the specified prior parameters. Additional parameters for each method can be found in the details below.

Usage

```
guessEnriched(object, method = "mle", ...)
```

Arguments

object a PhIPData object

method one of "mle" or "edgeR", specifying which method to use to identify clearly

enriched peptides

... additional parameters dependent on the method used. See details for more in-

formation

Details

edgeR. Identification of clearly enriched peptides relies on edgeR fold-change estimates, so edgeR must be run on the PhIPData object beforehand. Additional parameters for identifying clearly enriched peptides based on edgeR estimated fold-changes are listed below:

- object: a PhIPData object.
- threshold: minimum estimated fc for a peptide to be considered super-enriched. The default value is 15.
- fc.name: assay name corresponding to the assay that stores the edgeR estimated log2 foldchanges.

MLE. As the number of reads tends to be quite large, the estimates for the proportion of reads pulled are generally accurate. Clearly enriched peptides are identified by first comparing the observed read count to the expected read count based on the beads-only prior parameters. Peptides with observed read counts larger than 5 times the expected read counts are temporarily labeled as enriched, and attenuation constants are estimated by regressing the observed read counts on the expected read counts for all non-enriched peptides. Using this attenuation constant, peptides with fold-changes above some predefined threshold after adjusting for the attenuation constant are considered enriched. Parameters for identifying clearly enriched peptides using the MLE approach are listed below.

- object: a PhIPData object.
- threshold: minimum estimated fc for a peptide to be considered super-enriched.
- beads.prior: data.frame of prior parameters for beads-only samples.

22 guessInits

Value

a logical matrix of the with the same dimensions as object indicating which peptides are considered super-enriched.

Examples

```
sim_data <- readRDS(system.file("extdata", "sim_data.rds", package = "beer"))
edgeR_out <- runEdgeR(sim_data)

guessEnriched(edgeR_out, method = "edgeR", threshold = 15, fc.name = "logfc")
guessEnriched(edgeR_out,
    method = "mle",
    beads.prior = getAB(edgeR_out, method = "edgeR"),
    threshold = 15
)</pre>
```

guessInits

Derive initial estimates of unknown model parameters

Description

To reduce converge time and to reduce the likelihood of the slice sampler getting stuck, we use maximum likelihood to derive initial estimates for unknown model parameters.

Usage

```
guessInits(object, beads.prior)
```

Arguments

object a PhIPData object

beads.prior a data frame with two columns (named a_0, b_0) containing estimated shape

parameters from beads-only samples.

Details

Briefly initial values are defined as follows:

- 1. theta_guess[i, j] = Y[i, j]/n[j], or the the MLE for theta.
- 2. Z_guess[i, j] = 1 if j is a serum sample, and the observed read count is >2x the expected read count assuming c[j] = 1.
- 3. $pi_guess[j]$ is the mean of column j in Z_guess .
- 4. c_guess[j] is the estimated slope from regressing the observed read counts against the expected read counts (without adjusting for the attenuation constant) for non-enriched peptides only.
- 5. phi_guess[i,j] is the ratio of the observed read counts to the expected read counts multiplied by the attenuation constant.

runEdgeR 23

Value

a list of estimated initial values.

See Also

Methods in [Chen et. al 2022](https://www.biorxiv.org/content/10.1101/2022.01.19.476926v1)

runEdgeR

Run edgeR on PhIP-Seq data

Description

Run edgeR on PhIP-Seq data

Usage

```
runEdgeR(
  object,
  threshold.cpm = 0,
  threshold.prevalence = 0,
  assay.names = c(logfc = "logfc", prob = "prob"),
  beadsRR = FALSE,
  BPPARAM = BiocParallel::bpparam()
)
```

Arguments

object PhIPData object

threshold.cpm CPM threshold to be considered present in a sample

threshold.prevalence

proportion of beads-only samples that surpass threshold.cpm.

assay.names named vector specifying the assay names for the log2(fold-change) and exact

test p-values. If the vector is not names, the first and second entries are used as

defaults

beadsRR logical value specifying whether each beads-only sample should be compared

to all other beads-only samples.

BPPARAM [BiocParallel::BiocParallelParam] passed to BiocParallel functions.

Value

PhIPData object with log2 estimated fold-changes and p-values for enrichment stored in the assays specified by 'assay.names'.

See Also

```
[BiocParallel::BiocParallelParam], beadsRR
```

24 summarizeRun

Examples

```
sim_data <- readRDS(system.file("extdata", "sim_data.rds", package = "beer"))
## Default back-end evaluation
runEdgeR(sim_data)
## Serial
runEdgeR(sim_data, BPPARAM = BiocParallel::SerialParam())
## Snow
runEdgeR(sim_data, BPPARAM = BiocParallel::SnowParam())</pre>
```

summarizeRun

Summarize MCMC chain and return point estimates for BEER parameters

Description

Posterior means are used as point estimates for c,π,ϕ , and Z. As super-enriched peptides are tossed out before MCMC sampling, super-enriched peptides return NA for the ϕ and Z point estimates. Indices corresponding to a particular peptide in the MCMC sampler are mapped back to the original peptide names.

Usage

Arguments

```
object a PhIPData object

jags.files list of files containing MCMC sampling results

se.matrix logical matrix indicating which peptides were identified as super-enriched peptides

burn.in number of iterations to be burned

post.thin thinning parameter

assay.names named vector of specifying where to store point estimates

BPPARAM [BiocParallel::BiocParallelParam] passed to BiocParallel functions.
```

summarizeRunOne 25

Value

PhIPData object with point estimates stored in the assays specified by 'assay.names'.

summarizeRunOne	Derive point estimates for c, pi, phi, and Z for a particular sample

Description

Posterior means are used as point estimates for c,π,ϕ , and Z. As super-enriched peptides are tossed out before MCMC sampling, super-enriched peptides return NA for the ϕ and Z point estimates. Indices corresponding to a particular peptide in the MCMC sampler are mapped back to the original peptide names.

Usage

```
summarizeRunOne(object, file, se.matrix, burn.in = 0, post.thin = 1)
```

Arguments

object	a PhIPData object
file	path to rds file
se.matrix	logical matrix indicating which peptides were identified as super-enriched peptides
burn.in	number of iterations to be burned

post.thin thinning parameter

Value

list of point estimates for c, pi, phi and Z

Index

```
.beadsRRBeer, 2
.beadsRREdgeR, 3
. \, {\tt brewSamples}, \, {\tt 4}
.checkCounts, 5
.checkOverwrite, 6
.edgeRBeads, 6
.getABEdgeR, 7
. \mathtt{getABMLE}, 8
.getABMLEProp, 8
.getABMOM, 9
.getABMOMProp, 10
.guessEnrichedEdgeR, 11
.guessEnrichedMLE, 11
.\, \texttt{tidyAssayNames}, \, \textcolor{red}{12}
.tidyInputsJAGS, 12
.tidyInputsPrior, 13
.tidyInputsSE, 13
beadsRR, 14, 17, 23
brew, 14, 15, 17
brewOne, 17
edgeR, 14, 21
edgeR0ne, 18
getAB, 16, 17, 19
guessEnriched, 16, 17, 21
guessInits, 22
PhIPData, 6-8, 10, 11, 18, 19, 21-25
runEdgeR, 18, 23
summarizeRun, 24
summarizeRunOne, 25
```